Package insert and How to use guide

4N6FLOQSwabs™ Genetics

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4N6FLOQSwabs™ and NAO™ Basket are Copan Italia trademarks.

FLOQSwabs™ are PATENTED WORLDWIDE No. PCT/EP2004/003392 and in the following countries: European Patent #1608268, US Patent #US8114027B2, Canadian Patent #2515205, Japanese Patent #2007-523663, Australian Patent #2004226798, New Zealand Patent #541560, Chinese Patent #101103931.

EO TREATED

ENGLISH

4N6FLOQSwabs™ Genetics

Product Insert & How to Use Guide

INTENDED USE

The 4N6FLOQSwabs™ Genetics product line is intended for the collection and transport of swab samples for the detection and the analysis of human DNA for Genetics applications such as human identification testing and reference samples for paternity testing.

Not intended for diagnostic use.

SUMMARY AND PRINCIPLES

One of the routine procedures for genetics applications involves the collection and transportation of swab samples (for example buccal swab (1), vaginal swabs cells, saliva, semen (2), blood) for human DNA profiling. This can be accomplished using 4N6FLOQSwabs™ Genetics which are certified free of amplifiable human DNA, and detectable DNase and RNase.

4N6FLOQSwabs™ Genetics can be used to collect and transport specimens for DNA profiling and are available in two main formats:

- 4N6FLOQSwabs™ Genetics line, that can be used for the collection of specimens for different collection sites and sample types.

Samples collected with the 4N6FLOQSwabs™ Genetics product line without an Active Drying Agent should be transported, after proper drying, directly to the laboratory. Proper swab drying avoids contaminant overgrowth optimizing DNA stabilization. Swab samples should be stored at room temperature (20 – 25°C) or frozen for up to 6 month (3).

- 4N6FLOQSwabs™ with an Active Drying System

4N6FLOQSwabs™ Genetics line with an Active Drying System does not need drying and samples should be stored at room temperature (20 – 25°C) or frozen for up to 6 month (ref. 3).

REAGENTS

Some product codes are available with Active Drying System, which consists of a desiccant contained in the tube stopper to absorb the humidity of the collected sample.

LIMITATIONS

- 1) Not for diagnostic use.
- 2) 4N6FLOQSwabs™ Genetics product line is intended for the collection and transport of swab samples for the detection of human DNA, for Genetics applications such as human identification testing, and reference samples for paternity testing.
- 3) The amount of specimen collected for DNA analysis significantly varies from donor to donor and in relation to collection execution mode. For the proper collection procedure, please refer to Specimen Collection instructions.
- 4) The product codes supplied in peel pouch format must NOT be reinserted into the peel pouch or envelopes after sample collection. Appropriate transport systems should be adopted.
- 5) If codes W/O Active Drying System are used for collection, the swab must be air-dried, before reclosing into its own tube or cuvette (it should take approximately 1 hour). Reclosing a wet swab into a tube or a cuvette could result in nucleic acids degradation due to microbial overgrowth.
- 6) If the sample has been frozen bring it back at room temperature before use
- 7) The use of handling and testing procedures different from the ones described in this package insert may affect product performance and analysis results.

- 8) The use of the product in combination with extraction and amplification kits must be always qualified by the user.

PRODUCT STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 2 – 30°C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box, and on each individual pack.


PRODUCT DETERIORATION

4N6FLOQSwabs™ *Genetics* should not be used if 1) there is evidence of damage or contamination to the product, 2) the expiration date is passed, 3) the swab package is open and/or the security seal is compromised.



MATERIALS SUPPLIED

4N6FLOQSwabs™ *Genetics* line are available in different formats shown below.


Regular size tip 4N6FLOQSwabs™: individually wrapped in peelpouch (product code 4520CS01)


Regular size tip 4N6FLOQSwabs™ individually wrapped in peelpouch with a 2ml cuvette with and without evaporation duct (product codes 4520CF and 4520CA)


Regular size tip 4N6FLOQSwabs™ individually wrapped in peelpouch with a 2ml cuvette and NAO™ basket (product code 4103CS01)

Regular size tip 4N6FLOQSwabs™, single or dual, in a plain dry tube (product codes 4508C, 4511C and 40U001D)

Regular size tip 4N6FLOQSwabs™ or Buccal shape tip 4N6FLOQSwabs™, in a shorter tube format (tube length 109mm) with Active Drying System (product codes 4504C and 4500C), 4N6FLOQSwabs™ sublingual shape in a shorter tube format (tube length 109mm) with Active Drying System (product code 40U022D).

All 4N6FLOQSwabs™ *Genetics*, with the only exception of the buccal shape tip (code 4500C), are provided with a molded breakpoint which allows easy breakage of the swab tip end into a cuvette and/or NAO™ basket.

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for nucleic acids extractions, amplifications and DNA analysis. Refer to laboratory reference manuals for nucleic acid amplification and analysis.

DIRECTIONS FOR USE

Copan 4N6FLOQSwabs™ *Genetics* are available in the product configurations indicated in the table below:

Product Code	Description	Packaging	Sample
4520CS01	4N6FLOQSwabs™ <i>Genetics</i> , Regular size tip, with breakpoint at 20mm, individually wrapped in peelpouch	100 units/box 10x100 units/carton	Buccal Swabs (Saliva, cells, others)
4520CA	4N6FLOQSwabs™ <i>Genetics</i> , Regular size tip, with breakpoint at 20mm, plus 2 ml cuvette individually wrapped in peelpouch	100 units/box 6x100 units/carton	
4520CF	4N6FLOQSwabs™ Regular size tip breakpoint at 20mm plus 2 ml cuvette with evaporation duct individually wrapped in peelpouch	100 units/box 6x100 units/carton	
4103CS01	4N6FLOQSwabs™ <i>Genetics</i> , Regular size tip, with breakpoint at 20mm, plus NAO™ Basket and 2ml cuvette, individually wrapped in peelpouch	100 units/box 6x100 units/carton	Vaginal Swabs (Semen, Blood, others)
4508C	4N6FLOQSwabs™ Regular size tip, with breakpoint at 20mm, in dry tube	100 units/box 10x100 units/carton	
4511C	4N6FLOQSwabs™ Dual Swab Regular size tip, breakpoint at 20mm in dry tube	100 units/box 6x100 units/carton	Buccal Swabs (Saliva, cells, others)
40U001D	4N6FLOQSwabs™ <i>Genetics</i> , regular size tip with stopper without breakpoint, in dry tube	100 units/box 6x100 units/carton	
4504C	4N6FLOQSwabs™ Regular size tip, with breakpoint at 20mm, in 109mm long tube with Active Drying System	100 units/box 6x100 units/carton	
4500C	4N6FLOQSwabs™ Buccal Shape tip, unbreakable, in 109 mm long tube with Active Drying System	100 units/box 6x100 units/carton	
40U022D	4N6FLOQSwabs™ <i>Genetics</i> with Active Drying System, Subungual shape, breaking point at 20 mm, in 109mm long tube	100 units/box 6x100 units/carton	Subungual swab (blood, cells, others)

Additional product codes may be available. For updates please refer to our website: www.copangroup.com

Specimen Collection

Proper specimen collection from the donor is extremely critical for successful DNA analysis. For specific guidance regarding specimen collection procedures, follow procedures below or internal collecting procedures.

Gloves and protective clothing and eyewear should be worn when collecting and handling specimens.

During sample collection when handling the swab applicator, the operator must hold the swab from the handle or the upper part of the swab shaft.

NOTE: Do not use excessive force, pressure or bending when collecting swab samples from humans as this may result in accidental breakage of the swab shaft. Swab shafts often exhibit diameter changes to facilitate different sampling requirements. Swab shaft have a molded breakpoint point close to the tip (with the exception of product code 4500C) designed for intentional breakage of the swab into a transport tube.


NOTE: The donor must not eat, drink (apart from water) or brush teeth or chew gum one hour before sample collection. The swab must NOT be reinserted into the peelpouch after sample collection.

NOTE: In particular, for the product codes under the format "4N6FLOQSwabs™ with an Active Drying System", the sample collection must be performed immediately once the aluminum pouch is opened.


- 1) Wear gloves and open the peel pouch. While holding the swab from the handle or from the upper part of the shaft, remove it from the pouch.
- 2) Holding the swab by the handle, insert the swab inside the mouth over the tongue and wet the swab with saliva.

- 3) Move the swab to the side of the mouth and firmly roll the swab 10 times against the cheek, close to the gingival wall. Use the index and medium finger to press on the cheek while collecting the swab.
- 4) Repeat step n° 3 by transferring the swab on the other cheek


STEP 2

STEP 3

STEP 4

For Genetics 4N6FLOQSwabs™ codes 4520CS01, 4520CA, 4520CF, 4103CS01 (individually wrapped in peel pouch):

5) After collection, allow the swab to air dry for 1 hour, in an environment protected by contamination.


6/A) For Product Code 4520CS01: using the breaking point of the swab shaft, break the swab in a cuvette (not provided) to transport to the testing laboratory or process as per laboratory procedures.

6/B) For Product Codes 4520CA, 4520CF : using the breaking point close to the tip of the swab shaft, break the swab in the cuvette included in the pouch, close the cap of the cuvette to transport to the testing laboratory.

For Product Code 4103CS01: using the breaking point close to the tip of the swab shaft, break the swab in the NAO™ basket included in the pouch (previously inserted

into the relative provided cuvette), close the cap of the tube of NAO™ basket to transport to the testing laboratory.

STEP 6/B


For Genetics 4N6FLOQSwabs™ codes 4508C, 4511C and 40U001D (in dry tube)

5) After collection, allow the swab to air dry for 1 hour, in an environment protected by contamination.

6) Replace the swab into the plastic tube and close tightly.

7) Write the patient ID and date.


STEP 6

For Genetics 4N6FLOQSwabs™ codes 4504C, 4500C (with active drying system)

5) After collection, replace swab into the plastic tube and close tightly.


STEP 5

Vaginal swab specimen collection

For product codes 4520CS01, 4520CA, 4520CF, 4103CS01, 4508C, and 4511C

1. Open the peel pouch.
2. Remove the 4N6FLOQSwabs™ from the peel pouch. Collect the sample from the patient with the 4N6FLOQSwabs™ as per the following instructions.
NOTE: Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shafts often exhibit diameter changes to facilitate different sampling requirements. Swab shaft have a molded breakpoint close to the tip designed for intentional breakage of the swab into a transport tube.
3. Insert the swab about 5 cm (2 inches) into the opening of the vagina.
4. Gently rotate the swab 3 times against the sides of the vagina. Withdraw the swab carefully.


For Genetics 4N6FLOQSwabs™ codes 4520CS01:

- 5) using the breaking point of the swab shaft, break the swab in a cuvette (not provided) to transport to the testing laboratory or process as per laboratory procedures.


For Genetics 4N6FLOQSwabs™ codes 4520CA, 4520CF:

- 5) using the breaking point close to the tip of the swab shaft, break the swab in the cuvette included in the pouch, close the cap of the cuvette to transport to the testing laboratory. For Product Code 4103CS01: using the breaking point of the swab shaft, break the swab in the NAO™ basket included in the pouch (previously inserted into the relative provided cuvette), close the cap of the tube+NAO™ basket to transport to the testing laboratory.

For Genetics 4N6FLOQSwabs™ codes 4508C, 4511C (in dry tube)

- 5) After collection, allow the swab to air dry for 1 hour, in an environment protected by contamination.
- 6) Replace the swab into the plastic tube and close tightly


- 7) Write the patient ID and date.

8) Write patient information on the tube label or apply identification label. Send the sample to the test laboratory.

Subungueal swab specimen collection

Product code 40U022D:

- 1) Wear gloves and open the aluminum pouch. While holding tightly the cap, open the tube.
- 2) Place the tip of the swab under fingernail and gently move the swab along the length of the nail. Repeat the collection in the opposite direction as shown in the picture.
- 3) After collection, replace the swab into the plastic tube and close tightly.

Processing 4N6FLOQSwabs™ Genetics in the laboratory

Once the swab sample arrive to the laboratory, they can be processed immediately, or they can be kept at room temperature for short period of time (e.g. 1 month) or they can be frozen for long term storage (e.g. 6 month). The storage at +4°C of the swabs after collection is not recommended, because the humidity can give rise to accelerated DNA degradation during time.

Gloves and protective clothing and eyewear should be worn when collecting and handling forensic specimens.

When working with molecular methods care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are essential to prevent amplicon carry-over.

Nucleic acids extraction:

During sample processing when handling the 4N6FLOQSwab™ applicator, the operator must not touch the area below the breakpoint that is the area from the line to the tip of the nylon flocked swab, as this will lead to contamination of the applicator shaft invalidating the test results.

- 1) For Product Codes 4520CA, 4520CF: open the cuvette containing the 4N6FLOQSwabs™ tip, add directly into the cuvette the appropriate amount of lysing buffer, as per package insert of the extraction kit in use. Continue as per extraction and amplification kits procedures.
- 2) For Product Code 4103CS01: (In the case of swabs in NAO basket + cuvette): open the cap of the NAO basket containing the swabs, and add directly into the NAO basket the appropriate amount of lysing buffer, as per package insert of the extraction kit in use. The lysis step can proceed, keeping the swabs into NAO basket, as per package insert of the extraction kit in use. NOTE 1: the maximum volume of lysis buffer that can be added to NAO basket is 700 ul, the minimum volume recommended is 300 ul, in order to cover the major part of the swab surface. NOTE 2: for lysis temperatures of 70°C or more, some drops of liquid could come through the valves of NAO basket into the below cuvette. This doesn't represent a problem, if the swabs stay immersed into lysis buffer. At the end of the lysis step, centrifuge the NAO basket + cuvette containing the swab at 10000xg for 1 minute, in order to recover all the liquid containing the lysate sample in the below cuvette. Open the cap of the tube, discard the NAO™ basket containing the dry swab and proceed with the cuvette containing the lysate sample, as per package insert of the extraction kit in use.
- 3) For Product Codes 4508C, 4511C, 4504C, 4500C, 40U001D (4N6FLOQSwabs™ in tube) and 4520CS01: Using the breaking point of the swab shaft, transfer the swab to an extraction tube as per laboratory SOP. Add the appropriate amount of lysing buffer, as per package insert of the extraction kit in use. Continue as per extraction and amplification kits procedures. All Genetic swabs with the exception of code 4500C have a breaking point compatible with NAO basket so they can be used in association during the lysing step: in this case, follow the instructions given at previous point n°2.

4N6FLOQSwabs™ Genetics have been internally tested with the following extraction methods: ChargeSwitch Forensic DNA Purification Kit by Invitrogen, and Chelix 100 by Biorad (Ref. 2), PrepFiler Express by Life Technologies (Ref. 3); Nucleo Spin Tissue by Macherey Nagel, DNA IQ system by Promega and InstaGene matrix by Bio-Rad (Ref. 4).

Extracted DNA has been tested with the following PCR kits: and AmpFtSTR® Yfiler™ PCR Amplification Kit by Life Technologies (Ref. 2); Identifier Plus by Life Technologies and Quantifiler Trio Human DNA Quantification kit by Life Technologies (Ref. 3); Human DNA Quantification kit by Life Technologies (Ref. 4).

QUALITY CONTROL

Genetics 4N6FLOQSwabs™ are guaranteed for following features:

Analyte release: the swab is inoculated with a bacterial suspension, transferred into liquid phase and check for the complete analyte release.

≥70% of analyte must be released by the swab.

- **DNase free:** Following internal SOP, 4N6FLOQSwabs™ are saturated with a solution containing standard DNA ladder. DNA from solution exposed to the swabs must be intact.
- **RNase free:** Following internal SOP, 4N6FLOQSwabs™ are saturated with a solution containing standard RNA ladder. RNA from solution exposed to the swabs must be intact.
- **Human DNA free:** Following internal SOP 4N6FLOQSwabs™ are tested to be free of any amplifiable human DNA with the Identifier® Plus STR kit.
- **Active Drying System ability:** Following internal SOP, 4N6FLOQSwabs™ with Active Drying System are inoculated with a known amount of osmotic water. The swabs, reclosed in their own tubes, must be completely dry after 24 hours at room temperature.

WARNING AND PRECAUTIONS

1. Not for diagnostic use.
2. Not suitable for any other application than intended use.
3. To be used only by adequately trained or qualified personnel.
4. Directions for use must be followed carefully.

- Before transport, where there is no active drying system, make sure that proper drying procedure has been followed. Transport of a not dried sample may result in DNA degradation.
- Product is ready to use; it is not necessary to pretreat prior sample collection.
- All specimens and materials used to process them should be treated in human DNA free condition, to avoid post collection contamination by personnel, that could lead to false results: wear disposable gloves, mask, cap and other protection commensurate with universal precautions not to contaminate the samples.
- This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
- Do not re-pack.
- The use of this product in association with a rapid kit or with instrumentation should be previously validated by the user.
- Do not use if the swab is visibly damaged (i.e., if the swab tip or swab shaft is broken).
- Do not use excessive force or pressure when collecting swab samples from patients as this may result in breakage of the swab shaft.
- Dispose of used products, waste and specimens in accordance with local regulations.
- It must be assumed that all specimens may contain infectious micro-organisms; therefore all specimens must be handled and disposed following appropriate precautions.

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EO TREATED

BIBLIOGRAPHY

1. **VALUTAZIONE DELL'EFFICIENZA DELLA METODICA DI PRELIEVO SALIVARE MEDIANTE TAMPONI COPAN 4N6FLOQSwabs™ genetics** Corato S., Ponzano E., Giuliodori A., Rodriguez D., Caenazzo L. Presented at XXIV Convegno Nazionale Genetisti Forensi Italiani in Pavia on 20-22 Settembre 2012
2. **Enhanced semen elution from nylon 4N6FLOQSwabs™ prior DNA analysis** P. MarDnez, MG.Blanco, G.Vallejo National Institute of Toxicology & Forensic Science INTCF José Echegaray 4, 28232 Las Rozas Madrid SPAIN Presented at International Symposium on Human Identification in Phoenix on 29 September-2 October 2014.
3. **Genetics 4N6FLOQSwabs™ are supporting DNA preservation from buccal swabs for long term storage** A. Squassina S. Castriciano Copan Italia S.p.A. Brescia, Italy. Presented at International Society for Forensic Genetics in Krakow on 31 August – 5 September 2015
4. **FORENSIC FLOCKED COLLECTION DEVICES TO MAXIMIZE CRIME SCENE SAMPLE PROCUREMENT** Squassina, A., Poli, A., Gervasoni, A., Castriciano, S., Paroni, R., Triva, D. Copan Italia S.p.A. Brescia, Italy Presented at International Society for Forensic Genetics in Vienna on 29 August – 3 September 2011

Tabella dei Simboli / Table of Symbols / Tableau des pictogrammes / Tabelle der verwendeten Symbole / Tabla de símbolos / Tabela de Símbolos / Tabela simbola

Simbolo / Symbol / Pictogramme / Símbolo / Symbol / Símbol	Significato / Meaning / Signification / Bedeutung / Significado / Význam / Značenje
	Fabbricante / Manufacturer / Fabricant / Hersteller / Fabricante / Výrobce / Proizvođač
CE 0123	Numero di identificazione dell'organismo notificato/ Identification number of notified body/ Identification de l'organisme notifié/ Identifizierung der benannten Stelle/ Identificación del organismo notificado/ Identificação do organismo notificado / Identifikační číslo notifikované osoby / Identifikacioni broj notifikovanog tela
STERILE EO	Sterilizzato usando ossido di etilene/ Sterilized using ethylene oxide/ Stérilisé à l'aide d'oxyde d'éthylène/ Sterilisiert mit Äthylenoxid/ Esterilizado usando óxido de etileno/ Esterilizado por óxido de etileno / Sterilizováno ethylenoxidem / Sterilisano etilen oksidom
STERILE R	Sterilizzato usando radiazioni ionizzanti/ Sterilized using irradiation/ Stérilisé à l'aide de radiations ionisantes/Sterilisiert mit ionisierenden Strahlungen/ Esterilizado usando radiaciones ionizantes// Sterilisiert ved hjælp av ionisert stråling/Sterilizují ozařování / Sterilisano jonizujućim zračenjem
	Non riutilizzare / Do not reuse / Ne pas réutiliser / Nicht zur Wiederverwendung / No reutilizar / Não reutilizar / Nepoužívat znovu / Ne koristite ponovo
REF	Numero di catalogo / Catalogue number / Numéro de catalogue / Bestellnummer / Número de catálogo / Katalogové číslo / Broj kataloga
	Limiti di temperatura / Temperature limits / Limites de température / Temperaturbegrenzung / Limites de temperatura / Limites de temperatura / Teplotní limity / Ograničenja temperature
	Utilizzare entro / Use before / À utiliser avant / Verwendbar bis / Usar antes de / Usar antes de / Použit do / Upotrebiti do
	Consultare le istruzioni per l'uso / Consult the instructions for use / Consulter les instructions d'utilisation / Gebrauchsanweisung beachten / Consultar las instrucciones de uso / Consultar as instruções de utilização / Prostudovat návod k použití / Pročítajte uputstvo za upotrebu
	Strappare per aprire/ Peel/ Décoller/ Abziehen/ Desprender/ Destacável / Odtřhnout / Odlepite
LOT	Codice del lotto (partita) / Batch code (lot) / Code du lot (lot) / Chargenbezeichnung / Código del lote (Lote) / Código do lote (remessa) / Kód šarže (dávky) / Broj partije
	Contenuto sufficiente per <n> test / Contents sufficient for <n> tests / Contenu suffisant pour <n> tests / Inhalt ausreichend für <n> Prüfungen / Contenido suficiente para <n> pruebas / Conteúdo suficiente para <n> testes / Obsah stačí pro <n> testů / Sadržaj je dovoljan za <n> testova
	Non utilizzare in caso di confezionamento danneggiato / Do not use if the packaging is damaged / Ne pas utiliser si l'emballage est endommagé / Nicht verwenden, wenn Verpackung beschädigt / No utilizar si el envase está dañado / Não utilizar se a embalagem apresentar danos / Nepoužívat v případě poškozeného obalu/Ne koristite ukoliko je pakovanje oštećeno

Copan Flock Technologies Srl
Via F. Perotti, 18
25125 - Brescia, Italy

Copan Flock Technologies Srl
Via Perotti 18
25125 Brescia Italy
Tel +39 030 3666100
Fax +39 030 2659932
Email: info@copanflock.com
Website: www.copangroup.com

North American Distributor:
Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562 USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.com
Website: www.copanusa.com

This page is intentionally left blank

This page is intentionally left blank