

exposed to CyMol or UTM for 30 minutes, inoculated into shell vial culture

and stained after 48 hours incubation at 37°C

IMPROVING PRE-ANALYTIC COLLECTION SYSTEMS: INACTIVATION AND PRESERVATION OF INFLUENZA A FOR RAPID TESTING

<u>K. Luinstra¹</u>, S. Castriciano², M. Ackerman¹, A. Petrich¹, S. Chong¹, J.B.Mahony¹, and M. Smieja¹. St. Joseph's Healthcare, Hamilton, ON, Canada and Copan Italia



| Abstract Objective: An alternative sample collection system that would render a virus non-infectious but permit rapid identification by direct immunofluorescence (DFA) or molecular methods would be beneficial during a pandemic. In this study, the CVMol collection system was evaluated for its ability to inactivate influenza A and poreserve the | Methods Influenza A RNA Stability and Recovery •Viral lysate was adsorbed onto duplicate swabs and placed in CyMol and UTM Collection tubes •The stability and recovery of influenza A nucleic acid was assessed in duplicate after 1, 7, 14 and 21 days exposure to CyMol or UTM at -20°C, 4°C, room temperature and 37°C. •200 uL aliquots of mocked collected samples were extracted by easyMAG (bioMérieux) and eluted in 60 uL •5 uL of purified NA was tested by a quantitative Influenza A RT-PCR (Matrix gene, CDC) on the Roche LightCycler 2.0. | | Average RNA Copies Recovered/5uL | | | | | |
|--|---|--|---------------------------------------|--|----------------------------------|--|--|--|
| sample for molecular testing. Methods: (VMC) (Copan, Italia) is an alcohol-based medium that preserves cells for DFA. Flocked nacopharyngeal avaba (NFS) collected in universal transport media (UTM) were compared to NFS collected in CyMeI. Aliquots of Influenza A viral lystate were absorbed onto duplicate flocked avabas and placed into the UTM and CyMeI collection moched samples into R-mix shell vial culture followed by immunofluorescent staining at 48 hours. That Avabas assessed for 51 miluneza A subtypes. The stability and recovery of influenza A nuclei cald (NA) was also assessed after 1, 7, 14 and 21 dyns at 44°, 20°C, cron temperature (RY) and 37°C. Aliquots of mocked specimens were cartracted by searMAG (boldericas) and S ui ca purprised NA tested by quantitative R-FRC on the Roche Results: In contrast to the UTM, the influenza A virus was inactivated on exposure to the CyMeI collection medium and failed to grow in shell vial culture. Influenza A RNA levels were stable for up to 14 days in both the UTM and CyMeI collection systems at -20°C, 4°C and RT. The stability of Influenza A ANdA declined in both systems at 37°C. Conclusions: The Copan CyMeI andium inactivated influenza A rivel wise stabilizing RNA for molecular testing for up andemic stability. Cell and RT. CyMeI is a potential alternative for safe sample collection during an influenza pandemic stability. | | | Sample UTM-20 Cymol-20 UTM-4 | 1 day 2.13E+05 1.92E+05 2.09E+05 | 2.49E+05 | 14 days 2.62E+05 1.45E+05 2.29E+05 | 21 days 2.33E+05 2.59E+05 2.11E+05 | |
| | Results Inactivation of Infectivity | | Cymol-4 UTM-RT | 1.80E+05 2.05E+05 1.18E+05 | 2.07E+05 1.95E+05 1.10E+05 | 2.31E+05 1.86E+05 7.18E+04 | 2.47E+05 1.21E+04 4.24E+04 | |
| Objective | •In contrast to the UTM, Influenza A/ Victoria/3/75 H3N2 was in Medium at RT and failed to grow in shell vial culture. | . , | Cymol-RT UTM-37 Cymol-37 | 1.69E+05 5.80E+04 | 5.10E+04 | 1.08E+04 1.91E+03 | 1.19E+03 | |
| The ability of the CyMol collection system (Copan, Italy) to inactivate influenza A and preserve the sample for molecular testing was evaluated in this study. | 4 other influenza A subtypes (H1, H8, H10 and H15) were also inactivated by a 30 minute exposure to the CyMol Collection Medium at RT and failed to grow in shell vial culture. Influenza A virus inactivation was observed to occur after CyMol Collection Medium exposure at time zero. | | Average Crossing Points | | | | | |
| | Influenza A RNA Stability and Recovery | | Sample | 1 day | 7 days | 14 days | 21 days | |
| Methods | •Influenza A RNA levels were stable for up to 14 days in both the UTM and CyMol collection systems at -20°C, 4°C and RT. •The stability of Influenza A RNA declined in both systems at 37°C. | | UTM-20 Cymol-20 | 23.86 24.35 | 23.82 23.88 | 23.75 24.62 | 23.90 23.75 | |
| Influenza A Virus Inactivation | UTM vs CyMol RNA Recovery at -20 C | UTM vs CyMol RNA Recovery at 4 C | UTM-4 | 23.94 | 23.67 | 23.94 | 24.05 | |
| Viral lysate (Influenza A/Victoria/3/1975 (H3N2)) was diluted in pooled | OT M VS CYMOL RIVA RECOVERY at -20 C | UIM VS CYMOIRNA Recovery at 4 C | Cymol-4 | 24.44 | 24.07 | 23.93 | 23.83 | |
| negative respiratory samples to obtain a dilution representative of an | 1000000 | 1000000 | UTM-RT | 23.98 | 24.08 | 24.30 | 28.38 | |
| influenza A positive sample (10 ⁴ -10 ⁵ copies/5uL; Crossing Point 25-27) | 100000 | | Cymol-RT | 25.06 | 24.99 | 25.64 | 26.41 | |
| 50 uL of the diluted viral lysate was adsorbed onto a flocked nasopharyngeal swab and then placed in a CyMol or UTM collection tube | | 100000 <u>a Cymol-4</u> | UTM-37 | 24.25 | 25.98 | 28.40 | 31.55 | |
| •After 30 minutes exposure to the media, the mock sample collections were vortexed for 20 seconds, swabs removed and collection media spun for 10 | 1 day 7 days 14 days 21 days Storage Time | 1 day 7 days 14 days 21 days Storage Time | Cymol-37 | 26.09 | 27.28 | 31.47 | 36.13 | |
| minutes at 2000 rpm •200 uL of the neat supernatant and 200 uL of 1:10 diluted supernatant were inoculated in duplicate into R-Mix shell vials with coverslips | UTM vs CyMol RNA Recovery at RT UTM vs CyMol RNA Recovery at 37 C | | Conclusions | | | | | |
| •Immunofluorescent staining of the shell vials occurred after 48 hours incubation at 37°C •Virus inactivation was monitored by shell vial culture at time zero and after 5, 10 and 20 minute exposures to CyMol or UTM Collection Media •Viral lysate from 4 other influenza A subtypes (H1, H8, H10 and H15) was | | | | The Copan CyMol medium inactivates influenza A virus infectiv on exposure to the media. Influenza A RNA was stabilized for molecular testing for up to 3 days at -20°C, 4⁰C and RT. CyMol is a potential alternative for safe sample collection during | | | | |

1 day

7 days

Storage Time

14 days 21 days

Storage Time

•CyMol is a potential alternative for safe sample collection during an influenza pandemic situation.