

Background:




Automated Antimicrobial Susceptibility Systems provide reliable results improving clinical outcomes although the microbial suspension preparation is still completely manual. Colibrí® is a new automated system able to pick colonies from plates, re-suspend them in saline solution, check the turbidity of the suspension with the on-board nephelometer and adjust it if required.

Objective:

The aim of the study was to compare manual microbial suspensions preparation for Antimicrobial Susceptibility Testing (AST) to Colibrí® automated preparation.

Materials/methods:

Twenty-Two Gram negative and twenty-two Gram positive clinical strains were streaked on CPS®Elite (Biomérieux) agar plates using Copan WASP® and incubated in WASPLab® automation. Isolated colonies were manually picked and suspended in saline solution to obtain a 0.5 McFarland microbial concentration measured with the Biomérieux DensiCHECK Plus. The same plate was then processed by Colibrí®: selected bacteria colonies were automatically picked by the **pipettor** and transferred into a saline solution tube. The turbidity of the microbial suspension was verified by the **on-board nephelometer**. Manual and automated prepared microbial suspensions were used to perform the AST on the Biomérieux **Vitek2** (version 07.01) as per manufacturer instructions. Gram negative (*E. coli* (n=6), *E. cloacae* (n=5), *P. aeruginosa* (n=3), *S. marcescens* (n=2), *P. mirabilis* (n=3), *K. oxytoca* (n=1), *K. pneumoniae* (n=1), *C. freundii* (n=1)) were processed with AST-N204 card and Gram positive (*S. epidermidis* (n=6), *S. aureus* (n=8), *E. faecalis* (n=4), *E. faecium* (n=3), *E. gallinarum* (n=1)) were processed with AST-P592 card. The viable count of 0.5 McFarland suspensions was performed using the serial dilution method on both the manual and automated preparations.



Colibrí® is able to:

- Pick colonies using disposable tips (1)
- Release bacteria colonies in a saline solution tube (2)
- Check the microbial suspension turbidity using the on-board nephelometer (3)
- Adjust the microbial suspension adding colonies or saline to obtain a 0.5 MacFarland
- Prepare the microbial dilution needed for processing with AST system
- Prepare the purity plate from the microbial suspension (4)

Results:

The results were interpreted with Biomérieux Expert System EUCAST 2015:

- 318 antibiotic combinations were tested for Gram positive: 1 minor error (0.31%) and 1 major error (0.31%) were found.
- 306 antibiotic combinations were tested for Gram negative: 9 minor error (2.94%) and 1 very major error (0.32%) were found.
- Only 1 discrepant result could be related to the viable count where the manual preparation was 5-fold higher than Colibrí® microbial preparation.

Five out of the 44 strains tested were defined as MDR organisms according to Magiorakos *et al* 2011 ⁽¹⁾.

	Agreement	Minor error	Major error	Very major error
Gram positive (n=318)	n=316 (99.38%)	n=1 (0.31%)	n=1 (0.31%)	n=0
Gram negative (n=306)	n=296 (96.74%)	n=9 (2.94%)	n=0	n=1 (0.32%)

Conclusions:

The data generated in this study demonstrate Copan Colibrí®:

- is able to automatically prepare microbial suspensions from different bacteria colonies directly from agar plate culture.
- provides comparable results to those obtained with manual preparation
- standardizes the preparation of microbial suspensions for AST and purity plates eliminating the person to person variability.
- Reduce the exposure to pathogens improving staff safety

Reference: ⁽¹⁾ Magiorakos *et al* 2011; Multidrug-resistant, extensively drug-resistant and pandrug-resistant bacteria: an international expert proposal for interim standard definitions for acquired resistance; European Society of Clinical Microbiology and Infectious Diseases