Instructions for Use

Copan Universal Transport Medium (UTM-RT®) System Instructions for use

INTENDED USE

Copan Universal Transport Medium (UTM-RT®) System is intended for the collection and transport of clinical specimens containing viruses, chlamydiae, mycoplasma or ureaplasma from the collection site to the testing laboratory. UTM-RT® can be processed using standard clinical laboratory operating procedures for viral, chlamydial, mycoplasma and ureaplasma culture.

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of infections caused by viruses, chlamydiae, mycoplasmas or ureaplasmas involves collection and refrigerated transport of biological specimens. Using the UTM-RT® System, the collected specimen can be stored for up to 48 hours at 2-25 °C. The UTM-RT® consists of a Hanks' Balanced Salt Solution (HBSS) enriched with proteins and sugars with a neutral pH and pH indicator. The medium contains some antibiotics and antimycotics to inhibit overgrowth of bacteria and yeasts, maintain cellular integrity and encourage preservation of viruses and chlamydiae.

After collection, the specimen should be stored at 2-25 °C and processed within 48 hours.

If delivery and processing exceed 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70 °C or colder.

PRODUCT DESCRIPTION

UTM-RT® System is ready for use and requires no further preparation. It is available in the various configurations listed in **Table 1** and supplied in a labelled screw-cap test tube filled with different volumes of UTM-RT®. The packaging in kits also includes a sterile collection device.

REF	PRODUCT DESCRIPTION		PACKAGING
	TUBE	SWAB	
302C; 302C.LC.	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	Two regular size applicator swabs polyester tipped with breaking point	50 kits per package 6 x 50 kits per box
305C; 305CMH;	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
328C	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swabs polyester tipped with breaking point	50 kits per package 6 x 50 kits per box
306C; CA302MH; 346C	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
307C; CA303MH; 357C	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
321C	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
		One regular size applicator swab with flocked nylon fiber tip with breaking point	
323C, 3U031N	2 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
338C	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
		One regular size applicator swab with flocked nylon fiber tip with breaking point	
340C	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab polyester tipped with breaking point	50 kits per package 6 x 50 kits per box
		One Combo wire - plastic shaft minitip size applicator swab polyester tipped pre-scored	
347C	3 ml of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
355C; 355CW	3 ml of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
356C	3 ml of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab polyester tipped with breaking point	50 kits per package 6 x 50 kits per box
359C; 359CE.A	1 ml of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box

360C	1 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One flexible applicator swab with flocekd nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
361C	1 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocced nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
365C	2 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One flexible applicator swab with flocekd nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
366C; 3C004N	2 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocced nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
367C.TD	1 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One contoured size with low stopper applicator swab with flocced nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
368C	1 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocced nylon fiber tip with breaking point One flexible size applicator swab with flocekd nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
3U006N; 3U008N	3 ml of UTM-RT [®] medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible applicator swab with flocekd nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
3E010N.A	1 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with round bottom.	One regular size applicator swab with flocced nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
3E012N.A	1 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with round bottom.	One flexible size applicator swab with flocekd nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box
330C;	3 ml of UTM-RT [®] medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box
331C; CA339MH	10 ml of UTM-RT [®] medium in 25x90 mm screw-cap tube with internal shaped conical bottom.	NA	15 tubes per package 6 x 15 tubes per box
348C	2,5 ml of UTM-RT [®] medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box
349C;	1,5 ml of UTM-RT [®] medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box
363C	1,5 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box
350C; 350CV	1 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box
353C	3 ml of UTM-RT [®] medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box

Not all the product codes (REF) are salable in all countries. Please contact Copan Customer care service for product codes availability for a specific country.

Table 1: product description

For all the product codes, excluding 3U008N, 331C and 348C, the tube contains three glass beads that allow homogenizing the specimen by vortexing.

REAGENTS

The UTM-RT[®] formulation includes proteins for virus stabilization¹⁷, antibiotics and antimycotics to prevent overgrowth of bacterial and fungal flora and a buffer solution to maintain a neutral pH.

Components
Sucrose
HBSS solution
Bovine serum albumin
Buffered solution
Gelatin
Amino acids
Antibiotics
Phenol Red

pH 7,3 ± 0,2 a 2+25 °C

REQUIRED MATERIALS BUT NOT PROVIDED

Materials suitable for isolation, differentiation and culture of viruses, chlamydiae, mycoplasmas and ureaplasmas.

STORAGE

The product must be stored in its original packaging at a temperature between 2 and 25 °C until the time of use. Do not overheat or freeze prior to use.

LIMITATIONS

1. Because calcium alginate swabs are toxic for many enveloped viruses⁵ and may interfere with immunofluorescence tests², they should not be used for specimen collection.
2. Wooden shaft swabs may contain toxins and formaldehydes^{2,6} and should not be used.
3. The product codes without beads (REF 3U008N, 331C and 348C) are not suitable for use with mucous or particularly viscous specimens.
4. UTM-RT[®] kits are intended to be used with the medium tubes and swabs provided in the kit. The use of tubes of medium or swabs from any other source could affect the performance of the product.

WARNINGS AND PRECAUTIONS

1. Single-use device for professional in vitro diagnostic use.
2. Do not use beyond the expiry date.
3. Do not use the UTM-RT[®] medium for premoistening or prewetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites
4. Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.⁷
5. Do not re-sterilize unused swabs.
6. Do not re-pack
7. Not suitable to collect and transport microorganisms other than viruses, chlamydiae, mycoplasma and ureaplasma
8. Do not ingest the medium
9. Specimens for the search of viruses, chlamydiae, mycoplasmas and ureaplasmas must be collected and handled using personal protective equipment against biological risk according to published manuals and guidelines.^{1,4,6,7,9,17}
10. Repeated freezing and thawing of specimens may reduce the recovery of viable organisms.
11. Do not use UTM-RT[®] if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the color of the medium has changed from light orange-red, (4) the swab pouch is open, or (5) there are other signs of deterioration.
12. The use of this product in combination with diagnostic kits or instruments must be validated by the user prior to use.
13. Due to the design of the minitip flocked swab, the swab will coil when placed in the tube. Therefore, it is not recommended to remove the swab from the tube. To process the specimen, collect the liquid using a sterile pipet or loop. If the user must remove the swab, use caution and observe adequate biohazard precaution to protect the operator and the environment in case of splash.

INSTRUCTIONS FOR USE

Proper collection of the specimen from the patient is a crucial aspect for successful isolation and identification of infectious organisms. Specimens should be collected as soon as possible after the clinical onset of disease. Highest viral titers are present during the acute illness.

UTM-RT[®] in kit


1. Open the UTM-RT[®] kit package and remove the medium test tube and the internal bag containing the sterile swab.
2. Take the sterile swab out of its bag and collect the clinical specimen; to prevent the risk of contamination, make sure that the swab tip comes into contact with the collection site only.
3. After collecting the specimen, insert the swab into the test tube until the breakpoint is level with the test tube opening.
4. Bend the swab shaft at a 180 degrees angle to break it off at the breaking point. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
5. Discard the broken handle part of the swab shaft into an approved medical waste disposal container.
6. Screw the cap back onto the test tube and hermetically seal it.

UTM-RT[®] in bulk

1. Aseptically remove cap from tube taking care not to spill the medium.
2. Aseptically place vesicle aspirates²⁶, corneal or conjunctival scrapings²⁸, small pieces of tissue or stool³⁰ samples into the tube with UTM-RT[®] medium.
3. Screw the cap back onto the test tube and hermetically seal it.
4. Identify the tube containing the specimen.
5. Send to the laboratory for immediate analysis

If processing is delayed (over 48 hours), the specimens must be frozen at -70°C or colder.

Fig 1. Collection swab showing breakpoint indication line and area for holding the applicator

DISPOSAL

Waste must be disposed of in compliance with local legislation. Take the appropriate precautions for infected material if necessary.

QUALITY CONTROL

The UTM-RT[®] lots are tested for microbial contamination, toxicity for the cellular lines used for the viral cultures and the ability to maintain the viability of viral, chlamydia and mycoplasma strains for up to 48 hours at 2-2.5°C in accordance with the methods described in CLSI M40-A2⁵.

RESULTS AND PERFORMANCE

The results obtained largely depend on proper and adequate specimen collection as well as the promptness with which the specimens are transported to the laboratory and analysed.

Viability studies were performed using Copan UTM-RT[®] with a panel of representative strains of the various families supported by the UTM-RT[®]. The swabs that accompany each transport system were directly inoculated in triplicate with 100 µl of organism suspension. Subsequently, the swabs were inserted in the respective test tubes containing the transport medium and stored for 0 and 48 hours at 2-6°C and at controlled room temperature (20-25°C). At the time of processing, each swab was vortexed for 20 seconds and removed from its transport medium test tube, then, an aliquot of the suspension was inoculated into the cellular line (200 µl) or into the appropriate culture medium^{6,15}. All the cultures were processed using the standard laboratory culture technique^{6,15}. Organism viability was determined by fluorescent cell counting for viral and chlamydia strains and CFU counting for mycoplasma and ureoplasma strains. The acceptability limits for time zero and for 48 hours were defined in accordance with the regulations M40-A2⁵.

Copan UTM-RT[®] System preserved the viability of all the organisms tested for 48 hours at both controlled room temperature and in the refrigerator in the above described test conditions. The organisms evaluated and the results obtained are given in the table below.

Organism	ATCC number	% of reduction of fluorescing infected cells after 48 hours time 2-6 °C	% of reduction of fluorescing infected cells after 48 hours time 20-25 °C
Herpes Simplex Virus Type 1	ATCC VR-539	≤50%	≤50%
Herpes Simplex Virus Type 2	ATCC VR-734	≤50%	≤50%
Respiratory Syncytial Virus	ATCC VR-1580	≤50%	≤50%
Coxsackie B1 Virus	ATCC VR-28	≤50%	≤50%
Chlamydia trachomatis	ATCC VR-880	≤50%	≤50%
Influenza A	ATCC VR-1679	≤50%	≤50%
Cytomegalovirus	ATCC VR-977	≤50%	≤50%
Varicella-zoster virus	ATCC VR-1367	≤50%	≤50%
Chlamydia pneumoniae	ATCC VR-1360	≤50%	≤50%
Adenovirus	DHI 20-4740010	45%	47%
Parainfluenza virus type 3	DHI 20-4770010	0%	41%

Echovirus type 30	ATCC VR-322	13%	24%
Organism	ATCC number	% of reduction of CFUs after 48 hours time 2-6 °C	% of reduction of CFUs after 48 hours time 20-25 °C
Mycoplasma hominis	ATCC 23114	35%	25%
Ureaplasma urealyticum	ATCC 27618	0%	6%
Mycoplasma pneumoniae	ATCC 15531	50%	50%

TABLE OF SYMBOLS

See the table of symbols at the end of the instructions for use.

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UA.TR.039

BIBLIOGRAPHY

1. Gary W. Procop and Elmer W. Koneman, 2016. Color Atlas and Textbook of Diagnostic Microbiology, Seventh edition. Wolters Kluwer Health.
2. James H. Jorgensen, Michael A. Pfaller, Karen C. Carroll, Guido Funke, Marie Louise Landry, Sandra S. Richter, David W. Warnock, 2015. Manual of Clinical Microbiology, 11th Edition. ASM, Washington, DC.
3. James Versalovic, Karen C. Carroll, Guido Funke, James H. Jorgensen, Marie Louise Landry, David W. Warnock, 2011. Cumitech 15A. Laboratory Diagnosis of Viral Infections. ASM, Washington, DC.
4. Patricia Tille. 2014. Bailey & Scott's Diagnostic Microbiology, 13th Edition. Laboratory Medicine.
5. Clinical and Laboratory Standards Institute (CLSI), 2014. M40-A2 Quality Control of Microbiological Transport Systems; Approved Standard-Second Edition.
6. Clinical and Laboratory Standards Institute (CLSI), 2006. M41-A Viral Culture; Approved guidelines.
7. Wardford, A., M. Chernesky, and E. M. Peterson, 1999. Cumitech 19A, Laboratory Diagnosis of Chlamydia trachomatis Infections. ASM, Washington DC.
8. 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
9. J. Michael Miller, Shelley A. Miller, 2017. A Guide to Specimen Management in Clinical Microbiology, Third Edition. ASM, Washington DC.
10. Centers for Disease Control and Prevention (CDC), 2016. Guide for Shipping Infectious Substances.
11. Centers for Disease Control and Prevention (CDC), 2009. Biosafety in Microbiological and Biomedical Laboratories 5th Edition.
12. World Health Organization 2015. Guidance on regulations for the Transport of Infectious Substances 2015 – 2016.
13. Centers for Disease Control and Prevention (CDC), 2002. Screening Tests to Detect Chlamydia trachomatis and Neisseria gonorrhoeae Infections - 2002.
14. Centers for Disease Control and Prevention (CDC), 2014. Recommendations for the Laboratory-Based Detection of Chlamydia trachomatis and Neisseria gonorrhoeae - 2014.
15. European Collection of Authenticated Cell Cultures (ECACC). Fundamental Techniques in Cell Culture Laboratory Handbook, 3rd edition (Protocol 9 – Detection of Mycoplasma by Culture Isolation).
16. J.B. Mahony, M.A. Chernesky Effect of Swab Type and Storage Temperature on the Isolation of Chlamydia trachomatis from Clinical Specimens Journal of Clinical Microbiology, Nov. 1985, p. 865-867.
17. **S. Specter, R. L. Hodinka, S. A. Young. Clinical Virology Manual, fifth edition, 2016.**
18. Peter Daley, Santina Castriciano, Max Chernesky, Marek Smieja. Comparison of Flocked and Rayon Swabs for Collection of Respiratory Epithelial Cells from Uninfected Volunteers and Symptomatic Patients. Journal of Clinical Microbiology, June 2006.
19. K. Loens, L. Van Heirstraeten, S. Malhotra-Kumar, H. Goossens, and M. Ieven. Optimal Sampling Sites and Methods for Detection of Pathogens Possibly Causing Community-Acquired Lower Respiratory Tract Infections. Journal of Clinical Microbiology, Jan. 2009.
20. Richard Garceau MD, Danielle Leblanc RT, Louise Thibault MD, Gabriel Girouard MD, Manon Mallet PhD. Herpes simplex virus type 1 is the leading cause of genital herpes in New Brunswick. Can J Infect Dis Med Microbiol 2012.
21. Marek Smieja, Santina Castriciano, Susan Carruthers, Geoffrey So, Sylvia Chong, Kathy Luinstra, James B. Mahony, Astrid Petrich, Max Chernesky, Mario Savarese, and Daniele Triva. Development and Evaluation of a Flocked Nasal Midturbinate Swab for Self-Collection in Respiratory Virus Infection Diagnostic Testing. Journal of Clinical Microbiology, Sept. 2010.
22. Donghyok Kwon, Kyeoncheol Shin, Mihwa Kwon, Hee-Bok Oh, Chun Kang, Joo-Yeon Lee. Development and Evaluation of Rapid Influenza Diagnostic Test for the Pandemic (H1N1) 2009 Virus. JCM Accepts, published online ahead of print on 27 October 2010.
23. Bhupesh K. Prusty, Christine Siegl, Petra Hauck, Johannes Hain, Suvi J. Korhonen, Eija Hiltunen- Back, Mirja Puolakkainen, Thomas Rudel. Chlamydia trachomatis Infection Induces Replication of Latent HHV-6. Plos One, April 2013.
24. Ellen Vancutsem, Oriane Soetens, Maria Breugelmanns, Walter Foulon, Anne Naessens. Modified Real-Time PCR for Detecting, Differentiating, and Quantifying *Ureaplasma urealyticum* and *Ureaplasma parvum*. The Journal of Molecular Diagnostics, Vol. 13, No. 2, March 2011.
25. Turkiya Al-Siyabi, Khalifa Binkhamis, Melanie Wilcox, Sallene Wong, Kanti Pabbaraju, Raymond Tellier, Todd F. Hatchette, and Jason J. Le Blanc. A cost-effective real-time PCR for the detection of adenovirus from viral swabs. Al-Siyabi et al. Virology Journal 2013.
26. Eidgenössisches Departement für Verteidigung, Bevölkerungsschutz und Sport VBS; SAMPLING AND DISPATCH OF SAMPLES - VIROLOGY Ident-Nr./Vers. 10014156143/01 Aktenzeichen:323-FachinfoViro12016.
27. Matthew L. Faron et al (2016); Multicenter Evaluation of Meridian Bioscience HSV 1&2 Molecular Assay for Detection of Herpes Simplex Virus 1 and 2 from Clinical Cutaneous and Mucocutaneous Specimens; Journal of Clinical Microbiology Volume 54 Number 8
28. Al-Siyabi et al. (2013); A cost effective real-time PCR for the detection of adenovirus from viral swabs. Virology Journal 10:184.
29. Le Roy C, Pereyre S, Hénin N, Bébéar C. 2017. French prospective clinical evaluation of the Aptima Mycoplasma genitalium CE-IVD assay and macrolide resistance detection using three distinct assays. J Clin Microbiol 55:3194–3200. <https://doi.org/10.1128/JCM.00579-17>.
30. Minodier et al. (2017) Clinical and virological factors associated with gastrointestinal symptoms in patients with acute respiratory infection: a two-year prospective study in general practice medicine BMC Infectious Diseases 17:729 DOI 10.1186/s12879-017-2823-9

Index of Symbols / Tabella dei Simboli / Tabla de símbolos / Symboltabelle / Table des Symboles / Tabela de símbolos / Таблица символов

Symbol / Simbolo / Símbolo / Symbol / Símbolo / Symbole / Símbolos / Символ	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Значення
	Manufacturer / Fabbricante / Fabricante / Hersteller / Fabricant / Fabricante / Виробник
	In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro/ Dispositivo de diagnóstico in vitro / Виріб для діагностики in vitro
	CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE/ Marcação CE / Маркування CE
	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Ідентифікаційний номер уповноваженого органу
	Sterilized using ethylene oxide / Sterilizzato usando ossido di etilene / Esterilizado por óxido de etileno / Sterilisiert mit Äthylenoxid / Stérilisé à l'aide d'oxyde d'éthylène / Esterilizado utilizando óxido de etileno / Стерилізований завдяки етиленоксиду
	Sterilized using irradiation / Sterilizzato usando radiazioni ionizzanti / Esterilizado usando radiaciones ionizantes / Sterilisiert mit ionisierenden Strahlungen / Stérilisé à l'aide de radiations ionisantes / Esterilizado por radiación ionizante / Стерилізований завдяки опроміненню
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Не використовувати повторно
	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Número do catálogo / Номер за каталогом
	Temperature limitation / Limiti di temperatura / Límites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Температурні обмеження
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Використати до
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultar as instruções de utilização / Користуватись інструкціями по використанню
	Peel / Strappare per aprire / Desprender / Abziehen / Décoller / Destacar / Відкривати тут
	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Код серії (партії)
	Contains sufficient for <n> tests / Contenido suficiente per <n> test / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contém o suficiente para <n> testes / Придатний для проведення <n> кількості випробувань
	Do not use if package is damaged / Non utilizzare in caso di confezionamento danneggiato / No utilizar en caso de paquetedañado / Bei Beschädigung der Verpackung nicht verwenden / Ne pas utiliser si l'emballage est abîmé / Não utilizar se a embalagem estiver danificada / Не використовувати, якщо упаковка пошкоджена
	This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Isto aplica-se apenas aos EUA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Це стосується винятково США: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."

Copan Italia S.p.A.
Via F. Perotti, 10
25125 Brescia Italy
Tel +39 030 2687211
Fax +39 030 2687250

Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:
Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562 USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com

Copyright © 2020 Copan Italia S.p.A.
All rights reserved