ITALIANO

USO PREVISTO

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ENGLISH**INTENDED USE**

TODD HEWITT CNA BROTH (LIM) is a ready to use medium for selective enrichment for group-B Streptococci (*Streptococcus agalactiae*).

SUMMARY AND EXPLANATION

Streptococcus agalactiae (GBS) is one of the main causes of serious infections, such as sepsis, meningitis, and pneumonia in newborn babies, as well as an etiological agent in diabetic adults or chronic patients. It is now universally accepted that in most cases the transmission from a carrier mother to the fetus occurs in the uterus through ascending transmission right before or during delivery (1). The identification of maternal colonization during the last weeks of pregnancy is extremely important, as the administration of preventive treatment upon delivery, thus the fetus protection from a possible infection, will depend on such identification. The guidelines developed by the US Centers for Diseases Control provide indications for optimal treatment of samples to be analyzed for suspected GBS. Specifically, for the detection of anal-genital carrier status in pregnant women, samples should be taken on the 35th to 37th pregnancy week, with a vaginal swab and a second anal-rectal swab. The standard method for GBS diagnosis is based on culture in specific cultural media. The use of a selective broth inhibiting the growth of the other microorganisms in the swab increases the method selectivity by 60-90% (2).

PROCEDURE PRINCIPLES

Todd Hewitt forms the broth base and is a growth medium mainly used for culture of Streptococci, especially for serological studies (3). Dextrose stimulates haemolysin formation; at the same time, dextrose fermentation is the source of energy. The buffering system (sodium bicarbonate and sodium phosphate) contributes to maintaining the pH level, by countering the acidity formed by sugar fermentation. Through the addition of an antimicrobial supplement (Nalidixic Acid and Colistin) with inhibiting action against Gram-negative microorganisms, the cultural medium becomes selective.

REAGENTS**TODD HEWITT BROTH CNA (LIM) Composition**

Infusion of meat (fat-free)

Tryptone

Glucose

Sodium bicarbonate

Sodium chloride

Di-sodium phosphate

Colistin sulphate

Nalidixic acid

pH

7.7 ± 0.2 AT 25 °C

USE PRECAUTIONS AND WARNING NOTICES

1. TODD HEWITT CNA BROTH is a SINGLE USE device for diagnostic use in vitro. Its reuse may cause infections and/or inaccurate results.
2. Observe approved precautions for potentially hazardous biological materials, as well as techniques in asepsis.
3. Use must be limited to properly trained and qualified staff. The manufacturer shall not be liable for any use by non-qualified or non-authorized people.
4. Clinical samples may contain pathogens, including hepatitis or human immunodeficiency viruses. Handle all materials and items contaminated with blood and other biological fluids in accordance with the guidelines and "standard precautions" provided by institutional bodies.
5. All samples and materials used for product analysis shall be considered as potentially infected and shall be managed in a way as to avoid the risk of infection for laboratory staff. After use, sterilize all potentially hazardous waste, including samples, containers, and transport and culture media.
6. Observe the instructions rigorously.
7. Not appropriate for transport of microorganisms not included in the use destination.
8. Do not ingest the medium.
9. Do not use the medium to pre-moisten or pre-wet the swab before sample collection, or to rinse or irrigate the sampling sites.
10. The product cannot be reused; its reuse may cause inappropriate diagnosis.
11. Do not use in case of evident damage to the product.
12. Avoid any skin contact with the cultural medium.

In the case sample taking and transport have been performed using common taking and transport systems, and the sample should be processed with automatic plating systems, remove the swab from the tube containing the enriched medium before processing the sample through automatic systems.

STORAGE AND SHELF LIFE

TODD HEWITT CNA BROTH (LIM) is ready to use and does not require any other preparations. It shall be kept in its original package at 5-25 °C until it is used. Do not overheat. Do not incubate or freeze before use. Inappropriate storage causes a loss of effectiveness. Do not use after the expiry date shown on the external package and medium tube label.

PRODUCT DETERIORATION

Do not use TODD HEWITT CNA BROTH (LIM) if the product is clearly damaged or contaminated, and leakages or other deterioration signs are visible.

PROCEDURE

MATERIAL SUPPLIED with the number 476CE.A:

TODD HEWITT CNA BROTH (LIM)

MATERIALS REQUIRED BUT NOT SUPPLIED: Ancillary culture plates, reagents, organisms for quality assurance, various laboratory instruments, and Regular Flocked Swab with fracture point printed at 80 mm (53380CS01).

NOTE: In the case the swab is processed through automatic systems, the stem validated for use is the Regular Flocked Swab with fracture point printed at 80 mm (53380CS01).

SAMPLE MANAGEMENT: For screening of bacterial colonization by Streptococcus B in pregnant women, vaginal/rectal swabs are normally used. After samples have been taken, swabs kept in appropriate transport and storage medium shall be preferably kept in refrigerated conditions until tests are performed.

USE INSTRUCTIONS:

- Remove the swab(s) from the transport medium and introduce it in the selective medium tube, squeezing the tube walls to extract as much sample material as possible. If necessary, vortex for a few minutes to facilitate complete sample release from the swab. *In the case TODD HEWITT CNA BROTH (LIM) is used in combination with Regular Flocked Swab (53380CS01) as the maintenance system, homogenize the sample by shaking or vortexing the tube.*
- Incubate aerobically or in CO₂-enriched atmosphere for 18-24 hours at 35-37 °C.
- At the end of incubation, if turbidity is observed in the tube, spread on appropriate Blood Agar plate. Extend incubation in tube for 24 hours more before confirming any negatives. Spreading on plate can be performed through:
 - Plating with WASP automation; refer to the instruments manual for any additional explanations. Remember to remove the swab from the medium, if it is in the tube and if it is different from the Regular Flocked Swab (53380CS01).
 - Direct plating with swab.
 - Plating of 100 µl suspension taken with micropipette and sample distribution on plate with a spatula.
- Incubate the plates for 18-24 hours, or up to 48 hours, if necessary, at 35-37 °C in aerobic conditions.
- Observe the growth of typical grey-white beta-hemolytic or non-hemolytic colonies (GBS is a gram-positive catalysis-negative coccus). The incubation should be extended to 48 hours, if no group-B Streptococcus colonies are identified at the end of the 18-24 hour period.
- Use the isolated colonies for further confirmation through agglutination tests or recommended tests for detection of group-B Streptococcus antigens, referring to the manufacturer procedures for use instructions and reading of results.

RESULTS

The results obtained substantially depend on proper sampling, as well as on quick transport and treatment in laboratories.

QUALITY ASSURANCE AND PERFORMANCE MONITORING

Quality assurance procedures include system performance monitoring. TODD HEWITT CNA BROTH (LIM) is a selective growth medium for GBS. Therefore, the internal quality assurance procedures include the review of both characteristics. In particular, the medium maintenance and growth capacity is checked using a pure suspension of the target microorganism (GBS), whereas a pure suspension of a gram-negative microorganism (*Escherichia coli*) is used to check the inhibition capacity.

Procedure:

- Starting from a fresh culture on plate, prepare a 0.5 McFarland suspension of the microorganism in physiological solution.
- 0.1 ml dilutions 1:100 of the two microorganisms taken with graduated micropipette are used to incubate the medium tubes in aseptic conditions.
- The tubes are incubated in thermostat at the temperature of 35 °C for 24 hours; the inspection of tubes and subculture on Blood Agar Plate are used as acceptability parameters for batch release.

Microorganism	ATCC	Result
<i>Streptococcus agalactiae</i>	ATCC 12386	Growth
<i>Escherichia coli</i>	ATCC 25922	Inhibition (from partial to total)

AVAILABILITY

CODE	DESCRIPTION
476CE.A	6 boxes containing 50 tubes with 2 ml TODD HEWITT medium

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