Copan SELENITE- Product Insert & How to Use Guide

ENGLISH

PI24A REV.01 2012.10

INTENDED USE

Copan Selenite Broth is an Enrichment selective broth for clinical specimens. It can be used for research of *Salmonella* spp in stool samples or others biological materials . **Copan Selenite Broth** can be used for processing with automated systems for planting and streaking on culture media plates (Copan WASP included).

REAGENTS

Copan Selenite Broth formulation per liter :

Peptone 5g
Lactose 4g
Phosphate sodium salt 10 g
Sodium selenite salt 4 g

PRECAUTIONS

1. Observe approved biohazard precautions and aseptic techniques.
2. To be used only by adequately trained and qualified personnel.
3. All specimens and materials used to process them should be considered potentially infectious and handled in a manner which prevents infection of laboratory personnel. Sterilize all biohazard waste including specimens, containers and media after their use. Observe other CDC Biosafety Level 2 recommendations.
4. Directions should be read and followed carefully.

STORAGE

This product is ready to use and no further preparation is necessary. The product has to be stored in its original container at 2-25°C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box, on each individual kit and on the specimen transport tube label.

Copan Selenite Broth is supplied in boxes of 50 units individually packed and sterilize.

PRODUCT DETERIORATION

Do not use **Copan Selenite Broth** if: (1) the product shows visible marks of damage or contamination; (2) there is evidence of leakage; (3) the expiration date has passed; (4) there are other signs of deterioration.

MATERIALS SUPPLIED

Fifty (50) units are contained in a shelf pack and 6 x 50 units are contained in a box. Each unit consists of a sterile package containing a pre-labeled polypropylene screw-cap tube with conical shaped bottom filled with 2 ml of Selenite Broth .

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for isolation and for bacteria culture. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques.

Catalog No.	Copan Selenite Broth Product Descriptions	Pack Size	Suitable for Automation
475CE.A	2 ml of Selenite broth in 12X80 mm polypropylene screw-cap tube with internal conical shape	50 units per shelf pack 6x50 units per box	YES
475CE	2 ml of Selenite broth in 12X80 mm polypropylene screw-cap tube with internal conical shape	50 units per shelf pack 6x50 units per box	NO

INSTRUCTIONS FOR USE

Copan Selenite broth is used for the selective enrichment of *Salmonella* spp present in stool specimens or fecal samples collected by rectal swabbing. Copan Selenite is equipped with a special capture cap designed for use with Copan's regular size flock swabs. A rectal swab collected using a regular size flocked swabs, such as 552C, or transported using Copan FecalSwab Kit (470CE.A, 470CE) can be transferred into Copan Selenite Broth. The capture cap facilitates easy handling during planting and streaking in manual mode. In the case of 552C simply break the swab at the hemodend breakpoint inside the Selenite tube and close the cap; the swab will then remain attached to the cap. Alternatively a 502C or 552C flocked swab can be used for transferring approximately 0.2g of stool sample into Copan Selenite broth.

1. Unscrew the Selenite broth tube cap.
 - 2.. Add the specimen:
 - rectal swab: insert the swab in the tube and break it at the breaking point
 - feces sample: add about 0,2 g of stool by using a regular size flocked swab or a traditional paddle
- In order to have optimal medium performance the ratio between the medium and stool must be 10:1.
3. Replace the cap on the tube and screw the cap tightly.
 4. Write the patient or the specimen information on the tube label or apply patient/specimen identification label. Send immediately the specimen to the laboratory

SPECIMEN PROCESSING IN LABORATORY

PROCEDURE, MANUAL OPERATION

Samples collected in Selenite

- Vortex the Selenite tubes for 10sec.

Samples collected in FecalSwab

- Unscrew the Selenite broth tube cap.
 - Vortex the specimen FecalSwab tube for 10 sec.
 - Unscrew the FecalSwab tube cap and transfer the swab in the broth by using sterile forceps or alternatively transfer about 30µL of Cary-Blair medium by using a micropipette.
 - Re-cap the FecalSwab tube with the original cap
 - Re-cap the inoculated Selenite broth with the original cap.
 - Vortex the inoculated Selenite tubes for 10sec.
1. Incubate at 35°C ±2°C for 18-24 hours .
 2. After incubation , streak the sample on a surface of specific agar plate for isolation of *Salmonella* sp. The subculture of the specimen can be performed directly by using a swab or planting an aliquot of sample by an adequate instrument for liquid handling (micropipette, loop etc.) and following streaking on the plate. Plant at least 1µL of sample.
 3. Incubate the agar plates as per standard laboratory procedure.

SPECIMEN PROCESSING IN LABORATORY

PROCEDURE, AUTOMATIC OPERATION

The Selenite broth is manufactured in a container suitable to be processed by automatic system for planting and streaking (WASP). By using the WASP is possible to inoculate specimens collected by the FecalSwab system into the broth and to subculture the enriched Selenite broth sample.

BROTH INOCULATION BY THE WASP

1. Process the specimen collected in FecalSwab by a specific broth inoculation protocol by using the 30µL loop. Refer to the WASP User Manual for further information.
2. Examine the BHI broth tube for growth after 18-24 hours and, if negative, re-incubate for additional 24 hours.

SUBCULTURE OF THE ENRICHED SELENITE SAMPLE BY THE WASP

1. Process the enriched Selenite sample by a specific protocol for the streaking of an appropriate selective medium plate. The use of a 1 or 10µL loop and a multiple quadrant streaking pattern (3-4 Quadrants type) is suggested in order to improve the colonies isolation. It is recommended the laboratory to validate the streaking protocol. Refer to the WASP User Manual for further information.
2. Incubate the plates according to the standard laboratory procedure.

QUALITY CONTROL

The Copan Selenite Broth has been tested with the follow procedure :

PERFORMANCE TEST :

- From fresh culture of *Salmonella enterica* Serovar Typhimurium ATCC 14028, prepare a 0.5 Mac Farland in PBS.
- Prepare a serial dilution 10⁻⁵ of 0.5 Mac Farland
- Spike 100 µl of inoculum directly into Copan Selenite Broth
- Recap the tubes
- Vortex for 10 second at 2500 RPM
- Plate at time zero , 100 µl per dilution on a XLD agar plate.
- Incubate the tube at 35°C ±2°C for 24 hours .
- After incubation , plate 100 µl of Selenite broth , on a XLD agar plate for isolation of *Salmonella* .
- Incubate the agar plate as per standard laboratory procedure 35°C ±2°C for 24 hours.

Tine zero count should be 30-300 CFU/plate and the 24 hours time plate must show any grow.

INHIBITION TEST :

- From fresh culture of *E.coli* ATCC 25922, make a 0.5 Mac Farland in PBS.
- Prepare a serial dilution 10⁻¹ of 0.5 Mac Farland
- Spike 100 µl of inoculum directly into Copan Selenite Broth
- Recap the tubes
- Vortex for 10 second at 2500 RPM
- Plate at time zero , 100 µl per dilution on McConkey agar.
- Incubate at 35°C ±2°C for 24 hours .
- After incubation , plate 100 µl of Selenite broth , on a McConkey agar plate for isolation of *E.coli*
- Incubate the agar as per standard laboratory procedure

Tine zero count should show confluent grow and the 24 hours time plate must show partial to total inhibition.

Copan Italia S.p.A. Via Perotti, 10 - 25125 BRESCIA - Italy Tel.: +39 030 2687211 - Fax: +39 030 2687206 E-mail: info@copanitalia.com - Web site: www.copanitalia.com

Strains	Results cfu/plate Time Zero Average on 5 plates	Results cfu/plate Time 24 Hours Average on 5 plates
<i>Salmonella enterica</i> Serovar Typhimurium ATCC 14028	172	CONFLUENT GROW
<i>E.coli</i> ATCC 25922	CONFLUENT GROW	INIBITION
<i>Shigella sonnei</i> ATCC 9290	202	CONFLUENT GROW

LIMITATIONS

- Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.
- Performance testing with Copan BHI broth was conducted using laboratory strains spiked into the BHI broth tube and not using human specimens.
- Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals. Specimens should be collected as soon as possible after the clinical onset of disease. Highest bacterial titres are present during the acute illness.
- Selenite broth is an enrichment selective broth for isolation of *Salmonella* spp. The product is not intended to maintain the viability of viruses and Chlamydia
- Stool specimens require special attention. They should be collected early in the course of the disease and need to be cultured within 2 hours of collection. Optimally, the specimen should be transferred to transport media soon after collection (such as Cary-Blair or Enteric Pathogen Transport for specimens, or Stuarts for swabs) to maintain the viability of the organisms. Specimens should never be refrigerated.
- In Selenite Broth without dulcitol, coliform growth is suppressed for the first 8-12 hours; longer incubation can lead to overgrowth by the mixed flora.

WARNINGS

- For *in vitro* diagnostic use.
-  This product is for single use only; reuse may cause a risk of inaccurate results.
- Not suitable for any other application than intended use.
- Work under a biological safety cabinet, wear gloves and other protection commensurate with universal precautions when handling clinical specimens.
- The use of this product in association with any diagnostic assay or with any diagnostic instrumentation should be validated by the user before using.
- Do not use if the product is visibly damaged.
- The medium is TOXIC.Do not ingest the medium.
- The manufacturer cannot be held responsible for any inappropriate or unqualified use of the product.
- All specimens considered infectious and must be handled with appropriate precautions.
- Copan Selenite broth is for in-vitro diagnostics use only and is in no way intended for curative or prophylactic purposes.
- Observe approved biohazard precautions and aseptic techniques. Product to be used only by adequately trained and qualified personnel.
- Directions should be read and followed carefully.
- Caution: This medium is toxic if swallowed, inhaled or comes into contact with the skin. Wear gloves and eye/face protection.

BIBLIOGRAPHY

- Murray, P.R., E.J. Baron, M.A. Pfäller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Atlas, R.N.1993 Handbook of Microbiological Media p. 147-153. CRC Press, Boca Raton, FL

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BIBLIOGRAFIA

- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Atlas, R.N.1993 Handbook of Microbiological Media p. 147-153. CRC Press, Boca Raton, FL

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BIBLIOGRAPHIE

- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isernberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Atlas, R.N.1993 Handbook of Microbiological Media p. 147-153. CRC Press, Boca Raton, FL

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