THIOL Broth - Product Package Leaflet and Use Instructions

INTENDED USE

THIOL broth is an enrichment medium which facilitates the growth of facultative and obligate aerobic and anaerobic microorganisms. COPAN Thiol Broth is used for the enrichment of anaerobic bacteria isolated from cultures of clinical samples. THIOL Broth can be used with automatic plate culture systems (Copan WASP™ or Copan WASPLab™).

SUMMARY

THIOL broth is an enrichment medium consisting of dextrose, peptone and L-cystine solution and yeast extract, which provide the essential growth factors for bacterial reproduction.

The essential ions are provided by sodium chloride. Sodium thioglycollate and L-cystine act as reducing agents, which contribute to maintaining a low redox potential, supporting anaerobiosis.

Thiol broth also contains growth supplements such as Hemin and Vitamin K1 which facilitate the reproduction of some particularly demanding anaerobic bacterial strains.

REAGENTS

Components of THIOL Broth:

Names of components
L-cystine
Sodium chloride
Glucose
Yeast extract
Pancreatic digest of casein
Sodium thioglycollate
Hemin
Vitamin K1
Distilled water

STORAGE

Thiol Broth is a ready-made product and no additional preparation operations are required. The batch packed in its original package can be kept at 5-25°C until it is used or until the expiry date indicated on the tube and on the box. Do not overheat. Do not incubate or freeze before use. Inappropriate storage may cause a loss of product effectiveness.

PRODUCT DETERIORATION

Do not use THIOL Broth if: (1) the product shows visible signs of damage or contamination; (2) there are indications of leakages from the tube; (3) it is past its expiry date; (4) there are other visible deterioration signs.

SUPPLIED MATERIALS

Code	Product Description	Package Dimensions	Suitable for automation
4U002N	12x80 mm PET test tubes with PE perforated screw cap. One test tube contains 4 ml of THIOL Broth.	50 units per internal package. 50x6 units per box.	YES

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for culture and isolation of bacteria; refer to laboratory reference manuals for the protocols recommended for the culture and identification techniques.

INSTRUCTIONS FOR USE

For specific guidelines on sampling procedures, consult the reference manuals published. The Copan eSwab system can be used as a system for sampling, transport and preservation of the clinical sample. Refer to Copan for information on the most appropriate code, and refer to the product package leaflet for use instructions.

USER INSTRUCTIONS

MANUAL OPERATIONS

1. Unscrew the THIOL broth tube cap.
2. Inoculate the sample into the open tube.
For liquid samples, use a 10 µl loop or micropipette to transfer the sample to the tube.
For samples in eSwab:
a. Homogenize the tube for 10 sec with vortex at 2000/2500 RPM, unscrew the eSwab tube cap, immerse the swab in the Thiol Broth medium and shake vigorously for 5/10 seconds to release the sample into the enrichment medium, dispose of the swab and re-close the Thiol Broth tube.
b. Using a micropipette, transfer 30 µl of liquid Amies medium directly into the THIOL Broth. Fit the caps on the eSwab and Thiol Broth tubes.
3. Homogenize for 5-10 seconds at 2000/2500 rpm by a vortex.
4. Incubate the inoculated THIOL Broth tubes at 35°C ± 2°C.

5. Examine the THIOL Broth tube after 24-48 hours. The tube must be checked daily. Absence of turbidity in the tube does not indicate absence of the microbe growth.
6. Proceed with surface culture on appropriate medium in a plate. Surface culture on plate can be performed through surface culture of a sample portion taken by a micropipette or a loop, followed by distribution on plate. The minimum recommended quantity for surface culture in a plate after incubation of the Thiol Broth tube is 1-10 µl. Incubate the plates for 24-96 hours at 35°C ± 2°C in appropriate anaerobiosis conditions, or according to standard laboratory procedure.
7. At the end of incubation, proceed with plate reading.
8. You are advised to plate the inoculated Thiol Broth sample daily.
9. You are advised to check growth in the plate daily.

OPERATIONS WITH AUTOMATIC SYSTEM (WASP™)

Thiol Broth can be used on automatic culture systems like COPAN WASP™. For correct use, see the WASP™/WASPLab™ instruction manual.

LIMITATIONS

1. Some difficult samples with specific requirements may not grow in this medium.
2. The conditions and the volume of the sample collected for culture are significant variables for obtaining reliable culture results.
3. An appropriate sample collection from patients is an extremely critical factor for isolation and identification of the bacterial strains of interest.
4. Absence of turbidity in the Thiol Broth tube is not indicative of a negative result.
5. Check the tube daily until the end of incubation.
6. You are advised to plate both after incubation of 24 hours and after incubation of 48 hours in an appropriate growth medium.
7. Maintain the Thiol Broth tube in incubation for a maximum of 96 hours at 35°C ± 2°C.

NOTICES

1. For diagnostic use in vitro.
2. ⚠ This product is designed for single use; reuse may lead to unreliable results.
3. Not appropriate for any application other than its intended use.
4. Follow the recommended guidelines for collection of the samples.
5. Work in a laminar flow cabinet after inoculation of the product, wear gloves and other protections in line with universally recognized precautions for handling of clinical samples.
6. The use of this product with any diagnostic test or with any diagnostic tool shall be assessed and if necessary validated by the user before use.
7. Do not use the product if it is visibly damaged.
8. Do not ingest.
9. Do not inhale.
10. Do not inject.
11. The manufacturer shall not be considered liable for any inappropriate or non-qualified use of the product.
12. All clinical samples are considered potentially infective and shall be handled with appropriate precautions.
13. The THIOL Broth shall only be used for in vitro diagnostics. Its use for therapeutic or prevention purposes is not allowed under any circumstances.
14. Take appropriate precautions against biological hazards and apply aseptic techniques. The product may only be used by properly trained and qualified operators.
15. The instructions must be read and observed carefully.
16. Do not keep the tube in anaerobiosis conditions during incubation.
17. Do not place in an autoclave.
18. Do not filter.

WASTE DISPOSAL

Unused reagents can be considered as non-hazardous waste and disposed of accordingly.

For disposal, refer to the product MSDS.

Used reagents and any other contaminated waste material shall be disposed of according to the procedures established for infected or potentially infected products. The laboratory shall be responsible for managing waste materials by treatment and disposal, as established in all applicable regulations.

QUALITY CONTROL PROCEDURE

PERFORMANCE TEST

- > From a fresh culture, prepare a 0.5 McFarland suspension in PBS of every ATCC organism.
- > Prepare a serial 10 -4 dilution of 0.5 McFarland suspension.
- > Inoculate with 400 µl of the dilution in PBS directly in THIOL broth.
- > Put the caps back on the tubes.
- > Homogenize for 10 seconds at 2000/2500 RPM with vortex.
- > For the control at zero time, seed 100 µl of the bacterial suspension on nutrient medium (e.g. blood agar plate) and incubate in appropriate growth conditions.
- > Incubate THIOL broth tube at 35°C ± 2°C for 24-48 hours.
- > After incubation, seed 100 µl THIOL broth on appropriate nutrient medium (e.g. blood agar plate).
- > Incubate the plates at 35°C ± 2°C for 24-48 hours in appropriate growth conditions.

LIMITS OF ACCEPTABILITY: Time zero: CFU counting 30-300 CFU/100 µl

Time 24-48 h: growth

RESULTS

STRAIN	ZERO DEAD TIME COUNTING: CFU/PLATE	CFU COUNTING AFTER 18-24 HOURS AT 35°C±2°C; CFU/PLATE
Peptostreptococcus anaerobius ATCC 27337*	216	GROWTH
Bacteroides Vulgatus ATCC 8482*	198	GROWTH
Clostridium Perfringens ATCC 13124*	134	GROWTH
Prevotella Melaninogenica ATCC 25845	121	GROWTH
Propionibacterium acnes ATCC 6919	166	GROWTH
Bacteroides fragilis ATCC 25285	200	GROWTH
Staphylococcus aureus ATCC 25923	210	GROWTH

The results obtained and presented are representative of the values obtained during general evaluation of the product.
The results obtained and presented are representative of tests performed at the COPAN laboratory with ATCC strains.
The results obtained and presented were NOT obtained using clinical samples.

* Bacterial performance tests conducted using strains established by the CLSI-M22-A3 reference standard, specific for quality control of enrichment mediums for the growth of anaerobic bacteria.

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BIBLIOGRAPHIE

- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology. 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Atlas, Handbook of Microbiological Media 4th edition. CRC Press, Boca Raton FL

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Literaturverzeichnis:

- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Atlas, Handbook of Microbiological Media 4th edition 2010 CRC Press, Boca Raton, FL

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BIBLIOGRAFIA

- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Atlas, Handbook of Microbiological Media 4th edition 2010 CRC Press, Boca Raton, FL.

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BIBLIOGRAFIA

- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology, 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (NCCLS). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Isenberg, H. D., 2004. Clinical Microbiology Procedures Handbook, 2nd ed. ASM, Washington, DC.
- Atlas, Handbook of Microbiological Media 4th edition 2010 CRC Press, Boca Raton, FL.

Index of Symbols/Table des Symboles/Symbole/Tabla de Símbolos/Tabella dei Simboli

Symbol/Symbole/Simbolo	Meaning/Signification/Bedeutung/Significado/Significato
	Manufacturer/Fabricant/Hersteller/Fabricante/Fabbricante
	Diagnosi in vitro
	Do no reuse/Ne pas réutiliser/Nicht zur Wiederverwendung/No reutilizar/Non riutilizzare
	Catalogue number/Référence du catalogue/Bestellnummer/Número de catálogo/Numero di catalogo
	Temperature limitation/Limites de temperature/Temperaturbegrenzung/Límites de temperatura/Limiti di temperatura
	Use by/Utiliser jusque/Verwendbar bis/Fecha de caducidad/Utilizzare entro
	Consult Instructions for Use/Consulter les instructions d'utilisation/Gebrauchsanweisung beachten/Consulte las instrucciones de uso/Consultare le istruzioni per l'uso
	Batch code (Lot)/Code de lot (Lot)/Chargencode (Chagenbezeichnung)/Código de lote (Lote)/Codice del lotto (partita)
	Contains sufficient for <n> tests/Contenu suffisant pour <n> tests/Ausreichend für <n> Tests/Contenido suficiente para <n> pruebas/Contenuto sufficiente per <n> test

Copan Italia Spa
Via Perotti, 10
25125 Brescia, Italia

Tel: +39 030 2687211
Fax: +39 030 2687250
Email: info@copanitalia.com
www.copanitalia.com