MSwab™

**Collection, Preservation and Transport System  
Product Insert & How to Use Guide**

**Copan MSwab™ Collection, Preservation and Transport System - Product Insert & How To Use Guide**

See symbol glossary at end of insert

**INTENDED USE**

MSwab™ is a Collection, Transport and Preservation System intended for the collection and transport of clinical specimens containing Gram positive aerobic and facultative anaerobic bacteria, HSV 1 and HSV 2 from the collection site to the testing laboratory. In the laboratory, MSwab™ specimens are processed using standard clinical laboratory operating procedures for culture.

**SUMMARY AND PRINCIPLES**

One of the routine procedures in the diagnosis of bacteriological infections involves the collection and safe transportation of swab samples. This can be accomplished using the Copan MSwab™ that is a Collection, Transport and Preservation System. Copan MSwab™ incorporates a transport and preservation medium containing TRIS HCl, EDTA, TRIS Base, Dimethyl Sulfoxide (DMSO) and Bovine Serum Albumin. The medium is designed to maintain the viability of Gram positive aerobic and facultative anaerobic bacteria and HSV 1 and HSV 2 during transit to the testing laboratory.

Copan MSwab™ Collection, Transport and Preservation System is supplied in two different formats: a) Collection kit format. Each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 1 ml or 1.6 ml of MSwab™ transport and preservation medium and a small sterile peel pouch containing one specimen collection swab that has a tip flocced with soft nylon fiber. b) Tube only format. A plastic screw-cap tube with conical shaped bottom filled with 1 ml or 1.6 ml of MSwab™ transport and preservation medium.

Once a swab sample is collected, it should be placed immediately into the MSwab™ transport tube where it comes into contact with the transport medium. Swab specimens for bacterial or viral investigations collected using MSwab™ should be transported directly to the laboratory, preferably within 2 hours of collection (1, 2, 7) to maintain optimum organism viability. If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 hours. Bacterial viability studies for *Staphylococcus aureus*, ATCC 29213 and ATCC 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC 43300 and ATCC 700698 show viability of the tested organisms up to 14 days at refrigerated temperature (4 – 8°C) or 72 hours at room temperature (20 – 25°C). Independent scientific studies on swab transport systems have shown that for certain bacteria viability is superior at refrigerated temperatures compared with room temperature (12 – 21). If viral samples must be frozen, they should be frozen at -70°C.

**REAGENTS****Formulation of MSwab™ Transport Medium**

TRIS HCl  
EDTA  
TRIS Base  
Dimethyl Sulfoxide (DMSO)  
Bovine Serum Albumin  
Distilled water

**PRECAUTIONS**

1. This product is For In Vitro Diagnostic Use.
2. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
3. All specimens and materials used to process them should be considered potentially infectious and handled in a manner which prevents infection of laboratory personnel. Sterilize all biohazard waste including specimens, containers and media after their use. Observe other CDC Biosafety Level 2 recommendations (31, 32, 33, 34).
4. Directions should be read and followed carefully.

**STORAGE**

This product is ready for use and no further preparation is necessary. This product should be stored in its original container at 5 – 25°C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual collection unit and the specimen transport tube label.

**PRODUCT DETERIORATION**

Copan MSwab™ should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed, (4) the swab package is open, or (5) there are other signs of deterioration.

**SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION**

Specimens collected for microbiological investigations which comprise the isolation of bacteria or viruses should be collected and handled following published manuals and guidelines (7, 8, 4).

To maintain optimum organism viability, transport specimens collected using MSwab™ directly to the laboratory, preferably within 2 hours of collection (1, 2, 7). If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 hours. Bacterial viability studies for *Staphylococcus aureus*, ATCC 29213 and ATCC 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC 43300 and ATCC 700698 show viability of the tested organisms up to 14 days at refrigerated temperature (4 – 8°C) or 72 hours at room temperature (20 – 25°C). If viral samples must be frozen, they should be frozen at -70°C.

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations (34, 35, 36, 37). Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

**MATERIALS SUPPLIED**

Fifty (50) MSwab™ collection units are contained in a shelf pack and 6 x 50 units are contained in a box. Each collection unit consists of a package containing two components: a pre-labeled polypropylene screw-cap tube with conical shaped bottom filled with 1 ml or 1.6 ml of MSwab™ transport medium and a specimen collection swab which has a tip flocced with soft nylon fiber (see Fig 1). Two types of collection formats are available; both include a tube of medium, one includes also sterile swab applicator single wrapped in peel pouch.

Copan MSwab™ Collection, Transport and Preservation System is supplied in two different formats: a) Collection kit format. Each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 1 ml or 1.6 ml of MSwab™ transport and preservation medium and a small sterile peel pouch containing one specimen collection swab that has a tip flocked with soft nylon fiber. b) Tube only format. A plastic screw-cap tube with conical shaped bottom filled with 1 ml or 1.6 ml of MSwab™ transport and preservation medium.

#### MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for isolating and culturing aerobic and facultative anaerobic bacteria. These materials include culture media plates or tubes and incubation systems. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques for aerobic and facultative anaerobic from clinical swab samples (2, 4).

Appropriate materials for isolating, differentiating and culturing viruses. These materials include tissue culture cell lines, tissue culture medium, incubation systems and reading equipment. Refer to appropriate references for recommended protocols for isolation and identification of viruses (1, 7).

#### DIRECTIONS FOR USE

Copan MSwab™ Collection, Transport and Preservation System is available in product configurations indicated in the table below.

Catalog No.	Copan MSwab™ Product Descriptions	Pack Size	Capture Cap Feature*
403C	Single use sample collection pack containing: - Polypropylene screw-cap tube with internal conical shape filled with 1ml of MSwab™ Medium. - One regular size applicator swab with flocked nylon fiber tip, sterile and individually wrapped.	50 units per shelf pack 6x50 units per box	YES
404C 404C.R	Single use sample collection pack containing: - Polypropylene screw-cap tube with internal conical shape filled with 1.6ml of MSwab™ Medium. - One regular size applicator swab with flocked nylon fiber tip, sterile and individually wrapped.	50 units per shelf pack 6x50 units per box	YES
406C	Single transport and preservation tube: - Polypropylene screw-cap tube with internal conical shape filled with 1ml of MSwab™ Medium.	50 units per shelf pack 6x50 units per box	YES
407C	Single transport and preservation tube: - Polypropylene screw-cap tube with internal conical shape filled with 1.6 ml of MSwab™ Medium.	50 units per shelf pack 6x50 units per box	YES

Other product codes may be available. For updates please refer to our website: [www.copanflock.com](http://www.copanflock.com)

\*Capture Cap Feature is guaranteed only when using Copan Regular Size Flock Swab.

#### Specimen Collection

Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals (7, 2).



For MSwab™ codes 404C, 404C.R and 403C:

- Peel open the kit package and remove the tube of medium and inner pouch containing the sterile swab applicator (see Figure 1).
- Remove the swab applicator from its peel pouch and use to collect the clinical specimen. The operator must touch the swab applicator only above the colored breakpoint line, as illustrated in Figure 1, which is the opposite end to the nylon fiber tip. At all times when handling the swab applicator, the operator must not touch the area below the colored breakpoint line (the area from the line to the tip of the nylon flocked swab) as this will lead to contamination of the applicator shaft and the subsequent culture. During sample collection when handling the swab applicator, the operator must not touch the area below the pink breakpoint indication line; that is the area from the line to the tip of the nylon flocked swab (see Figure 1), as this will lead to contamination of the applicator shaft and the culture thus invalidating the test results.
- After the swab sample is taken from the patient, break off the swab applicator shaft at the colored breakpoint line into the MSwab™ tube containing the MSwab™ transport medium.
- Screw the cap and replace and secure tightly onto the tube (see Figure 1). Write patient's name and demographics on the tube label and send the sample to the laboratory.

For MSwab™ codes 406C and 407C:

- After the swab sample is taken from the patient, break off the swab applicator shaft at the colored breakpoint line, if present, into the MSwab™ tube containing the MSwab™ transport medium.
- Screw the cap and replace and secure tightly onto the tube (see Figure 1). Write patient's name and demographics on the tube label and send the sample to the laboratory.

**Fig 1. Collection swab showing breakpoint indication line and area for holding the applicator**

The operator must only handle the part of the swab applicator shaft above the breakpoint indication line as shown in Fig 1. After the swab sample is taken from the patient, the swab applicator shaft is broken off at the colored breakpoint indication line into the MSwab™ tube of transport medium. The operator then discards the handle part of the swab into an approved medical waste disposal container. The tube's screw cap is then replaced and secured tightly. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a funnel shaped molded docking receptacle in the cap (see Fig 2). This molded funnel shape captures the end of the broken applicator shaft and secures it firmly in the dock by friction grip.

**Fig 2. Capture of broken swab applicator stick by MSwab™ tube cap**

In the testing laboratory when the MSwab™ cap is unscrewed and removed, the swab applicator stick is securely attached to the cap. This feature allows the operator to conveniently remove the swab and perform various microbiology analyses using the tube cap as a handle to hold and manipulate the swab.

**Processing MSwab™ Specimens in the Laboratory – Bacteriology**

MSwab™ samples should be processed for bacteriological culture using recommended culture media and laboratory techniques which will depend on the specimen type and the organism under investigation. For recommended culture media and techniques for the isolation and identification of bacteria from clinical swab specimens refer to published microbiology manuals and guidelines (1-6).

Culture investigations of swab specimens for the presence of bacteria routinely involve the use of solid agar culture medium in Petri dish plates. The procedure for inoculation of MSwab™ samples onto solid agar in Petri dishes is as follows.


Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (31, 32, 33, 34).

Vortex mix the MSwab™ tube containing the swab sample for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the medium

1. Unscrew the MSwab™ cap and remove the swab applicator.
2. Roll the tip of the MSwab™ applicator onto the surface of one quadrant of the culture media plate to provide the primary inoculum.
3. If it is necessary to culture the swab specimen onto a second culture media plate, return the MSwab™ applicator to the transport medium tube for two seconds to absorb and recharge the applicator tip with transport medium/patient sample suspension then repeat Step No. 3.
4. If it is necessary to inoculate additional culture media plates, return the MSwab™ applicator to the transport medium tube and recharge the swab applicator tip with the transport medium/patient sample suspension before inoculating each additional plate.


The procedure described above utilizes the MSwab™ applicator like an inoculation wand to transfer the suspension of patient sample in transport medium onto the surface of a culture plate creating the primary inoculum (see Fig 3). Alternatively, the operator can vortex mix the MSwab™ tube with the swab inside for 5 seconds and then transfer 100µl volumes of the suspension onto each culture plate using a volumetric pipeter and sterile pipet tips. Standard laboratory techniques should then be used to streak the primary inoculum of patient sample across the surface of the culture plate (see Fig 4).

**Fig 3. Procedures for inoculation of MSwab™ specimens onto solid agar in Petri dishes**

1. Using swab to inoculate specimen
2. Using pipetor and sterile pipet tips to inoculate 100µl of specimen

**Fig 4. Procedure for streaking MSwab™ specimens on agar Petri dishes for primary isolation (33)**

Seed a primary inoculum of MSwab™ specimen onto the surface of an appropriate agar culture plate in the first quadrant.

Use a sterile bacteriology loop to streak the primary inoculum across the surface of the second, third and fourth quadrants of the agar culture plate.

#### **Preparation of Gram Stain Smears of MSwab™ Specimens**

Laboratory analysis of clinical swab samples collected from certain sites on the patient can routinely include microscopic examination of stained preparations (“direct smears”) using the Gram stain procedure. This can provide valuable information to physicians who are managing patients with infectious diseases (22). There are many instances in which a Gram stain can assist in making a diagnosis (23, 27).

The Gram stain can also help to judge specimen quality and contribute to the selection of culture media especially with mixed flora.

Microscope slides of patient specimens transported in Copan MSwab™ transport system can be prepared for Gram stain analysis, as describe below, by sampling an aliquot of vortexed suspension of the swab (3, 4). Sample transported in MSwab™ elution medium represent a homogeneous suspension in liquid phase. It can be uniformly smeared allowing clear and easy reading.

Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (31, 32, 33, 34).

1. Take a clean glass microscope slide, place it on a flat surface and inscribe an area using a diamond-tipped or similar glass marker to identify the location of the specimen inoculum. Note: a slide with a pre-marked 20 mm well can be used.
2. Vortex mix the MSwab™ tube containing the swab sample for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the medium.
3. Unscrew the MSwab™ cap and using a sterile pipet, transfer 1 – 2 drops of sample suspension to the inscribed area on the glass slide. Note: about 30µl would be a suitable amount of liquid for a pre-marked 20 mm diameter well slide.

In case of bloody or thicker specimens particular care should be taken to thinly spread the sample on the slide. Bacteria are difficult to detect if the sample shows many red cells and debris.

4. Allow the specimen on the slide to air dry at room temperature or place the slide in an electric slide warmer or incubator set at a temperature not exceeding 42°C.
5. Fix smears using methanol. Methanol fixation is recommended as it prevents lysis of Red Blood Cells, avoids damage to all host cells and results in a cleaner background (3, 4, 22).
6. Follow published laboratory reference manuals and guidelines for performing the Gram stain. If commercial Gram stain reagents are used, it is important to comply with instructions in the manufacturer’s product insert for performance test procedure.

For further information or guidance on the preparation of specimen slides for microscopic analysis, for information on Gram staining procedures and the interpretation and reporting of microscopic analysis, consult published laboratory reference manuals (1 - 5, 22 - 27).

#### **Processing MSwab™ Specimens in the Laboratory – Virology**

The survival of HSV 1 and HSV 2 depends on many factors including the type and concentration of the microorganism, duration of transport and storage temperature. To maintain optimum viability, specimens should be transported directly to the laboratory preferably within 2 hours of collection (1, 2, 7, 29). If immediate delivery or processing is delayed, then specimens collected using Copan MSwab™ Collection, Transport and Preservation System should be refrigerated at 4-8°C or stored at room temperature (20-25°C) and processed within 48 hours. If samples must be frozen, they should be frozen placed at -70°C.

In simulated transport and storage studies, the Copan MSwab™ System was shown to be capable of maintaining the viability of HSV 1 and HSV 2 at refrigerated (4-8°C) and room temperature (20-25°C) conditions for up to 48 hours. Based on performance studies conducted by Copan and independent scientific publications, viability of certain microorganisms is superior at refrigerated temperatures compared with room temperature (12 – 21, 29).

MSwab™ specimens should be processed for virology culture using recommended cells line and laboratory techniques that will depend on the specimen type and the organism under investigation. For recommended shell vials and techniques for the isolation and identification of HSV 1 and HSV 2 from clinical swab specimens, refer to published virology manuals and guidelines (1 – 6, 29, 30).

Culture investigations of swab specimens for the presence of HSV 1 and HSV 2 routinely involve the use of cells culture in shell vials. The procedure for inoculation of MSwab™ specimens onto shell vials is described below.

1. Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other BSL 2 recommendations.
2. Shake using a vortex mixer for 5 seconds the MSwab™ tube containing the swab specimen to release the specimen from the swab tip and evenly disperse and suspend the patient specimen in the liquid transport medium.
3. Unscrew the MSwab™ cap and remove the swab applicator.
4. Transfer 200 µl volumes of the suspension into the shell vial and proceed as per laboratory internal procedure.
- NOTE: Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.
5. Proceed with the appropriate techniques for virus detection.

#### QUALITY CONTROL

MSwab™ applicators are tested to ensure they are non-toxic to bacteria. MSwab™ medium and applicators are tested to ensure they are non-toxic to cell lines used for HSV 1 and HSV 2 cultivation. MSwab™ transport medium is tested for pH stability (9). MSwab™ is quality control tested before release for ability to maintain viable Gram positive aerobic and facultative anaerobic bacteria and HSV viruses at room temperature (20 – 25°C) for specified time points. Procedures for quality control of microbiology transport devices should be conducted using testing methods described in Clinical and Laboratory Standards Institute M40-A and other publications (9). If aberrant quality control results are noted, patient results should not be reported.

#### LIMITATIONS

1. In the laboratory, wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (31, 32, 33, 34) when handling or analyzing patient samples.
2. Condition, timing and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection (7, 8, 4).
3. MSwab™ is intended for use as a collection and transport medium for Gram positive aerobic and facultative anaerobic bacteria, HSV 1 and HSV 2 viruses. MSwab™ cannot be used as enrichment, selective or differential medium.
4. MSwab™ is an antibiotics free medium. Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.
5. Performance testing with Copan MSwab™ was conducted using laboratory strains spiked onto a swab following the test protocols described in Clinical and Laboratory Standards Institute M40-A Approved Standard (9). Performance testing was not conducted using human specimens.
6. Performance testing with Copan MSwab™ was conducted using Copan flocked swabs.

#### WARNINGS

1. Do not re-sterilize unused swabs.
2.  This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
3. Do not re-pack.
4. Not suitable to collect and transport microorganisms other than Gram positive aerobic and facultative anaerobic bacteria, HSV 1 and HSV 2 viruses.
5. Not suitable to collect and transport fastidious or anaerobic bacteria
6. Not suitable for any other application than intended use.
7. The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
8. Do not use if the swab is visibly damaged (i.e., if the swab tip or swab shaft is broken).
9. Do not use excessive force or pressure when collecting swab samples from patients as this may result in breakage of the swab shaft.
10. Applicator swab is qualified as Class IIa Medical Device according to European Medical Device Directive 93/42/EEC - Surgically Invasive Transient Use.
11. Class IIa means swabs can be used for sampling body surfaces, body orifices (e.g., for example nose, throat, vagina, wounds, groin or skin)).
12. Do not ingest the medium.
13. Directions for use must be followed carefully. The manufacturer cannot be held responsible for any unauthorized or unqualified use of the product.
14. To be handled by trained personnel only.
15. It must be assumed that all specimens contain infectious micro-organisms; therefore all specimens must be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste. Observe CDC Biosafety Level 2 recommendations (31, 32, 33, 34).
16. Do not use the MSwab™ medium for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.

#### RESULTS

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

#### PERFORMANCE CHARACTERISTICS

The test procedures employed for determining bacterial viability performance were based upon the quality control methods described in Clinical and Laboratory Standards Institute M40-A (9).

MSwab™ system has an intended use limited to Gram positive aerobic and facultative anaerobic bacteria and HSV 1 and HSV 2, therefore its field application is more restricted than some other devices. For this reason the bacterial recovery studies were conducted under the simulated transport and storage conditions as described and defined in CLSI M40-A, Quality Control of Microbiological Transport Systems: Approved Standard and included only the Gram positive aerobic and facultative anaerobic strains from the Group 1 of paragraph 7.11.1 of the CLSI M40-A document, in particular:

Streptococcus pyogenes	ATCC® 19615
Streptococcus pneumoniae	ATCC® 6305

*Pseudomonas aeruginosa*

ATCC® BAA-427

In addition, Copan included testing of additional Gram positive aerobic and facultative anaerobic microorganisms, clinically relevant, not required by CLSI M40-A. The specific bacterial strains used in these studies are here listed:

Enterococcus faecalis	ATCC® 29212
Staphylococcus epidermidis	ATCC® 12228
Staphylococcus aureus	ATCC® 29213
Streptococcus agalactiae (Group B Strep)	ATCC® 13813
Kocuria rhizophila	ATCC® 9341
Listeria monocytogenes	ATCC® 19114
Bacillus cereus	ATCC® 10876
Staphylococcus aureus (Methicillin resistant)	ATCC® 43300
Staphylococcus aureus	ATCC® 6538
Staphylococcus aureus (Methicillin resistant)	ATCC® 700698
Streptococcus pneumoniae	ATCC® 49136

All cultures of bacteria were ATCC (American Type Culture Collection) and were obtained commercially.

The selection of these organisms also reflects those Gram positive aerobic and facultative anaerobic bacteria that would normally be encountered in specimens collected and analyzed in a typical clinical microbiology laboratory.

Bacterial viability studies were performed on the Copan MSwab™ at two different ranges of temperature, 4 – 8°C and 20 – 25°C, corresponding to refrigerator and room temperature, respectively. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs, 24 hrs and 48 hrs. At the appropriate time intervals, each swab was processed according to the Roll-Plate or Swab Elution Method.

Additional bacterial viability studies for *Staphylococcus aureus*, ATCC 29213 and ATCC 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC 43300 and ATCC 700698 were performed on the Copan MSwab™ at two different ranges of temperature, 4 – 8°C and 20 – 25°C, corresponding to refrigerator and room temperature, respectively.

Swab accompanying the transport system was inoculated in triplicate with 100µl of specific concentrations of organism suspension.

Swabs were then placed in their respective transport medium tubes and:

For studies performed at 4 – 8°C inoculated MSwab™ tubes were held for 0 hrs, 10 days and 14 days. At the appropriate time intervals, each MSwab™ was processed according to the Roll-Plate Method.

For studies performed at 20 – 25°C inoculated MSwab™ tubes were held for 0 hrs and 72 hrs. At the appropriate time intervals, each MSwab™ was processed according to the Roll-Plate Method.

Bacterial overgrowth studies were performed on the Copan MSwab™ at 4 – 8°C, corresponding to refrigerator temperature. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs and 48 hrs. At the appropriate time intervals, each swab was processed according to the Roll-Plate Method.

Bacterial overgrowth studies were performed using *Pseudomonas aeruginosa*.

Viral viability studies were performed using HSV 1 and HSV 2. Swabs accompanying each transport system were directly inoculated in triplicate with 100µl of organisms suspension. Swabs were then placed in their respective transport medium tubes and were held for 0, 24 and 48 hours at both 4°C and room temperature (20-25°C). At the appropriate time interval, each swab was vortexed, removed from its transport medium tube and then 200µl aliquots of this suspension was inoculated into shell vials. All cultures were processed by standard laboratory culture technique and examined after a specified incubation time. Organism viability was determined by fluorescing foci counts.

Organisms evaluated were:

Herpes Simplex Virus Type 1 (HSV 1)	ATCC VR-539
Herpes Simplex Virus Type 2 (HSV 2)	ATCC VR-734.

**TEST RESULTS**

**\*Table 1 - SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES SWAB ELUTION METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 24 hrs	Average of CFUs recovered at time 48 hrs	Log <sub>10</sub> decline	Interpretation
<i>Streptococcus pyogenes</i> ATCC 19615	diluted 1:10	MSwab™	2045	3.9 x 10 <sup>5</sup>	2.7 x 10 <sup>5</sup>	2.3 x 10 <sup>5</sup>	-0.23	Acceptable Recovery
			2045/1	4.8 x 10 <sup>5</sup>	2.5 x 10 <sup>5</sup>	2.1 x 10 <sup>5</sup>	-0.36	Acceptable Recovery
			2045/2	4.2 x 10 <sup>5</sup>	2.3 x 10 <sup>5</sup>	2.0 x 10 <sup>5</sup>	-0.32	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC 6305	diluted 1:10	MSwab™	2045	1.2 x 10 <sup>6</sup>	1.8 x 10 <sup>4</sup>	2.1 x 10 <sup>3</sup>	-1.76	Acceptable Recovery
			2045/1	1.2 x 10 <sup>6</sup>	1.9 x 10 <sup>4</sup>	1.7 x 10 <sup>3</sup>	-1.85	Acceptable Recovery
			2045/2	1.2 x 10 <sup>5</sup>	2.0 x 10 <sup>4</sup>	1.9 x 10 <sup>3</sup>	-1.80	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC 49136	diluted 1:10	MSwab™	2045	2.0 x 10 <sup>5</sup>	2.0 x 10 <sup>3</sup>	2.1 x 10 <sup>2</sup>	-1.98	Acceptable Recovery
			2045/1	1.8 x 10 <sup>5</sup>	2.1 x 10 <sup>3</sup>	1.8 x 10 <sup>2</sup>	-2.00	Acceptable Recovery
			2045/2	1.9 x 10 <sup>5</sup>	1.8 x 10 <sup>3</sup>	1.7 x 10 <sup>2</sup>	-2.05	Acceptable Recovery
<i>Enterococcus faecalis</i> ATCC 29212	diluted 1:10	MSwab™	2045	1.1 x 10 <sup>6</sup>	1.0 x 10 <sup>6</sup>	9.4 x 10 <sup>5</sup>	-0.07	Acceptable Recovery
			2045/1	1.1 x 10 <sup>6</sup>	1.0 x 10 <sup>6</sup>	9.2 x 10 <sup>5</sup>	-0.08	Acceptable Recovery
			2045/2	1.1 x 10 <sup>6</sup>	1.0 x 10 <sup>6</sup>	9.5 x 10 <sup>5</sup>	-0.06	Acceptable Recovery
<i>Staphylococcus epidermidis</i> ATCC 12228	diluted 1:10	MSwab™	2045	1.6 x 10 <sup>6</sup>	1.2 x 10 <sup>6</sup>	9.2 x 10 <sup>5</sup>	-0.24	Acceptable Recovery
			2045/1	1.7 x 10 <sup>6</sup>	1.3 x 10 <sup>6</sup>	1.1 x 10 <sup>6</sup>	-0.19	Acceptable Recovery
			2045/2	1.4 x 10 <sup>6</sup>	1.2 x 10 <sup>6</sup>	9.5 x 10 <sup>5</sup>	-0.17	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC 29213	diluted 1:10	MSwab™	2045	1.3 x 10 <sup>6</sup>	1.0 x 10 <sup>6</sup>	8.8 x 10 <sup>5</sup>	-0.17	Acceptable Recovery
			2045/1	1.0 x 10 <sup>6</sup>	9.5 x 10 <sup>5</sup>	7.4 x 10 <sup>5</sup>	-0.13	Acceptable Recovery
			2045/2	1.0 x 10 <sup>6</sup>	9.5 x 10 <sup>5</sup>	8.0 x 10 <sup>5</sup>	-0.10	Acceptable Recovery
<i>Streptococcus agalactiae</i> ATCC 13813	diluted 1:10	MSwab™	2045	2.1 x 10 <sup>6</sup>	1.8 x 10 <sup>6</sup>	1.5 x 10 <sup>6</sup>	-0.15	Acceptable Recovery
			2045/1	1.9 x 10 <sup>6</sup>	1.8 x 10 <sup>6</sup>	1.5 x 10 <sup>6</sup>	-0.10	Acceptable Recovery
			2045/2	2.2 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	1.3 x 10 <sup>6</sup>	-0.23	Acceptable Recovery
<i>Kocuria rhizophila</i> ATCC 9341	diluted 1:10	MSwab™	2045	2.0 x 10 <sup>6</sup>	1.5 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	-0.15	Acceptable Recovery
			2045/1	2.1 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	-0.18	Acceptable Recovery
			2045/2	2.2 x 10 <sup>6</sup>	1.6 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	-0.19	Acceptable Recovery
<i>Listeria monocytogenes</i> ATCC 19114	diluted 1:10	MSwab™	2045	1.4 x 10 <sup>6</sup>	8.1 x 10 <sup>5</sup>	3.0 x 10 <sup>5</sup>	-0.67	Acceptable Recovery
			2045/1	1.3 x 10 <sup>6</sup>	9.0 x 10 <sup>5</sup>	3.4 x 10 <sup>5</sup>	-0.58	Acceptable Recovery
			2045/2	1.5 x 10 <sup>6</sup>	8.8 x 10 <sup>5</sup>	4.2 x 10 <sup>5</sup>	-0.55	Acceptable Recovery
<i>Bacillus cereus</i> ATCC 10876	diluted 1:10	MSwab™	2045	1.5 x 10 <sup>6</sup>	8.9 x 10 <sup>5</sup>	1.6 x 10 <sup>5</sup>	-0.97	Acceptable Recovery
			2045/1	1.6 x 10 <sup>6</sup>	7.6 x 10 <sup>5</sup>	1.8 x 10 <sup>5</sup>	-0.95	Acceptable Recovery
			2045/2	1.6 x 10 <sup>6</sup>	8.3 x 10 <sup>5</sup>	1.4 x 10 <sup>5</sup>	-1.06	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	diluted 1:10	MSwab™	2045	2.1 x 10 <sup>6</sup>	1.2 x 10 <sup>6</sup>	8.0 x 10 <sup>5</sup>	-0.42	Acceptable Recovery
			2045/1	1.9 x 10 <sup>6</sup>	1.3 x 10 <sup>6</sup>	8.9 x 10 <sup>5</sup>	-0.33	Acceptable Recovery
			2045/2	1.9 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	9.4 x 10 <sup>5</sup>	-0.30	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC 6538	diluted 1:10	MSwab™	2045	1.5 x 10 <sup>6</sup>	9.2 x 10 <sup>5</sup>	6.9 x 10 <sup>5</sup>	-0.34	Acceptable Recovery
			2045/1	1.6 x 10 <sup>6</sup>	1.0 x 10 <sup>6</sup>	8.3 x 10 <sup>5</sup>	-0.29	Acceptable Recovery
			2045/2	1.7 x 10 <sup>6</sup>	1.3 x 10 <sup>6</sup>	9.5 x 10 <sup>5</sup>	-0.25	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 700698	diluted 1:10	MSwab™	2045	2.0 x 10 <sup>6</sup>	1.6 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	-0.15	Acceptable Recovery
			2045/1	2.1 x 10 <sup>6</sup>	1.6 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	-0.18	Acceptable Recovery
			2045/2	1.9 x 10 <sup>6</sup>	1.5 x 10 <sup>6</sup>	9.5 x 10 <sup>5</sup>	-0.30	Acceptable Recovery

**\*Table 2 - SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES SWAB ELUTION METHOD, 20-25°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 24 hrs	Average of CFUs recovered at time 48 hrs	Log <sub>10</sub> decline	Interpretation
<i>Streptococcus pyogenes</i> ATCC 19615	diluted 1:10	MSwab™	2045	3.9 x 10 <sup>5</sup>	2.2 x 10 <sup>5</sup>	1.5 x 10 <sup>5</sup>	-0.41	Acceptable Recovery
			2045/1	4.8 x 10 <sup>5</sup>	2.1 x 10 <sup>5</sup>	1.5 x 10 <sup>5</sup>	-0.51	Acceptable Recovery
			2045/2	4.2 x 10 <sup>5</sup>	2.0 x 10 <sup>5</sup>	1.7 x 10 <sup>5</sup>	-0.39	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC 6305	diluted 1:10	MSwab™	2045	1.2 x 10 <sup>6</sup>	1.0 x 10 <sup>4</sup>	9.4 x 10 <sup>2</sup>	-2.11	Acceptable Recovery
			2045/1	1.2 x 10 <sup>6</sup>	1.2 x 10 <sup>4</sup>	1.1 x 10 <sup>3</sup>	-2.04	Acceptable Recovery
			2045/2	1.2 x 10 <sup>6</sup>	1.2 x 10 <sup>4</sup>	1.1 x 10 <sup>3</sup>	-2.04	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC 49136	diluted 1:10	MSwab™	2045	2.0 x 10 <sup>6</sup>	1.5 x 10 <sup>4</sup>	1.6 x 10 <sup>3</sup>	-2.10	Acceptable Recovery
			2045/1	1.8 x 10 <sup>6</sup>	1.5 x 10 <sup>4</sup>	1.1 x 10 <sup>3</sup>	-2.21	Acceptable Recovery
			2045/2	1.9 x 10 <sup>6</sup>	1.7 x 10 <sup>4</sup>	1.4 x 10 <sup>3</sup>	-2.13	Acceptable Recovery
<i>Enterococcus faecalis</i> ATCC 29212	diluted 1:10	MSwab™	2045	1.1 x 10 <sup>6</sup>	1.1 x 10 <sup>6</sup>	7.7 x 10 <sup>5</sup>	-0.15	Acceptable Recovery
			2045/1	1.1 x 10 <sup>6</sup>	9.6 x 10 <sup>5</sup>	8.1 x 10 <sup>5</sup>	-0.13	Acceptable Recovery
			2045/2	1.1 x 10 <sup>6</sup>	9.3 x 10 <sup>5</sup>	8.6 x 10 <sup>5</sup>	-0.11	Acceptable Recovery
<i>Staphylococcus epidermidis</i> ATCC 12228	diluted 1:10	MSwab™	2045	1.6 x 10 <sup>6</sup>	1.1 x 10 <sup>6</sup>	6.2 x 10 <sup>5</sup>	-0.49	Acceptable Recovery
			2045/1	1.7 x 10 <sup>6</sup>	1.1 x 10 <sup>6</sup>	6.8 x 10 <sup>5</sup>	-0.40	Acceptable Recovery
			2045/2	1.4 x 10 <sup>6</sup>	9.0 x 10 <sup>5</sup>	5.5 x 10 <sup>5</sup>	-0.41	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC 29213	diluted 1:10	MSwab™	2045	1.3 x 10 <sup>6</sup>	8.7 x 10 <sup>5</sup>	5.2 x 10 <sup>5</sup>	-0.40	Acceptable Recovery
			2045/1	1.0 x 10 <sup>6</sup>	8.7 x 10 <sup>5</sup>	5.8 x 10 <sup>5</sup>	-0.24	Acceptable Recovery
			2045/2	1.0 x 10 <sup>6</sup>	9.9 x 10 <sup>5</sup>	6.0 x 10 <sup>5</sup>	-0.22	Acceptable Recovery
<i>Streptococcus agalactiae</i> ATCC 13813	diluted 1:10	MSwab™	2045	2.1 x 10 <sup>6</sup>	1.5 x 10 <sup>6</sup>	1.3 x 10 <sup>6</sup>	-0.21	Acceptable Recovery
			2045/1	1.9 x 10 <sup>6</sup>	1.7 x 10 <sup>6</sup>	1.3 x 10 <sup>6</sup>	-0.16	Acceptable Recovery
			2045/2	2.2 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	1.0 x 10 <sup>6</sup>	-0.34	Acceptable Recovery
<i>Kocuria rhizophila</i> ATCC 9341	diluted 1:10	MSwab™	2045	2.0 x 10 <sup>6</sup>	1.1 x 10 <sup>6</sup>	6.0 x 10 <sup>5</sup>	-0.52	Acceptable Recovery
			2045/1	2.1 x 10 <sup>6</sup>	1.5 x 10 <sup>6</sup>	5.4 x 10 <sup>5</sup>	-0.59	Acceptable Recovery
			2045/2	2.2 x 10 <sup>6</sup>	1.3 x 10 <sup>6</sup>	6.5 x 10 <sup>5</sup>	-0.53	Acceptable Recovery
<i>Listeria monocytogenes</i>	diluted	MSwab™	2045	1.4 x 10 <sup>6</sup>	3.8 x 10 <sup>5</sup>	2.0 x 10 <sup>5</sup>	-0.84	Acceptable Recovery

ATCC 19114	1:10		2045/1	1.3 x 10 <sup>6</sup>	5.3 x 10 <sup>5</sup>	1.8 x 10 <sup>7</sup>	-0.86	Acceptable Recovery
			2045/2	1.5 x 10 <sup>6</sup>	5.3 x 10 <sup>5</sup>	1.7 x 10 <sup>7</sup>	-0.94	Acceptable Recovery
Bacillus cereus ATCC 10876	diluted 1:10	MSwab™	2045	1.5 x 10 <sup>6</sup>	5.8 x 10 <sup>4</sup>	2.0 x 10 <sup>7</sup>	-1.87	Acceptable Recovery
			2045/1	1.6 x 10 <sup>6</sup>	7.8 x 10 <sup>4</sup>	2.1 x 10 <sup>7</sup>	-1.88	Acceptable Recovery
			2045/2	1.6 x 10 <sup>6</sup>	8.1 x 10 <sup>4</sup>	3.8 x 10 <sup>7</sup>	-1.62	Acceptable Recovery
Staphylococcus aureus(MRSA) ATCC 43300	diluted 1:10	MSwab™	2045	2.1 x 10 <sup>6</sup>	1.0 x 10 <sup>6</sup>	5.0 x 10 <sup>7</sup>	-0.62	Acceptable Recovery
			2045/1	1.9 x 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	4.2 x 10 <sup>7</sup>	-0.65	Acceptable Recovery
			2045/2	1.9 x 10 <sup>6</sup>	9.1 x 10 <sup>5</sup>	4.3 x 10 <sup>7</sup>	-0.65	Acceptable Recovery
Staphylococcus aureus ATCC 6538	diluted 1:10	MSwab™	2045	1.5 x 10 <sup>6</sup>	7.8 x 10 <sup>5</sup>	5.1 x 10 <sup>7</sup>	-0.47	Acceptable Recovery
			2045/1	1.6 x 10 <sup>6</sup>	1.1 x 10 <sup>6</sup>	6.4 x 10 <sup>7</sup>	-0.40	Acceptable Recovery
			2045/2	1.7 x 10 <sup>6</sup>	1.0 x 10 <sup>6</sup>	7.9 x 10 <sup>7</sup>	-0.33	Acceptable Recovery
Staphylococcus aureus (MRSA) ATCC 700698	diluted 1:10	MSwab™	2045	2.0 x 10 <sup>6</sup>	1.6 x 10 <sup>6</sup>	1.1 x 10 <sup>7</sup>	-0.26	Acceptable Recovery
			2045/1	2.1 X 10 <sup>6</sup>	1.5 x 10 <sup>6</sup>	1.3 x 10 <sup>7</sup>	-0.21	Acceptable Recovery
			2045/2	1.9 X 10 <sup>6</sup>	1.4 x 10 <sup>6</sup>	8.2 x 10 <sup>7</sup>	-0.36	Acceptable Recovery

**\*Table 3 - SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES ROLL-PLATE METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 24 hrs	Average of CFUs recovered at time 48 hrs	Interpretation
Streptococcus pyogenes ATCC 19615	diluted 10 <sup>-3</sup>	MSwab™	2045	283.3	172.3	153.0	Acceptable Recovery
			2045/1	266.7	220.0	169.0	Acceptable Recovery
			2045/2	275.7	240.7	160.3	Acceptable Recovery
Streptococcus pneumoniae ATCC 6305	diluted 10 <sup>-1.5</sup>	MSwab™	2045	227.3	153.0	111.7	Acceptable Recovery
			2045/1	215.7	171.0	109.0	Acceptable Recovery
			2045/2	217.3	159.0	104.3	Acceptable Recovery
Streptococcus pneumoniae ATCC 49136	diluted 10 <sup>-1.5</sup>	MSwab™	2045	279.7	206.3	123.7	Acceptable Recovery
			2045/1	256.7	207.7	113.0	Acceptable Recovery
			2045/2	270.0	195.7	123.0	Acceptable Recovery
Enterococcus faecalis ATCC 29212	diluted 10 <sup>-1.5</sup>	MSwab™	2045	207.3	172.3	158.3	Acceptable Recovery
			2045/1	214.0	167.7	158.3	Acceptable Recovery
			2045/2	223.3	183.3	140.0	Acceptable Recovery
Staphylococcus epidermidis ATCC 12228	diluted 10 <sup>-1.5</sup>	MSwab™	2045	280.7	238.7	123.0	Acceptable Recovery
			2045/1	277.0	238.0	113.0	Acceptable Recovery
			2045/2	273.3	151.3	117.0	Acceptable Recovery
Staphylococcus aureus ATCC 29213	diluted 10 <sup>-1.5</sup>	MSwab™	2045	252.7	209.3	179.0	Acceptable Recovery
			2045/1	212.7	188.3	160.0	Acceptable Recovery
			2045/2	217.0	188.7	157.3	Acceptable Recovery
Streptococcus agalactiae ATCC 13813	diluted 10 <sup>-3</sup>	MSwab™	2045	174.0	121.0	117.7	Acceptable Recovery
			2045/1	162.3	127.0	100.3	Acceptable Recovery
			2045/2	186.0	148.3	124.7	Acceptable Recovery
Kocuria rhizophila ATCC 9341	diluted 10 <sup>-2</sup>	MSwab™	2045	271.7	237.0	212.0	Acceptable Recovery
			2045/1	267.7	223.3	202.0	Acceptable Recovery
			2045/2	269.3	235.0	179.0	Acceptable Recovery
Listeria monocytogenes ATCC 19114	diluted 10 <sup>-3</sup>	MSwab™	2045	267.3	243.3	209.7	Acceptable Recovery
			2045/1	275.0	242.0	192.7	Acceptable Recovery
			2045/2	274.3	226.7	196.3	Acceptable Recovery
Bacillus cereus ATCC 10876	diluted 10 <sup>-1.5</sup>	MSwab™	2045	234.7	185.3	167.0	Acceptable Recovery
			2045/1	246.0	182.7	151.7	Acceptable Recovery
			2045/2	231.0	179.0	144.7	Acceptable Recovery
Staphylococcus aureus(MRSA) ATCC 43300	diluted 10 <sup>-1.5</sup>	MSwab™	2045	204.7	185.0	172.7	Acceptable Recovery
			2045/1	211.7	185.7	177.3	Acceptable Recovery
			2045/2	196.7	183.7	161.0	Acceptable Recovery
Staphylococcus aureus ATCC 6538	diluted 10 <sup>-1.5</sup>	MSwab™	2045	149.0	123.0	98.7	Acceptable Recovery
			2045/1	172.0	138.0	108.3	Acceptable Recovery
			2045/2	167.3	142.3	98.3	Acceptable Recovery
Staphylococcus aureus (MRSA) ATCC 700698	diluted 10 <sup>-1.5</sup>	MSwab™	2045	257.0	223.3	204.0	Acceptable Recovery
			2045/1	254.3	229.0	203.3	Acceptable Recovery
			2045/2	257.7	238.0	202.7	Acceptable Recovery

**\*Table 4 - SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES ROLL-PLATE METHOD, 20-25°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 24 hrs	Average of CFUs recovered at time 48 hrs	Interpretation
Streptococcus pyogenes ATCC 19615	diluted 10 <sup>-3</sup>	MSwab™	2045	283.3	184.0	113.3	Acceptable Recovery
			2045/1	266.7	196.0	123.0	Acceptable Recovery
			2045/2	275.7	211.0	130.0	Acceptable Recovery
Streptococcus pneumoniae ATCC 6305	diluted 10 <sup>-1.5</sup>	MSwab™	2045	227.3	131.7	83.0	Acceptable Recovery
			2045/1	215.7	149.3	73.0	Acceptable Recovery
			2045/2	217.3	137.3	80.0	Acceptable Recovery
Streptococcus pneumoniae ATCC 49136	diluted 10 <sup>-1.5</sup>	MSwab™	2045	279.7	148.3	71.0	Acceptable Recovery
			2045/1	256.7	171.0	81.7	Acceptable Recovery
			2045/2	270.0	171.0	71.0	Acceptable Recovery
Enterococcus faecalis ATCC 29212	diluted 10 <sup>-1.5</sup>	MSwab™	2045	207.3	140.0	117.7	Acceptable Recovery
			2045/1	214.0	140.0	117.7	Acceptable Recovery
			2045/2	223.3	134.0	96.0	Acceptable Recovery
Staphylococcus epidermidis ATCC 12228	diluted 10 <sup>-1.5</sup>	MSwab™	2045	280.7	147.0	89.7	Acceptable Recovery
			2045/1	277.0	154.0	93.3	Acceptable Recovery
			2045/2	273.3	167.7	88.0	Acceptable Recovery
Staphylococcus aureus	diluted	MSwab™	2045	252.7	176.7	112.0	Acceptable Recovery

ATCC 29213	$10^{3.5}$		2045/1	212.7	178.3	126.7	Acceptable Recovery
			2045/2	217.0	131.7	107.3	Acceptable Recovery
			2045	174.0	112.7	80.3	Acceptable Recovery
<i>Streptococcus agalactiae</i> ATCC 9341	diluted $10^{-3}$	MSwab™	2045/1	162.3	119.7	88.0	Acceptable Recovery
			2045/2	186.0	112.7	84.0	Acceptable Recovery
			2045	271.7	198.3	165.3	Acceptable Recovery
<i>Kocuria rhizophila</i> ATCC 9341	diluted $10^{2.2}$	MSwab™	2045/1	267.7	196.3	145.7	Acceptable Recovery
			2045/2	269.3	204.7	142.3	Acceptable Recovery
			2045	267.3	221.3	167.0	Acceptable Recovery
<i>Listeria monocytogenes</i> ATCC 19114	diluted $10^{-3}$	MSwab™	2045/1	275.0	221.0	162.0	Acceptable Recovery
			2045/2	274.3	211.3	172.0	Acceptable Recovery
			2045	234.7	156.0	116.0	Acceptable Recovery
<i>Bacillus cereus</i> ATCC 10876	diluted $10^{1.5}$	MSwab™	2045/1	246.0	167.0	120.3	Acceptable Recovery
			2045/2	231.0	161.7	117.3	Acceptable Recovery
			2045	204.7	173.7	153.0	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	diluted $10^{3.5}$	MSwab™	2045/1	211.7	174.3	142.3	Acceptable Recovery
			2045/2	196.7	174.3	137.3	Acceptable Recovery
			2045	149.0	117.0	72.7	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC 6538	diluted $10^{3.5}$	MSwab™	2045/1	172.0	113.3	82.0	Acceptable Recovery
			2045/2	167.3	133.3	86.3	Acceptable Recovery
			2045	257.0	205.0	171.7	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 700698	diluted $10^{3.5}$	MSwab™	2045/1	254.3	203.0	176.3	Acceptable Recovery
			2045/2	257.7	205.3	170.3	Acceptable Recovery

**\*Table 5 - SUMMARY OF RESULTS FOR ADDITIONAL BACTERIAL RECOVERY STUDIES ON SPECIFIC STRAINS ROLL-PLATE METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 10 days	Average of CFUs recovered at time 14 days	Interpretation
<i>Staphylococcus aureus</i> ATCC 29213	diluted $10^{3.5}$	MSwab™	1210	248	349	401	Acceptable Recovery
			25C080	232	321	497	Acceptable Recovery
			43A133	240	300	521	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	diluted $10^{3.5}$	MSwab™	1210	231	264	305	Acceptable Recovery
			25C080	220	287	289	Acceptable Recovery
			43A133	212	292	215	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC 6538	diluted $10^{3.5}$	MSwab™	1210	152	175	156	Acceptable Recovery
			25C080	88	139	160	Acceptable Recovery
			43A133	256	147	165	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 700698	diluted $10^{3.5}$	MSwab™	1210	168	184	225	Acceptable Recovery
			25C080	176	240	209	Acceptable Recovery
			43A133	180	327	247	Acceptable Recovery

**\*Table 6 - SUMMARY OF RESULTS FOR ADDITIONAL BACTERIAL RECOVERY STUDIES ON SPECIFIC STRAINS ROLL-PLATE METHOD, 20-25°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 72 hrs	Interpretation
<i>Staphylococcus aureus</i> ATCC 29213	diluted $10^{3.5}$	MSwab™	1210	248	233	Acceptable Recovery
			25C080	232	303	Acceptable Recovery
			43A133	240	327	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	diluted $10^{3.5}$	MSwab™	1210	231	60	Acceptable Recovery
			25C080	220	80	Acceptable Recovery
			43A133	212	69	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC 6538	diluted $10^{3.5}$	MSwab™	1210	152	112	Acceptable Recovery
			25C080	88	194	Acceptable Recovery
			43A133	256	125	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 700698	diluted $10^{3.5}$	MSwab™	1210	168	67	Acceptable Recovery
			25C080	176	84	Acceptable Recovery
			43A133	180	75	Acceptable Recovery

**\*Table 7 - SUMMARY OF RESULTS FOR BACTERIAL OVERGROWTH STUDIES, ROLL-PLATE METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension in saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 48 hrs	Interpretation
<i>Pseudomonas aeruginosa</i> ATCC BAA-427	diluted $10^{-4}$	MSwab™	2045	39.0	0	Acceptable Recovery
			2045/1	40.3	0	Acceptable Recovery
			2045/2	36.3	0	Acceptable Recovery

**\*Table 8 - SUMMARY OF RESULTS FOR VIRAL RECOVERY STUDIES, 4-8°C**

Organism	Dilution	Product	Lot Number	Average of foci of infected cells at time 0 hrs	Average of foci of infected cells at time 24 hrs	Average of foci of infected cells at time 48 hrs	Log <sub>10</sub> decline	Interpretation
HSV 1 ATCC VR539	diluted 10 <sup>-2</sup>	MSwab™	2045	535.3	302.6	175.6	-0.48	Acceptable Recovery
			2045/ 1	505.7	354.7	195.0	-0.41	Acceptable Recovery
			2045/ 2	501.7	369.3	213.7	-0.37	Acceptable Recovery
HSV 2 ATCC VR734	diluted 10 <sup>-1</sup>	MSwab™	2045	902.0	526.0	209.7	-0.63	Acceptable Recovery
			2045/ 1	889.7	419.0	245.7	-0.56	Acceptable Recovery
			2045/ 2	954.0	486.0	275.7	-0.54	Acceptable Recovery

**\*Table 9 - SUMMARY OF RESULTS FOR VIRAL RECOVERY STUDIES, 20-25°C**

Organism	Dilution	Product	Lot Number	Average of foci of infected cells at time 0 hrs	Average of foci of infected cells at time 24 hrs	Average of foci of infected cells at time 48 hrs	Log <sub>10</sub> decline	Interpretation
HSV 1 ATCC VR539	diluted 10 <sup>-2</sup>	MSwab™	2045	535.3	199.3	142.6	-0.57	Acceptable Recovery
			2045/ 1	505.7	195.7	140.7	-0.55	Acceptable Recovery
			2045/ 2	501.7	152.0	103.3	-0.69	Acceptable Recovery
HSV 2 ATCC VR734	diluted 10 <sup>-1</sup>	MSwab™	2045	902.0	387.7	184.7	-0.69	Acceptable Recovery
			2045/ 1	889.7	414.7	203.7	-0.64	Acceptable Recovery
			2045/ 2	954.0	375.3	187.0	-0.71	Acceptable Recovery

\*Translation in other languages of words/sentences used in Table 1, 2, 3, 4, 5, 6, 7, 8, 9:

**Organism:** Organismo/Organisme/ Organismus/ Organismo/ Organisme

**Dilution:** 0,5 McFarland bacterial suspension with saline/ Diluizione Sospensione batterica 0,5 McFarland con sol. Salina/ Dilution: Suspension Bactérienne 0,5 McFarland avec sol. Saline/ Verdünnung: 0,5-McFarland-Bakteriensuspension mit Salzlösung/ Dilución: Suspensión bacteriana

**McFarland 0,5 con sol. Salina/ Diluição:** Suspensão bacterica McFarland 0,5 com sol. Salina/ Fortynding: McFarland bakteriesuspension 0,5 m. saltvandsopl / Fortynning: Organisme 0,5 McFarland bakteriell suspensjon i saltvann/ Utspådnig:0,5 McFarland bakteriell koksaltløsning

**Product:** Prodotto/ Produit/ Produkt/ Product/ Produtol Produkt

**Lot Number:** Numero di Lotto/ Numéro Lot/ Chargen-nummer/ Número Lote/ Número Lote/ Parti nr./Batchnr / Parti-nummer

**Average of CFUs recovered at time 0 hrs, 24hrs, 48hrs, 72hrs, 10 days, 14 days :** UFC medie recuperate a 0 ore, 24 ore, 48 ore, 72 ore, 10 giorni/14 giorni / Moyenne d'UFC récupérées à 0 heures, 24 heures, 48 heures, 72 heures, 10 journées, 14 journées / Durchschnitt der nach 0, 24, 48, 72 Stunden 10, 14 tage wiedergefundenen KBE/ Promedio de UFC recuperadas a las 0 horas, 24 horas, 48 horas, 72 horas, 10, 14 días / Média de UFC recuperadas a 0 hora 24 horas, 48 horas, 72 hora, 10, 14 días / Gennemsnitligt UFC genvundet efter 0, 24, 28, 72 timer, 10, 14 dage / Gjennem-snitlig antall CFU (koloni-dannende enheter) gjen-opprettet etter 0, 24, 48, 72 time, 10, 14 dager / Medelvärde av CFUs återhämtningstid vid 0, 24, 48, 72 timma, 10, 14 dagar.

**Average of foci of infected cells at time 0, 24, 48 hrs:** Média di foci di cellule infettate a , 24, 48 ore/ Moyenne de foyers de cellules infectés à 0, 24, 48 heures/ Durchschnitt der nach 0, 24, 48 Stunden infizierten Bakterienherde/ Promedio de focos de células infectadas a las 0, 24, 48 horas/ Média de fozes de células infectadas a 0, 24, 48 horas/ Gennemsnit af inficerede cellekolonier efter 0, 24, 48 timer/ Gjennomsnittlig antall knuter av infiserte celler etter 0, 24, 48 timer/ Focimedelvärde av infekterade celler vid 0, 24, 48 timmar.

**Log<sub>10</sub> decline:** Riduzione Log<sub>10</sub> / Réduction Log<sub>10</sub>/ Abnahme Log<sub>10</sub>/ Decremento Log<sub>10</sub>/ Redução Log<sub>10</sub>/ Tab Log<sub>10</sub> / Log<sub>10</sub> reduksjon/ Log<sub>10</sub> minskning

**Interpretation:** Interpretazione/ Interprétation/ Interpretation/ Interpretación/ Interpretação/ Fortolkning/ Tolkning/ Tolkning

In accordance with Clinical and Laboratory Standards Institute M40-A, viability performance is measured for each test organism at the 48 hrs time point and compared with the acceptance criteria.

In both the Roll-Plate and Swab Elution viability performance studies, Copan MSwab™ System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8°C) and room temperature (20 – 25°C). Acceptable recovery for the Roll-Plate Method is defined as  $\geq 5$  CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU. Acceptable recovery for the Swab Elution Method is defined as no more than a  $3 \log_{10}$  ( $1 \times 10^3 \pm 10\%$ ) decline in CFU between the zero-time CFU count and the CFU of the swabs after the specified holding time.

Additional time-points were tested for *Staphylococcus aureus*, ATCC 29213 and ATCC 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC 43300 and ATCC 700698.

In the Roll-Plate viability performance studies, Copan MSwab™ System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8°C) for 14 days and room temperature (20 – 25°C) for 72 hrs. Acceptable recovery for the Roll-Plate Method is defined as  $\geq 5$  CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU.

Viability performance studies also include an assessment of bacterial overgrowth at refrigerated temperatures (4 – 8°C). For the Swab Elution Method, an overgrowth assessment is made on all bacteria species tested at the 48 hrs holding time point. Overgrowth assessment using the Swab Elution Method is defined as greater than 1 log<sub>10</sub> increase in CFU between the zero-time CFU count and the holding time point. For the Roll-Plate Method, an overgrowth assessment is made with a separate analysis in which swabs are dosed with 100µl containing  $10^2$  CFU of *Pseudomonas aeruginosa* culture. Overgrowth under these conditions is defined as greater than 1 log<sub>10</sub> increase in CFU between zero-time CFU and the 48 hrs holding time point.

Copan MSwab™ System demonstrated no overgrowth based on the acceptance criteria described in Clinical and Laboratory Standards Institute M40-A.

Copan MSwab™ System was able to maintain the viability of the following organisms for at least 48 hours at both room temperature (20-25°C) and in the refrigerator (2-8°C) under the test conditions described above: Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2.

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**BILIOGRAFI**

1. Murray PR, Baron EJ, Pfaller MA, Tenoer FC, Tenover FC, Tenover FC, eds. *Manual of Clinical Microbiology*. 7th edition. Washington, DC: ASM; 1999.
2. Isenberg HD, Schoenkencht FD, Von Graeventiz A, Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, S. J. Rubin. *American Society for Microbiology*, Washington, DC, 1979.
3. Isenberg HD. (Editor in Chief). *Clinical Microbiology Procedures Handbook*. American Society for Microbiology, Washington, DC, 1992.
4. Isenberg HD. (Editor in Chief). *Essential Procedures for Clinical Microbiology*. American Society for Microbiology, Washington DC, 1998.
5. Summanen P, Baron EJ, Citron D, Strong C, Wexler HM, Finegold SM. (1993). *Wadsworth Anaerobic Bacteriology Manual*, 5th edn. Star Publishing Company, Belmont, CA.
6. ASM Cumitech publications number 10A, 13A, 17A. *American Society for Microbiology*, Washington, DC.
7. Miller JM. *A Guide to Specimen Management in Clinical Microbiology*. Second Edition. American Society for Microbiology, Washington, DC, 1999.
8. Miller JM, Holmes HT. Specimen collection, transport, and storage. In: *Manual of Clinical Microbiology*. 6th ed. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. Washington, DC: ASM; 1995:19-20.
9. Quality Control of Microbiology Transport Systems. M40-A Vol 23 No. 34. December 2003. Clinical and Laboratory Standards Institute (formerly NCCLS).
10. Sng E-H, Rajan VS, Teo K-L, Goh A-J. The recovery of *Neisseria gonorrhoeae* from clinical specimens: effects of different temperatures, transport times, and media. *Sex Trans Dis*. 1982; 9:74-78.
11. Sun Y, Taylor T, Williams L, Sautter RL. Comparison of bacterial viability using both the EZ brand collection and transport system with the Difco swab transport pack. Presented at: 96th ASM General Meeting. 1996; Washington DC. Abstract C35.
12. Arbique JC, Forward KR, LeBlanc J. Evaluation of four commercial transport media for the survival of *Neisseria gonorrhoeae*. *Diagnostic Microbiology and Infectious Disease*. 2000; 36:163-168.
13. Perry JL. Effects of temperature on fastidious organism viability during swab transport. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-55.
14. Wilson DA, Tuohy MS, Procop GW, Hall GS. Effects of storage on the recovery of bacteria from three swab transport systems: BD CultureSwab, BD Culturette and Starplex StarSwab II. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-61.
15. Arbique J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies described in NCCLS document M40-P Quality Control of Microbiology Transport Devices. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-40.
16. Mitchell E, Berman M, Ginocchio CC. Evaluation of two new Liquid Stuart transport systems: Platinum StarSwab II (Starplex Scientific) and BBL CultureSwab (Becton Dickinson). 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-74.
17. Perry JL, Matthews JS. Compliance of two popular swab transport systems with performance standards detailed by the new NCCLS Proposed Standard, M40-P. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-42.
18. Robinson A, Gruver ML. Comparison of bacterial survival in two transport systems stored at room temperature and refrigerator temperatures. 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-69.
19. Human RP, Jones GA. Evaluation of 4 transport systems against a published standard. 104th General Meeting of the American Society for Microbiology. 2004; New Orleans, LA. Abstract C-161.
20. Human RP, Jones GA. Evaluation of swab transport systems against a published standard. *J Clin Pathol* 2004; 57:762-763.
21. Arbique J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies for anaerobic organisms described in NCCLS document M40-P, Quality Control of Microbiology Transport Devices. 13th European Congress of Clinical Microbiology and Infectious Disease (ECCMID). 2003; Glasgow, UK. Abstract P-652.
22. Marler LM, Siders JA, Allen SD. *Direct Smear Atlas, A Monograph of Gram-Stained Preparations of Clinical Specimens*. Lippincott Williams and Wilkins, 2001.
23. Rotimi VO, Yakubu Z, Abudu OO, Banjo TO. Direct Gram's stain of vaginal discharge as a means of diagnosing bacterial vaginosis. *Journal of Medical Microbiology*. 1991 Vol 35, Issue 2 103-106.
24. Spiegel CA, Amsel R, Holmes KK. Diagnosis of bacterial vaginosis by direct gram stain of vaginal fluid. *J. Clin Microbiol.* 1983 Jul;18 (1):170-177.
25. Benavides MI, Moncada X, Rodriguez B, Castillo C. Gonococcal urethritis in men: clinical experience in 1978-1988. *Rev Med Chil*. 1992 Oct;120(10):1140-3.
26. Mayaud P, Msuya W, Todd J, Kaatano G, West B, Begkoyian G, Grosskurth H, Mabey D. Rapid assessment in Rwandan refugee camps in Tanzania. *Genitourin Med*, 1997 Feb;73 (1):33-8
27. Lynne S. *Garcia Clinical Microbiology Procedures Handbook, Third Edition* 2010 ASM Press 3.2.1
28. AE Greenberg, LS Clesceri, and AD Eaton. 9215 heterotrophic plate count. In: *Standard Methods for the Examination of Water and Waste Water*. 18th ed. Washington, DC APHA; 1992: 9-33-9-34
29. Henry D. *Essential procedures for Clinical Microbiology* 1998 ASM PRESS 477-478
30. Steven Specter, Richard L. Hodinka, Stephen A. Young *Clinical Virology Manual* Third edition, 384-409
31. Fleming D. *Biological Safety: Principles and Practices*. January 2000. ASM, Washington DC.
32. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA. <http://www.cdc.gov/od/ohs/symp5/jyrtxt.htm>.
33. Richardson JH. *Biosafety in Microbiological and Biomedical Laboratories*. December 1994. Diane Publishing Company.
34. Hansen DJ. *Healthcare, Laboratories and Biosafety*. Vol 2., 1992. CRC Press.

This page is intentionally blank

This page is intentionally blank

Index of Symbols/Table des Symboles/Symbole/Tabla de Simbolos/Tabella dei Simboli/Quadro de Simbolos/Symboltabell

Symbol/Symbole/Simbolo	Meaning/Signification/Bedeutung/Significado/Significato/Significato/Significato/Betydning
	Manufacturer/Fabricant/Hersteller/Fabricante/Fabbricante/Fabricante/ Fabrikant/Tillverkare
	In vitro diagnostic device /Dispositivo diagnostico in vitro/ dispositif de diagnostic in vitro/ Diagnosegerät in vitro/ dispositivo de diagnóstico in vitro/ diagnostisk anordning in vitro/ diagnostisk utstyr in vitro
	Identification number of notified body/Identification de l'organisme notifié/Identifizierung der benannten Stelle/Identificación del organismo notificado/Identificazione dell'organismo notificato/Identificação do organismo notificado/ Identifikasjonsnummer ved påvist organisme/Identifiseringsnummer av sertifisert myndighet
	Sterilized using ethylene oxide/Sterilizzato usando ossido di etilene/ Sterilisiert mit Äthylenoxid/ Esterilizado usando óxido de etileno/ Esterilizado por óxido de etileno/ stérilisé à l'aide d'oxyde d'éthylène/ Sterilisert ved hjelp av etylenoksid/ Steriliserad med etylenoxid/ Steriliseret med ethylenoxid
	Do not reuse/Ne pas réutiliser/Nicht zur Wiederverwendung/No reutilizar/Non riutilizzare/Não voltar a usar/ Må ikke brukes på nytt/Får inte återanvändas
	Catalogue number/Référence du catalogue/Bestellnummer/Número de catálogo/Numero di catalogo/Referência do catálogo/ Katalognummer/Katalognummer
	Temperature limitation/Limites de temperature/Temperaturbegrenzung/Limites de temperatura/Limiti di temperatura/Limites de temperatura/Temperaturgreenser/Temperaturgränser
	Use by/Utiliser jusque/Verwendbar bis/Fecha de caducidad/Utilizzare entro/Prazo de validade/ Må brukes innen/ Ska användas innan
	Consult Instructions for Use/Consulter les instructions d'utilisation/Gebrauchsanweisung beachten/Consulte las instrucciones de uso/Consultare le istruzioni per l'uso/Consultar as instruções de utilização/ Se instruksjoner for bruk/Se bruksanvisningen
	Peel/Décoller/Abziehen/Desprender/Strappare per aprire/Destacável/Dra/Dra för att öppna
	Batch code (Lot)/Code de lot (Lot)/Chargencode (Chagenbezeichnung)/Código de lote (Lote)/Codice del lotto (partita)/ Código do lote (Lote)/ Lot nummer (parti)/Serienummer(parti)
	Contains sufficient for <n> tests/Contenu suffisant pour <n> tests/Ausreichend für <n> Tests/Contenido suficiente para <n> pruebas/Contenuto sufficiente per <n> test/Contémto suficiente para <n> testes/ Innhold tilstrekkelig for <n> test/Innehåller tillräckligt för<n> tester

Copan Flock Technologies Srl  
Via F. Perotti, 18  
25125 - Brescia, Italy

Copan Flock Technologies Srl  
Via F. Perotti, 18  
25125 - Brescia, Italy  
Tel +39 030 3666100  
Fax +39 030 2659932  
Email: info@copanflock.com  
Website: www.copaninnovation.com

North American Distributor:  
Copan Diagnostics Inc.  
26055 Jefferson Avenue  
Murrieta, CA 92562 USA  
Tel: 951-696-6957  
Fax: 951-600-1832

E-mail: copanswabs@aol.com  
Website: www.copanusa.com