UTM[®]

Instructions for Use

CE **IVD**

English	3
Italiano	9
Español	16
Deutsch	22
Français	28
Português	34
БЪЛГАРСКИ	40
Česky	48
Dansk	55
ΕΛΛΗΝΙΚΑ	63
Eesti Keel	71
Hrvatski	78
Latviešu	85
Lietuvių K.	92
Magyar	100
Nederlands	106
Norsk	114
Polski	120
Română	128
Slovenčina	135
Slovenščina	142
Srpski	148
Suomi	155
Svenska	162
Türkçe	169
УКРАЇНСЬКА	176

Copan Universal Transport Medium (UTM-RT®) System Instructions for use

INTENDED USE

Copan Universal Transport Medium (UTM-RT®) System is intended for the collection and transport of clinical specimens containing viruses, chlamydiae, mycoplasma or ureaplasma from the collection site to the testing laboratory. UTM-RT® can be processed using standard clinical laboratory operating procedures for viral, chlamydial, mycoplasma and ureaplasma culture.

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of infections caused by viruses, chlamydiae, mycoplasmas or ureaplasmas involves collection and refrigerated transport of biological specimens. Using the UTM-RT® System, the collected specimen can be stored for up to 48 hours at 2-25°C. The UTM-RT® consists of a Hanks' Balanced Salt Solution (HBSS) enriched with proteins and sugars with a neutral pH and pH indicator. The medium contains some antibiotics and antimycotics to inhibit overgrowth of bacteria and yeasts, maintain cellular integrity and encourage preservation of viruses and chlamydiae.

Copan UTM-RT® System medium is provided in labeled screw-cap tubes designed for transport of the clinical sample. Copan UTM-RT® System is also supplied as a sample collection kit that comprises a package which contains one screw-cap tube of UTM-RT® medium and a peel pouch incorporating one or two sterile specimen collection swabs. A range of UTM-RT® sample collection kits are available which incorporate different types of shaft swabs which facilitate the collection of specimens from different sites of the patient as described below in the Directions for Use section. Once a swab sample is collected it should be placed immediately into the transport tube where it comes into contact with transport medium. To maintain optimum microorganism viability, transport the specimens to the laboratory as soon as possible.

After collection, the specimen should be stored at 2-25°C and processed within 48 hours.

If delivery and processing exceed 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder.

REAGENTS

The UTM-RT® formulation includes proteins for virus stabilization¹⁷, antibiotics and antimycotics to prevent overgrowth of bacterial and fungal flora and a buffer solution to maintain a neutral pH.

Components
Sucrose
HBSS solution
Bovine serum albumin
Buffered solution
Gelatin
Amino acids
Antibiotics
Phenol Red

pH 7,3 ± 0,2 a 2+25°C

REQUIRED MATERIALS BUT NOT PROVIDED

Materials suitable for isolation, differentiation and culture of viruses, chlamydiae, mycoplasmas and ureaplasmas.

STORAGE

The product must be stored in its original packaging at a temperature between 2 and 25°C until the time of use. Do not overheat or freeze prior to use.

SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION

Specimens for virus, chlamydia, mycoplasma or ureaplasma investigation should be collected and handled following published manuals and guidelines^{2,3,4,7,9,10,11}. To maintain optimum viability, transport the specimen to the laboratory as soon as possible. After the collection, the specimen should be stored at 2-25°C and processed within 48 hours. If there will be a long delay exceeding 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder.

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations^{8,11,12}. Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

MATERIALS SUPPLIED

Copan UTM-RT® System is available in the format of kit (tube with specimen collection swabs) or tube only, the system is contained in box and boxes are contained in an external carton (Pack size are provided in table one).


Copan UTM-RT® System includes a screw-cap tube with conical or round bottom containing 1.5ml, 1ml, 2ml, 2.5ml, 3ml or 10ml of light orange-red transport medium. For all the product codes, excluding 3U008N, 331C, 348C and 3E056N, the tube contains three glass beads that allow homogenizing the specimen by vortexing.

UTM-RT® System tubes of transport medium are supplied alone or in a kit format. Each kit unit consists of a package containing: a pre-labeled screw-cap tube and a peel pouch incorporating one or two sterile specimen collection swabs (e.g. applicator swabs polyester tipped, flexible size applicator swab with flocked nylon fiber). For more details on available configurations please refer to Table 1.

The collection swab applicators FLOQSwabs® and CLASSIQSwabs™ provided with UTM-RT® have a molded breakpoint in the shaft of the applicator. For the line of FLOQSwabs® applicator the molded breakpoint show a colored mark.

UTM-RT® tube capture caps have an internal moulded design that is able to capture the swab shaft when it is broken off into the tube and the cap is closed. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a moulded docking receptacle in the cap (Fig.1). In the testing laboratory when the cap is unscrewed and removed, the swab applicator is attached to the cap. This feature allows the operator to conveniently remove the swab from the transport tube.

Fig 1. Capture of broken swab applicator stick by UTM-RT® tube cap

LIMITATIONS

1. Because calcium alginate swabs are toxic for many enveloped viruses⁵ and may interfere with immunofluorescence tests², they should not be used for specimen collection.
2. Wooden shaft swabs may contain toxins and formaldehydes^{2,6} and should not be used.
3. The product codes without beads (REF 3U008N, 331C, 348C and 3E056N) are not suitable for use with mucous or particularly viscous specimens.
4. UTM-RT® kits are intended to be used with the medium tubes and swabs provided in the kit. The use of tubes of medium or swabs from another source could affect the performance of the product.

WARNINGS AND PRECAUTIONS

1. Single-use device for professional in vitro diagnostic use.
2. UTM-RT® System is ready for use and requires no further preparation.
3. Do not use beyond the expiry date.
4. Do not use the UTM-RT® medium for premoistening or prewetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
5. Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.⁷
6. Do not re-sterilize unused swabs.
7. Do not re-pack.
8. Not suitable to collect and transport microorganisms other than viruses, chlamydiae, mycoplasma and ureaplasma.
9. Do not ingest the medium.
10. Specimens for the search of viruses, chlamydiae, mycoplasmas and ureaplasmas must be collected and handled using personal protective equipment against biological risk according to published manuals and guidelines^{1,4,6,7,8,17}.
11. Repeated freezing and thawing of specimens may reduce the recovery of viable organisms.
12. Do not use UTM-RT® if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the color of the medium has changed from light orange-red, (4) the swab pouch is open, or (5) there are other signs of deterioration.
13. The use of this product in combination with diagnostic kits or instruments must be validated by the user prior to use.
14. Do not bend or shape the swab before the collection of the specimen. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
15. Due to the design of the flexible minitip, the swab will coil when placed in the tube. Therefore, if necessary, remove the swab from the tube, use caution and observe adequate biohazard precaution to protect the operator and the environment in case of splash.
16. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.

INSTRUCTION FOR USE

UTM-RT® System is ready for use and requires no further preparation. It is available in the various configurations listed in **Table 1**.

REF	PRODUCT DESCRIPTION		PACK SIZE	SAMPLING SITES*	PRODUCT DESIGNED FOR HAVING CAPTURE CAP FEATURE (FOR UTM KIT ONLY)
	TUBE	SWAB			
302C; 302C.LC	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	Two regular size applicator swabs polyester tipped with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
305C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO
305CMH	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*

328C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swabs polyester tipped with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
306C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	NO
CA302MH; 346C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	NO*
307C; CA303MH	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO
357C	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
321C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip size applicator swab with flocked nylon fiber tip with breaking point One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Flexible Minitip: Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites Regular: Nose, throat, vagina, rectum, faeces and wounds	NO*
323C	2 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	NO*
338C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Minitip: eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites Regular: Nose, throat, vagina, rectum, faeces and wounds	NO*
347C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One mini flexible applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
355C; 355CW	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
356C	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab polyester tipped with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	NO*

358C	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
359C; 359CE.A	1 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
360C	1 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
361C	1 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
365C	2 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
366C; 3C004N	2 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
3U006N	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
3U008N	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom – No Beads.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
330C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box		NOT APPLICABLE Capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service
331C; CA339MH	10 mL of UTM-RT® medium in 25x90 mm screw-cap tube with internal shaped conical bottom.	NA	15 tubes per package 6 x 15 tubes per box		NO
348C	2.5 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom – No Beads.	NA	50 tubes per package 6 x 50 tubes per box		NO

349C	1,5 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box	NOT APPLICABLE Capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service
3U044N	3 mL of UTM-RT® medium in 16x100 mm screw- blue cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box	NO
3E056N	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom – No Beads.	NA	50 tubes per package 6 x 50 tubes per box	NO
350C; 350CV	1 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box	NOT APPLICABLE Capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service
353C	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box	NOT APPLICABLE Capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service

Table 1: product description

Not all the product codes (REF) are salable in all countries. Please contact Copan Customer care service for product codes availability for a specific country.
NO* = The swab could be accidentally and not firmly captured and might fall down.

The sampling sites Table is just a suggested table. Performance testing with Copan UTM-RT® System was conducted using laboratory strains. Performance testing was not conducted using human specimens. Please refer to your internal procedures to choose the most appropriate device for the specific sampling site.


Specimen Collection

Proper collection of the specimen from the patient is a crucial aspect for successful isolation and identification of infectious organisms. Specimens should be collected as soon as possible after the clinical onset of disease. Highest viral titers are present during the acute illness.


UTM-RT® in kit

1. Open the UTM-RT® kit package and remove the medium test tube and the internal bag containing the sterile swab.
2. Take the sterile swab out of its bag and collect the clinical specimen; to prevent the risk of contamination, make sure that the swab tip comes into contact with the collection site only. **NOTE:** Do not bend the swab before the collection of the specimen. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
3. After collecting the specimen, insert the swab into the test tube until the breakpoint is level with the test tube opening.
4. Bend the swab shaft at a 180 degrees angle to break it off at the breaking point. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
5. Discard the broken handle part of the swab shaft into an approved medical waste disposal container.
6. Screw the cap back onto the test tube and hermetically seal it.

Fig 1. Collection swab showing breakpoint indication line and area for holding the applicator

Sterile gloves and protective clothing and eyewear should be worn when collecting and handling microbiology specimens and care should be taken to avoid splashes and aerosols when breaking the swab stick into the tube of medium. During sample collection when handling the swab applicator, the operator must not touch the area below the colored breakpoint indication line; that is the area from the line to the tip of the nylon flocked swab, as this will lead to contamination of the applicator shaft and the culture thus invalidating the test results.

UTM-RT® in bulk

1. Aseptically remove cap from tube taking care not to spill the medium.
2. Aseptically place vesicle aspirates²⁶, corneal or conjunctival scrapings²⁸, small pieces of tissue or stool³⁰ samples into the tube with UTM-RT® medium.
3. Screw the cap back onto the test tube and hermetically seal it.
4. Identify the tube containing the specimen.
5. Send to the laboratory for immediate analysis.

If processing is delayed (over 48 hours), the specimens must be frozen at -70°C or colder.

DISPOSAL

Waste must be disposed of in compliance with local legislation.

QUALITY CONTROL

The UTM-RT® lots are tested for microbial contamination, toxicity for the cellular lines used for the viral cultures and the ability to maintain the viability of viral, chlamydia and mycoplasma strains for up to 48 hours at 2-25°C in accordance with the methods described in CLSI M40-A2⁵.

RESULTS AND PERFORMANCE

The results obtained largely depend on proper and adequate specimen collection as well as the promptness with which the specimens are transported to the laboratory and analysed.

Viability studies were performed using Copan UTM-RT® with a panel of representative strains of the various families supported by the UTM-RT®. The swabs that accompany each transport system were directly inoculated in triplicate with 100 µl of organism suspension using laboratory strains. Performance testing was not conducted using human specimens. Subsequently, the swabs were inserted in the respective test tubes containing the transport medium and stored for 0 and 48 hours at 2-6°C and at controlled room temperature (20-25°C). At the time of processing, each swab was vortexed for 20 seconds and removed from its transport medium test tube, then, an aliquot of the suspension was inoculated into the cellular line (200 µl) or into the appropriate culture medium^{6,15}. All the cultures were processed using the standard laboratory culture technique^{6,15}. Organism viability was determined by fluorescent cell counting for viral and chlamydia strains and CFU counting for mycoplasma and ureoplasma strains. The acceptability limits for time zero and for 48 hours were defined in accordance with the regulations M40-A2⁵.

Copan UTM-RT® System preserved the viability of all the organisms tested for 48 hours at both controlled room temperature and in the refrigerator in the above described test conditions. The organisms evaluated and the results obtained are given in the table below.

Organism	ATCC® number	% of reduction of fluorescing infected cells after 48 hours time 2-6°C	% of reduction of fluorescing infected cells after 48 hours time 20-25°C
Herpes Simplex Virus Type 1	ATCC® VR-539	≤50%	≤50%
Herpes Simplex Virus Type 2	ATCC® VR-734	≤50%	≤50%
Respiratory Syncytial Virus	ATCC® VR-1580	≤50%	≤50%
Coxsackie B1 Virus	ATCC® VR-28	≤50%	≤50%
Chlamydia trachomatis	ATCC® VR-880	≤50%	≤50%
Influenza A	ATCC® VR-1679	≤50%	≤50%

Cytomegalovirus	ATCC® VR-977	≤50%	≤50%
Varicella-zoster virus	ATCC® VR-1367	≤50%	≤50%
Chlamydia pneumoniae	ATCC® VR-1360	≤50%	≤50%
Adenovirus	DHI 20-4740010	45%	47%
Parainfluenza virus type 3	DHI 20-4770010	0%	41%
Echovirus type 30	ATCC® VR-322	13%	24%
Organism	ATCC® number	% of reduction of CFUs after 48 hours time 2-6°C	% of reduction of CFUs after 48 hours time 20-25°C
Mycoplasma hominis	ATCC® 23114	35%	25%
Ureaplasma urealyticum	ATCC® 27618	0%	6%
Mycoplasma pneumoniae	ATCC® 15531	50%	50%

TABLE OF SYMBOLS
See the table of symbols at the end of the instructions for use.

NOTES FOR THE PROFESSIONAL USER
In the event that a serious incident occurs in relation to this device, it must be reported to the manufacturer (see the contacts at the end of the Instructions for Use) and the competent authority in the country where the user and/or patient is located.

REVISION HISTORY

Last Revision No.*	Release date	Changes made
03	03-2025	Addition of CH-REP symbol in section TABLE OF SYMBOLS and update of MATERIALS SUPPLIED section (product codes discontinued)

* Should you need earlier revisions, contact Copan Customer Service.

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BIBLIOGRAPHY

1. Gary W. Procop and Elmer W. Koneman, 2016. Color Atlas and Textbook of Diagnostic Microbiology, Seventh edition. Wolters Kluwer Health.
2. James H. Jorgensen, Michael A. Pfaller, Karen C. Carroll, Guido Funke, Marie Louise Landry, Sandra S. Richter, David W. Warnock, 2015. Manual of Clinical Microbiology, 11th Edition. ASM, Washington, DC.
3. James Versalovic, Karen C. Carroll, Guido Funke, James H. Jorgensen, Marie Louise Landry, David W. Warnock, 2011. Cumitech 15A. Laboratory Diagnosis of Viral Infections. ASM, Washington, DC.
4. Patricia Tille. 2014. Bailey & Scott's Diagnostic Microbiology, 13th Edition. Laboratory Medicine.
5. Clinical and Laboratory Standards Institute (CLSI), 2014. M40-A2 Quality Control of Microbiological Transport Systems; Approved Standard-Second Edition.
6. Clinical and Laboratory Standards Institute (CLSI), 2006. M41-A Viral Culture; Approved guidelines.
7. Wardford, A., M. Chernesky, and E. M. Peterson, 1999. Cumitech 19A, Laboratory Diagnosis of Chlamydia trachomatis Infections. ASM, Washington DC.
8. 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
9. J. Michael Miller, Shelley A. Miller, 2017. A Guide to Specimen Management in Clinical Microbiology, Third Edition. ASM, Washington DC.
10. Centers for Disease Control and Prevention (CDC), 2016. Guide for Shipping Infectious Substances.
11. Centers for Disease Control and Prevention (CDC), 2009. Biosafety in Microbiological and Biomedical Laboratories 5th Edition.
12. World Health Organization 2015. Guidance on regulations for the Transport of Infectious Substances 2015 – 2016.
13. Centers for Disease Control and Prevention (CDC), 2002. Screening Tests to Detect Chlamydia trachomatis and Neisseria gonorrhoeae Infections - 2002.
14. Centers for Disease Control and Prevention (CDC), 2014. Recommendations for the Laboratory-Based Detection of Chlamydia trachomatis and Neisseria gonorrhoeae - 2014.
15. European Collection of Authenticated Cell Cultures (ECACC). Fundamental Techniques in Cell Culture Laboratory Handbook, 3rd edition (Protocol 9 – Detection of Mycoplasma by Culture Isolation).
16. J.B. Mahony, M.A. Chernesky Effect of Swab Type and Storage Temperature on the Isolation of Chlamydia trachomatis from Clinical Specimens Journal of Clinical Microbiology, Nov. 1985, p. 865-867.
17. S. Spector, R. L. Hodinka, S. A. Young. Clinical Virology Manual, fifth edition, 2016.
18. Peter Daley, Santina Castriaciano, Max Chernesky, Marek Smieja. Comparison of Flocked and Rayon Swabs for Collection of Respiratory Epithelial Cells from Uninfected Volunteers and Symptomatic Patients. Journal of Clinical Microbiology, June 2006.
19. K. Loens, L. Van Heirstraeten, S. Malhotra-Kumar, H. Goossens, and M. Ieven. Optimal Sampling Sites and Methods for Detection of Pathogens Possibly Causing Community-Acquired Lower Respiratory Tract Infections. Journal of Clinical Microbiology, Jan. 2009.
20. Richard Garceau MD, Danielle Leblanc RT, Louise Thibault MD, Gabriel Girouard MD, Manon Mallet PhD. Herpes simplex virus type 1 is the leading cause of genital herpes in New Brunswick. Can J Infect Dis Med Microbiol 2012.
21. Marek Smieja, Santina Castriaciano, Susan Carruthers, Geoffrey So, Sylvia Chong, Kathy Luinstra, James B. Mahony, Astrid Petrich, Max Chernesky, Mario Savarese, and Daniele Triva. Development and Evaluation of a Flocked Nasal Midturbinate Swab for Self-Collection in Respiratory Virus Infection Diagnostic Testing. Journal of Clinical Microbiology, Sept. 2010.
22. Donghyok Kwon, Kyeongcheol Shin, Mihwa Kwon, Hee-Bok Oh, Chun Kang, Joo-Yeon Lee. Development and Evaluation of Rapid Influenza Diagnostic Test for the Pandemic (H1N1) 2009 Virus. JCM Accepts, published online ahead of print on 27 October 2010.
23. Bhupesh K. Prusty, Christine Siegl, Petra Hauck, Johannes Hain, Suvi J. Korhonen, Eija Hiltunen- Back, Mirja Puolakkainen, Thomas Rudel. Chlamydia trachomatis Infection Induces Replication of Latent HHV-6. Plos One, April 2013.
24. Ellen Vancutsem, Oriane Soetens, Maria Breugelmans, Walter Foulon, Anne Naessens. Modified Real-Time PCR for Detecting, Differentiating, and Quantifying *Ureaplasma urealyticum* and *Ureaplasma parvum*. The Journal of Molecular Diagnostics, Vol. 13, No. 2, March 2011.
25. Turkiya Al-Siyabi, Khalifa Binkhamis, Melanie Wilcox, Sallene Wong, Kanti Pabbaraju, Raymond Tellier, Todd F. Hatchette, and Jason J. Le Blanc. A cost-effective real-time PCR for the detection of adenovirus from viral swabs. Al-Siyabi et al. Virology Journal 2013.
26. Eidgenössisches Departement für Verteidigung, Bevölkerungsschutz und Sport VBS; SAMPLING AND DISPATCH OF SAMPLES - VIROLOGY Ident-Nr/Vers. 10014156143/01 Aktenzeichen:323-FachinfoViro/2016.
27. Matthew L. Faron et al (2016); Multicenter Evaluation of Meridian Bioscience HSV 1&2 Molecular Assay for Detection of Herpes Simplex Virus 1 and 2 from Clinical Cutaneous and Mucocutaneous Specimens; Journal of Clinical Microbiology Volume 54 Number 8
28. Al-Siyabi et al. (2013); A cost effective real-time PCR for the detection of adenovirus from viral swabs. Virology Journal 10:184.
29. Le Roy C, Pereyre S, Hénin N, Bébéar C. 2017. French prospective clinical evaluation of the Aptima Mycoplasma genitalium CE-IVD assay and macrolide resistance detection using three distinct assays. J Clin Microbiol 55:3194–3200. <https://doi.org/10.1128/JCM.00579-17>.
30. Minodier et al. (2017) Clinical and virological factors associated with gastrointestinal symptoms in patients with acute respiratory infection: a two-year prospective study in general practice medicine BMC Infectious Diseases 17:729 DOI 10.1186/s12879-017-2823-9

Table of symbols / Tabella dei simboli / Tabla de símbolos / Tabelle der verwendeten Symbole / Tableau des symboles / Tabela de símbolos / Таблица със символи / Tabulka symbolů / Tabel over symboler / Πίνακας συμβόλων / Tingmärkide tabel / Tablica simbolă / Simbolu tabula / Simbolių lentelė / Szimbólumok táblázat / Tabel met symbolen / Symboltabel / Tabela simboli / Tabelul simbolurilor / Tabuľka symbolov / Tabela simbolov / Tabela simbola / Symbolitaulukko / Symbolföreteckning / Semboller tablosu / Таблица символів

Symbol / Simbolo / Simbolo / Symbol / Symbole / Simbolos / Символ / Symbol / Symbol / Σύμβολο / Tingmärk / Symboli / Simbol / Szimbólum / Simbolis / Simbols / Symbool / Symbol / Symbol / Simbol / Symbol / Symbol / Symboll / Символ	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Значение / Význam / Betyder / Σημασία / Tähendus / Merkitys / Značenje / Jelentés/ Reikšmė / Nozīme / Betekenis / Betydning / Znaczenie / Semnificație / Význam / Pomen / Značenje / Betydelse / Anlam / Значення
	Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante /Производител / Výrobce / Producent / Κατασκευαστής / Tootja / Valmistaja / Proizvođač / Gyártó/ Gamintojas / Ražotājs / Fabrikant / Produsent / Producent / Producător / Výrobca / Proizvajalec / Proizvodač / Tilverkare / Üretici / Виробник
IVD	In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro/ Dispositivo de diagnóstico in vitro / Инвитро диагностично изделие / Diagnostický prostředek pro in vitro diagnostiku / Udstyr til in vitro-diagnostik / Διαγνωστικό προϊόν in vitro / In vitro diagnostika seade / In vitro -diagnostikkalaita / In vitro -diagnostiikki proizvod / In vitro diagnostikai eszköz / „In vitro“ diagnostikos prietaisas / In vitro diagnostikas ierīce / Hulpmiddel voor in-vitrodiagnostiek / In vitro-diagnostisk enhet / Urządzenie do diagnostyki in vitro / Dispozitiv de diagnostic in vitro / Diagnostická pomůcka in vitro / In vitro diagnostiční připoček // In vitro dijagnostičko medicinsko sredstvo / Enhet för in vitro-diagnostik / In vitro tani cihazı / Виріб для діагностики in vitro
MD	Medical Device / Dispositivo Medico / Producto sanitario / Medizinprodukt / Dispositif médical / Dispositivo Médico / Медицинско изделие / Zdravotnický prostředek / Medicinsk udstyr / Ιατροτεχνολογικό προϊόν / Meditsiinisade / Lääkinläännelinen laite / Medicinski proizvod / Orvostechnikai eszköz / Medicinos prietaisas / Mediciniska ierīce / Medisch hulpmiddel / Medisinsk utstyr / Wyrób medyczny / Dispozitiv medical / Zdravotnícka pomůcka / Medicinski pripomoček / Medicinski aparat / Medicinteknisk enhet / Tibbi Cihaz / Медицин пристрій
UDI	Unique Device Identifier / Identificativo unico del dispositivo / Identificador único de dispositivo / Einmalige Produktkennung / Identifiant dispositif univoque / Identificador de dispositivo único / Унікален ідентифікатор на изделието / Jedinečný identifikátor zařízení / Unik enheds-id / Αποκλειστική ταυτοποίησης ιατροτεχνολογικού προϊόντος / Seadme kordumatu identifitseerimistunnus / Yksilöivä laitetunniste / Jedinstveni identifikator proizvoda / Egyedi eszközazonosító / Unikalus jreginio identifikatorius / Unikālais ierices identifikators / Unieke hulpmiddelidentificatie / Identifikator for Unique-enhet / Unikatowy numer urządzenia / Identificator unic dispozitiv / Unikátny identifikátor pomůcky / Edinstveni identifikator naprave / Jedinstveni identifikator uređaja / Unik enhetsidentifiare / Benzersiz Cihaz Tanımlayıcı / Унікальний ідентифікатор пристрою
CE	CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE/ Marcação CE / CE маркировка / Označení CE / CE-mærkning / Σήμανση CE / CE-märgis / CE-merkintä / Oznaka CE / CE-jelölés/ CE ženklas / CE markējums / CE-markering / CE-merking / Oznakowanie CE / Marcaj CE / Označenje CE / Oznaka CE / CE oznaka / CE-märkning / CE işareti / Маркування CE
CE 0123	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Идентификационный номер на нотифицирующа орган / Identifikační číslo notifikovaného subjektu / Identifikationsnummer af bemyndiget organ / Αριθμός ταυτοποίησης κοινοποιημένου οργάνου / Teavitatud asutuse identifitseerimisnumber / Ilmoitetun laitoksen tunnusnumero / Identifikacijski broj obavještenog tijela / A bejelentett szervezet azonosító száma/ Notifikuotosios įstaigos identifikavimo numeris / Pilnvarotās iestādes identifikācijas numurs / Identificatienummer van de aangemelde instantie / Identifikationsnummeret til det tekniske kontrolorganet / Numer identyfikacyjny jednostki notyfikowanej / Numărul de identificare al organismului notificat / Identifikačné číslo notifikovanej osoby / Identifikacijska številka prilaščenega organa / Identifikacioni broj nadležnog tela / Identifikationsnummer för anmält organ / Onaylı kuruluş tanımlama numarası / Ідентифікаційний номер уповноваженого органу
STERILE EO	Sterilized using ethylene oxide / Sterilizzato usando ossido di etilene / Esterilizado por óxido de etileno / Sterilisiert mit Äthylenoxid / Stérilisé à l'aide d'oxyde d'éthylène / Esterilizado utilizando óxido de etileno / Стерилизовано с етилен оксид / Sterilizováno pomocí ethylenoxidu / Steriliseret med ethylenoxid / Αποστειρωμένο με τη χρήση οξειδίου του αιθυλίου / Steriliseeritud, kasutades etüleenoksiidi / Steriloitu etyleenioksidilla / Sterilizirano etilen-oksidad / Etilén-oxidall sterilizálva / Sterilizuota etileno oksidu / Sterilizēts, izmantojot etilēnoksidu / Gesteriliseerd met ethyleenoxide / Steriliseret med etylenoksid / Sterylizowany tlenkiem etylenu / Sterilizat cu oxid de etilenă / Sterilizované etýlenoxidom / Sterilizirano z etilenoksidom / Sterilisano etilen-oksidad / Steriliserad med etylenoxid / Etilen oksit kullanılarak sterilize edilmiştir / Стерилизований завдяки етиленоксиду

STERILE R	STERILIZED using irradiation / Sterilized using radiation / Sterilized using radiations ionizantes / Esterilizado por radiación ionizante / Стерилизовано чрез облъчване / Sterilizováno pomocí ozařování / Steriliseret med stråling / Αποστειρωμένο με τη χρήση ακτινοβολίας / Steriliseeritud, kasutades kiintust / Steriloitu säteilyttämällä / Sterilizirano zračenjem / Sugárzással sterilizálva / Sterilizuota švintināt / Sterilizēts, izmantojot apstarošanu / Gesteriliseerd met bestraling / Steriliseret ved stråling / Sterylizowany promieniowaniem / Sterilizat prin iradiere / Sterilizované ožiením / Sterilizirano z ionizirajočim sevanjem / Sterilisano zračenjem / Steriliserad med strålning / İşinlama kullanılarak sterilize edilmiştir / Стерилизований завдяки опроміненню
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Не използвайте повторно / Výrobek nepoužívejte opakovaně / Må ikke genbruges / Μην επαναχρησιμοποιείτε / Ühekordne kasutus / Ei saa kasutada uudeleen / Ne upotrebjavajte ponovno / N használja újra / Nenaudoti pakartotinai / Nelietot atkārtoti / Niet hergebruiken / Skal ikke brukes flere ganger / Nie używać ponownie / A nu se reutiliza / Nepoužívať opakovane / Ne uporabljajte ponovno / Ne koristiti ponovo / Får ej återanvändas / Yeniden kullanmayın / Не використовувати повторно
	Single sterile barrier system / Sistema a barriera sterile singola / Sistema de barrera estéril único / Einfaches Sterilbarriersystem / Système de barrière stérile unique / Sistema de barreira estéril simples / Единична стерилна бариерна система / Systém jedné sterilní bariéry / Enkelt sterilt barriersystem / Movó ósúttmá stéirle / Ühekordne steriilne barjäärsüsteem / Yksittäinen steriili estejärjestelmä / Jednostruki sterilni sustav barijera / Egyszeres steril gátrendszer / Vienkartinė steriliaus barjero sistema / Vienas sterilas barjeras sistēma / Enkel steril barrièresysteem / Sterilt barriersystem til engangsbruk / System pojedynczej bariery sterylnej / Sistem cu barieră sterilă unică / Systém s jednou sterilnou bariérou / Enojni sterilni pregradni sistem / Sistem jednostruktne sterilne barijere / Enkelt sterilt barriärsystem / Tek kullanımlık steril ambalaj sistemi / Система з одним стерильним бар'єром
	Do not resterilize / Non risterilizzare / No reesterilizar / Nicht resterilisieren / Ne pas restériliser / Não reesterilizar / Не стерилизуйте отново / Výrobek opakovaně nesterilizujte / Må ikke gensteriliseres / Μην επαναποστειρώνετε / Mitte resterilisierida / Ei saa steriloida uudeleen / Ne sterilizirajte ponovno / Ne sterilizálja újra / Nesterilizuokite pakartotinai / Nesterilizēt atkārtoti / Niet opnieuw steriliseren / Må ikke steriliseres på nytt / Nie sterylizować ponownie / A nu se reesteriliza / Nesterilizujte znovu / Ne sterilizirajte ponovno / Nemojte ponovo sterilisati / Får ej omsteriliseras / Tekrar sterilize etmeyin / Не стерилізувати повторно
REF	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Número do catálogo / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Katalooginumber / Luettelonumero / Kataloški broj / Katalógusszám / Katalogo numeris / Kataloga numurs / Catalogusnummer / Katalognummer / Numer katalogowy / Număr catalog / Katalogové číslo / Številka kataloga / Kataloški broj / Katalognummer / Katalog numarası / Номер за каталогом
	Temperature limitation / Limiti di temperatura / Límites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Температурно ограничение / Teplotní omezení / Temperaturbegrænsning / Περιορισμός θερμοκρασίας / Temperaturuppirang / Lämpötilaraja / Temperaturu ograničenje / Hörmersäktet-korlátozás / Temperatūros apribojimas / Temperatūras ierobežojums / Temperaturuulimiet / Temperaturbegrensning / Ograniczenie temperatury / Limite de temperatură / Teplotné limity / Temperatura omejetev / Temperaturu ograničenje / Temperaturbegrænsning / Sicaklık kisitlamaları / Температурні обмеження
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Използвайте преди / Spotřebujte do / Udløbsdato / Χρήση έως / Kasuta kuni kuupäevani / Viimeinen käyttöpäivä / Upotrijebite do / Felhasználhatósági idő / Tinka naudoti iki / Izlietot līdz / Te gebruiken voor / Utløpsdato / Užyc do / A se utiliza înainte de / Použit do / Rok uporabe / Upotrebiti do / Används före / Son kullanim tarihi / Використати до
	Peel / Strappare per aprire / Desprender / Abziehen / Décoller / Destacar / Отлепете тук / Odleupněte / Riv op / Avoyõõvno / Eemaldamine / Kuori / Oguliti / Húzza le / Atpleşti / Nolobit / Open trekken / Rides av / Peel / Desprindeți / Odtrhnúť / Sterilizacijska vrečka / Odlepti / Riv av / Soylarak Açın / Відкрити тут
LOT	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Код на партидата (Серия) / Kód dávky (šarže) / Batchcode (lot) / Кодикós партидас (Партида) / Partii tähis (Lot) / Eräkoodi (erä) / Broj serije (Lot) / Tételkód (tétel) / Partijos kodas (partija) / Partijas kods (partija) / Batchcode (partij) / Batchcode (lot) / Kod partii (Seña) / Cod lot (lot) / Kód dávky (šarža) / Serijasa oznaka (šarža) / Šifra serije / Satskod (Lot) / Parti kodu (Lot) / Код партии (партия)
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> test / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contém o suficiente para <n> testes / Достаточна кількість, достатньо для <n> тестів / Obsah postačuje pro <n> testů / Indhold tilstrækkeligt til <n> prøver / Περιέχει επαρκή ποσότητα για <n> εξετάσεις / Sisust piisab <n> analüüsi tegemiseks / Sisältää riittävästi <n> testiin / Sadržaj dovoljno za <n> ispitivanja / <n> testhez elegendő anyagot tartalmaz / Pakanka <n> tyrimų / Satura pietiek <n> testiem / Bevat voldoende voor <n> testen / Inneholder tilstrækkelig til <n> tester / Zawiera materiał pozwalający na wykonanie <n> testów / Conținut suficient pentru <n> teste / Obsahuje dostatočné množstvo pre <n> testov / Vsebuje dovolj za <n> testov / Sadržaj dovoljno za <n> testova / Inneholder tilrækkelig mykdet for <n> tester / <n> test için yeterli içerik / Придатний для проведення <n> кількості випробувань

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Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:
Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562, USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com

Copan

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Copan Italia S.p.A.