MSwab[®]

Instructions for Use

English	3
Italiano	13
Español	20
Deutsch	26
Français	33
Português	40
БЪЛГАРСКИ	46
Česky	53
Dansk	59
ΕΛΛΗΝΙΚΑ	66
Eesti Keel	73
Hrvatski	79
Latviešu	85
Lietuvių K.	91
Magyar	97
Nederlands	103
Norsk	110
Polski	116
Română	122
Slovenčina	128
Slovenščina	134
Suomi	140
Svenska	146
Türkçe	152

Copan MSwab® Collection, Preservation and Transport System

Instruction for Use

INTENDED USE

The MSwab® system is used for the collection, transport and preservation of clinical specimens containing Gram positive aerobic and facultative anaerobic bacteria, HSV 1 and HSV 2 from the collection site to the testing laboratory. In the laboratory, MSwab® specimens are processed using standard clinical laboratory operating procedures for culture.

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of bacteriological infections involves the collection and safe transportation of swab samples. This can be accomplished by using the Copan MSwab® collection, transport and preservation system. Copan MSwab® incorporates a transport and preservation medium containing Organic solvent, Buffer, Bovine Serum Albumin and Distilled water. The medium is designed to maintain the viability of Gram positive aerobic and facultative anaerobic bacteria and HSV 1 and HSV 2 during transit to the testing laboratory.

Copan MSwab® Collection, Transport and Preservation System is supplied in collection kit format. Each collection kit consists of a package containing a plastic screw-cap tube filled with 1.6 mL of MSwab® transport and preservation medium and a small sterile peel pouch containing one specimen collection swab that has a tip flocced with soft nylon fiber.

Once a swab sample is collected, it should be placed immediately into the MSwab® transport tube where it comes into contact with the transport medium. Swab specimens for bacterial or viral investigations collected using MSwab® should be transported directly to the laboratory, preferably within 2 hours of collection ^(1, 2, 7) in order to maintain optimum organism viability. If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 hours. Bacterial viability studies for *Staphylococcus aureus*, ATCC® 29213 and ATCC® 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC® 43300 and ATCC® 700698 show viability of the tested organisms up to 14 days at refrigerated temperature (4 – 8°C) or 72 hours at room temperature (20 – 25°C). Independent scientific studies on swab transport systems have shown that for certain bacteria viability is superior at refrigerated temperatures compared with room temperature ^(12–21). If viral samples must be frozen, they should be frozen at -70 °C.

REAGENTS

Formulation of MSwab® Transport Medium

- Organic solvent
- Buffer
- Bovine Serum Albumin
- Distilled water

pH: 8.5 ± 0.20

PRODUCT DESCRIPTION

Copan MSwab® Collection, Transport and Preservation System is supplied the in product configurations indicated in the table below.

Catalog No.	Copan MSwab® Product Descriptions	Pack Size	Capture Cap Feature
404C 404C.R	Single use sample collection pack containing: - Polypropylene screw-cap tube with internal conical shape filled with 1.6mL of MSwab® Medium. - One regular size applicator swab with flocced nylon fiber tip and breaking point, sterile and individually wrapped.	50 units per shelf pack 6x50 units per box	YES

Copan MSwab® Collection, Transport and Preservation System is supplied in Collection kit format. The collection kit consists of a package containing a tube filled with MSwab® medium and a small sterile peel pouch containing one regular swab that has a tip flocced with nylon fiber intended for the collection of specimens from anatomical sites such as throat, vagina, wounds, rectum and faeces. The swab has a breakpoint in the shaft which is highlighted with a colored indication line marked. After the sample is collected from the patient, the breakpoint facilitates easy breakage of the swab applicator into the tube. The tube of both formats has a plastic screw capture cap and a conical shaped bottom filled with MSwab® medium.

The MSwab® tube capture caps have an internal moulded design that is able to capture the swab shaft when it is broken off into the tube and the cap is closed. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a moulded docking receptacle in the cap (Fig.1). In the testing laboratory when the cap is unscrewed and removed, the swab applicator is attached to the cap. This feature allows the operator to conveniently remove the swab from the transport tube.

Fig 1. Capture of broken swab applicator stick by MSwab® tube cap

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for isolating and culturing aerobic and facultative anaerobic bacteria. These materials include culture media plates or tubes and incubation systems. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques for aerobic and facultative anaerobic from clinical swab samples^(2, 4).

Appropriate materials for isolating, differentiating and culturing viruses. These materials include tissue culture cell lines, tissue culture medium, incubation systems and reading equipment. Refer to appropriate references for recommended protocols for isolation and identification of viruses^(1, 7).

PRODUCT STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 5 – 25 °C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual collection unit and the specimen transport tube label.

SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION

Specimens collected for microbiological investigations which comprise the isolation of bacteria or viruses should be collected and handled following published manuals and guidelines^(7, 8, 4).

To maintain optimum organism viability, transport specimens collected using MSwab® directly to the laboratory, preferably within 2 hours of collection^(1, 2, 7). If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8 °C or stored at room temperature (20 – 25 °C) and processed within 48 hours. Bacterial viability studies for *Staphylococcus aureus* ATCC® 29213 and ATCC® 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC® 43300 and ATCC® 700698 show viability of the tested organisms up to 14 days at refrigerated temperature (4 – 8 °C) or 72 hours at room temperature (20 – 25 °C). If viral samples must be frozen, they should be frozen at -70 °C.

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations^(34, 35, 36, 37). Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

LIMITATIONS

1. Condition, timing and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection^(7, 8, 4).
2. MSwab® is intended for use as a collection and transport medium for Gram positive aerobic and facultative anaerobic bacteria, HSV 1 and HSV 2 viruses. MSwab® cannot be used as enrichment, selective or differential medium.
3. Not suitable to collect and transport fastidious or anaerobic bacteria.
4. MSwab® is an antibiotics free medium. Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.
5. Performance testing with Copan MSwab® was conducted using laboratory strains spiked onto a swab following the test protocols described in Clinical and Laboratory Standards Institute M40-A2 Approved Standard⁽⁹⁾. Performance testing was not conducted using human specimens.
6. Performance testing with Copan MSwab® was conducted using Copan flocked swabs.

WARNINGS

1. For in vitro diagnostics use.
2. Do not re-sterilize unused swabs.
3. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
4. Do not re-pack.
5. Not suitable for any other application than intended use.
6. The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
7. Do not use if the swab is visibly damaged (i.e. if the swab tip or swab shaft is broken).
8. Do not use the same test tube for more than one patient. This will result in incorrect diagnosis.
9. Do not bend or shape the swab before the collection of the specimen. Do not use excessive force or pressure when collecting swab samples from patients as this may result in breakage of the swab shaft.
10. Do not ingest the medium.
11. To be handled by trained personnel only.
12. It must be assumed that all specimens contain infectious micro-organisms; therefore all specimens must be handled with appropriate precautions against biological risk and aseptic techniques should be used. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste. Observe CDC Biosafety Level 2 recommendations^(31, 32, 33, 34).
13. Copan MSwab® should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed, (4) the swab package is open, or (5) there are other signs of deterioration.
14. Do not use the MSwab® medium for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
15. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.
16. Repeated freezing and thawing of specimens may reduce the recovery of viable organisms^(8, 35).

INSTRUCTION FOR USE

Specimen Collection


Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals^(7, 2).

For MSwab® codes 404C and 404C.R:

1. Peel open the kit package and remove the tube of medium and inner pouch containing the sterile swab applicator (see Figure 2).

- Remove the swab applicator from its peel pouch (see Figure 2) and collect the clinical specimen. The operator must touch the swab applicator only above the colored breakpoint line, as illustrated in Figure 3, which is the opposite end to the nylon fiber tip. At all times when handling the swab applicator, the operator must not touch the area below the colored breakpoint line (the area from the line to the tip of the nylon flocked swab) as this will lead to contamination of the applicator shaft and the subsequent culture.
- Collect the sample from the patient.
- Unscrew and remove the cap from MSwab® tube making sure not to spill the medium.
- After the swab sample is taken from the patient, insert the swab into the tube until the breakpoint is level with the test tube opening.
- Bend the swab shaft at a 180 degrees angle to break it off at the breaking point. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
- Discard the broken off part of the swab shaft in an approved medical waste disposal container.
- Screw the cap and secure tightly onto the tube (see Figure 2).
- Write patient's name and data on the tube label.
- Send the sample to the laboratory.


Fig 2. Collection swab showing breakpoint indication line and area for holding the applicator

Sterile gloves and protective clothing and eyewear should be worn when collecting and handling microbiology specimens and care should be taken to avoid splashes and aerosols when breaking the swab stick into the tube of medium. During sample collection when handling the swab applicator, the operator must not touch the area below the colored breakpoint indication line; that is the area from the line to the tip of the nylon flocked swab (see Fig. 3), as this will lead to contamination of the applicator shaft and the culture thus invalidating the test results.

Fig. 3 Collection swab showing breakpoint indication line and area for holding the applicator

Do not touch the applicator in the area below the breakpoint indication line.
During specimen collection hold applicator above the breakpoint indication line, in this area

Molded breakpoint with colored indication line
The operator must only handle the part of the swab applicator shaft above the breakpoint indication line.

Processing MSwab® Specimens in the Laboratory – Bacteriology

MSwab® samples should be processed for bacteriological culture using recommended culture media and laboratory techniques which will depend on the specimen type and the organism under investigation. For recommended culture media and techniques for the isolation and identification of bacteria from clinical swab specimens refer to published microbiology manuals and guidelines⁽¹⁻⁶⁾.

Culture investigations of swab specimens for the presence of bacteria routinely involve the use of solid agar culture medium in Petri dish plates. The procedure for inoculation of MSwab® samples onto solid agar in Petri dishes is as follows.

Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations^(31, 32, 33, 34).

Vortex mix the MSwab® tube containing the swab sample for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the medium

- Unscrew the MSwab® cap and remove the swab applicator.
- Roll the tip of the MSwab® applicator onto the surface of one quadrant of the culture media plate to provide the primary inoculum.
- If it is necessary to culture the swab specimen onto a second culture media plate, return the MSwab® applicator to the transport medium tube for two seconds to absorb and recharge the applicator tip with transport medium/patient sample suspension then repeat Step No. 3.
- If it is necessary to inoculate additional culture media plates, return the MSwab® applicator to the transport medium tube and recharge the swab applicator tip with the transport medium/patient sample suspension before inoculating each additional plate.

The procedure described above utilizes the MSwab® applicator like an inoculation wand to transfer the suspension of patient sample in transport medium onto the surface of a culture plate creating the primary inoculum (see Fig 4).


Alternatively, the operator can vortex mix the MSwab® tube with the swab inside for 5 seconds and then transfer 100µl volumes of the suspension onto each culture plate using a volumetric pipetor and sterile pipet tips. Standard laboratory techniques should then be used to streak the primary inoculum of patient sample across the surface of the culture plate (see Fig 5).

Fig 4. Procedures for inoculation of MSwab® specimens onto solid agar in Petri dishes

1. Using swab to inoculate specimen
2. Using pipettor and sterile pipet tips to inoculate 100µl of specimen

Fig 5. Procedure for streaking MSwab® specimens on agar Petri dishes for primary isolation ⁽³³⁾

Seed a primary inoculum of MSwab® specimen onto the surface of an appropriate agar culture plate in the first quadrant.

Use a sterile bacteriology loop to streak the primary inoculum across the surface of the second, third and fourth quadrants of the agar culture plate.

Preparation of Gram Stain Smears of MSwab® Specimens

Laboratory analysis of clinical swab samples collected from certain sites on the patient can routinely include microscopic examination of stained preparations ("direct smears") using the Gram stain procedure. This can provide valuable information to physicians who are managing patients with infectious diseases ⁽²²⁾. There are many instances in which a Gram stain can assist in making a diagnosis ^(23, 27).

The Gram stain can also help to judge specimen quality and contribute to the selection of culture media especially with mixed flora.

Microscope slides of patient specimens transported in Copan MSwab® transport system can be prepared for Gram stain analysis, as describe below, by sampling an aliquot of vortexed suspension of the swab ^(3, 4). Sample transported in MSwab® elution medium represent a homogeneous suspension in liquid phase. It can be uniformly smeared allowing clear and easy reading.

Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations ^(31, 32, 33, 34).

1. Take a clean glass microscope slide, place it on a flat surface and inscribe an area using a diamond-tipped or similar glass marker to identify the location of the specimen inoculum. Note: a slide with a pre-marked 20 mm well can be used.
2. Vortex mix the MSwab® tube containing the swab sample for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the medium.
3. Unscrew the MSwab® cap and using a sterile pipet, transfer 1 – 2 drops of sample suspension to the inscribed area on the glass slide. Note: about 30µl would be a suitable amount of liquid for a pre-marked 20 mm diameter well slide.

In case of bloody or thicker specimens particular care should be taken to thinly spread the sample on the slide. Bacteria are difficult to detect if the sample shows many red cells and debris.

4. Allow the specimen on the slide to air dry at room temperature or place the slide in an electric slide warmer or incubator set at a temperature not exceeding 42 °C.
5. Fix smears using methanol. Methanol fixation is recommended as it prevents lysis of Red Blood Cells, avoids damage to all host cells and results in a cleaner background ^(3, 4, 22).
6. Follow published laboratory reference manuals and guidelines for performing the Gram stain. If commercial Gram stain reagents are used, it is important to comply with instructions in the manufacturer's product insert for performance test procedure.

For further information or guidance on the preparation of specimen slides for microscopic analysis, for information on Gram staining procedures and the interpretation and reporting of microscopic analysis, consult published laboratory reference manuals ^(1 - 5, 22 - 27).

Processing MSwab® Specimens in the Laboratory – Virology

The survival of HSV 1 and HSV 2 depends on many factors including the type and concentration of the microorganism, duration of transport and storage temperature. To maintain optimum viability, specimens should be transported directly to the laboratory preferably within 2 hours of collection ^(1, 2, 7, 29). If immediate delivery or processing is delayed, then specimens collected using Copan MSwab® Collection, Transport and Preservation System should be refrigerated at 4-8 °C or stored at room temperature (20-25°C) and processed within 48 hours. If samples must be frozen, they should be frozen placed at -70°C.

In simulated transport and storage studies, the Copan MSwab® System was shown to be capable of maintaining the viability of HSV 1 and HSV 2 at refrigerated (4-8°C) and room temperature (20-25°C) conditions for up to 48 hours. Based on performance studies conducted by Copan and independent scientific publications, viability of certain microorganisms is superior at refrigerated temperatures compared with room temperature ^(12 - 21, 29).

MSwab® specimens should be processed for virology culture using recommended cells line and laboratory techniques that will depend on the specimen type and the organism under investigation. For recommended shell vials and techniques for the isolation and identification of HSV 1 and HSV 2 from clinical swab specimens, refer to published virology manuals and guidelines ^(1, 6, 29, 30).

Culture investigations of swab specimens for the presence of HSV 1 and HSV 2 routinely involve the use of cells culture in shell vials. The procedure for inoculation of MSwab® specimens onto shell vials is described below.

1. Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other BSL 2 recommendations.
 2. Shake using a vortex mixer for 5 seconds the MSwab® tube containing the swab specimen to release the specimen from the swab tip and evenly disperse and suspend the patient specimen in the liquid transport medium.
 3. Unscrew the MSwab® cap and remove the swab applicator.
 4. Transfer 200 µl volumes of the suspension into the shell vial and proceed as per laboratory internal procedure.
- NOTE: Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.
5. Proceed with the appropriate techniques for virus detection.

QUALITY CONTROL

MSwab® applicators are tested to ensure they are non-toxic to bacteria. MSwab® medium and applicators are tested to ensure they are non-toxic to cell lines used for HSV 1 and HSV 2 cultivation. MSwab® transport medium is tested for pH stability⁽⁹⁾. MSwab® is quality control tested before release for ability to maintain viable Gram positive aerobic and facultative anaerobic bacteria and HSV viruses at room temperature (20 – 25°C) for specified time points. Procedures for quality control of microbiology transport devices should be conducted using testing methods described in Clinical and Laboratory Standards Institute M40-A2 and other publications⁽⁹⁾. If aberrant quality control results are noted, patient results should not be reported.

RESULTS

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

PERFORMANCE CHARACTERISTICS

The test procedures employed for determining bacterial viability performance were based upon the quality control methods described in Clinical and Laboratory Standards Institute M40-A2⁽⁹⁾.

MSwab® system has an intended use limited to Gram positive aerobic and facultative anaerobic bacteria and HSV 1 and HSV 2, therefore its field application is more restricted than some other devices. For this reason the bacterial recovery studies were conducted under the simulated transport and storage conditions as described and defined in CLSI M40-A2, Quality Control of Microbiological Transport Systems: Approved Standard and included only the Gram positive aerobic and facultative anaerobic strains from the Group 1 of paragraph 7.11.1 of the CLSI M40-A2 document, in particular:

<i>Streptococcus pyogenes</i>	ATCC® 19615
<i>Streptococcus pneumoniae</i>	ATCC® 6305
<i>Pseudomonas aeruginosa</i>	ATCC® BAA-427

In addition, Copan included testing of additional Gram positive aerobic and facultative anaerobic microorganisms, clinically relevant, not required by CLSI M40-A2. The specific bacterial strains used in these studies are here listed:

<i>Enterococcus faecalis</i>	ATCC® 29212
<i>Staphylococcus epidermidis</i>	ATCC® 12228
<i>Staphylococcus aureus</i>	ATCC® 29213
<i>Streptococcus agalactiae</i> (Group B Strep)	ATCC® 13813
<i>Listeria monocytogenes</i>	ATCC® 19114
<i>Bacillus cereus</i>	ATCC® 10876
<i>Staphylococcus aureus</i> (Methicillin resistant)	ATCC® 43300
<i>Staphylococcus aureus</i>	ATCC® 6538
<i>Staphylococcus aureus</i> (Methicillin resistant)	ATCC® 700698

All cultures of bacteria were ATCC® (American Type Culture Collection) and were obtained commercially.

The selection of these organisms also reflects those Gram positive aerobic and facultative anaerobic bacteria that would normally be encountered in specimens collected and analyzed in a typical clinical microbiology laboratory.

Bacterial viability studies were performed on the Copan MSwab® at two different ranges of temperature, 4 – 8 °C and 20 – 25 °C, corresponding to refrigerator and room temperature, respectively. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs, 24 hrs and 48 hrs. At the appropriate time intervals, each swab was processed according to the Roll-Plate or Swab Elution Method.

Additional bacterial viability studies for *Staphylococcus aureus*, ATCC® 29213 and ATCC® 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC® 43300 and ATCC® 700698 were performed on the Copan MSwab® at two different ranges of temperature, 4 – 8 °C and 20 – 25 °C, corresponding to refrigerator and room temperature, respectively.

Swab accompanying the transport system was inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and:

For studies performed at 4 – 8 °C inoculated MSwab® tubes were held for 0 hrs, 10 days and 14 days. At the appropriate time intervals, each MSwab® was processed according to the Roll-Plate Method.

For studies performed at 20 – 25 °C inoculated MSwab® tubes were held for 0 hrs and 72 hrs. At the appropriate time intervals, each MSwab® was processed according to the Roll-Plate Method.

Bacterial overgrowth studies were performed on the Copan MSwab® at 4 – 8 °C, corresponding to refrigerator temperature. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs and 48 hrs. At the appropriate time intervals, each swab was processed according to the Roll-Plate Method.

Bacterial overgrowth studies were performed using *Pseudomonas aeruginosa*.

Viral viability studies were performed using HSV 1 and HSV 2. Swabs accompanying each transport system were directly inoculated in triplicate with 100µl of organisms suspension.

Swabs were then placed in their respective transport medium tubes and were held for 0, 24 and 48 hours at both 4 °C and room temperature (20-25 °C). At the appropriate time interval, each swab was vortexed, removed from its transport medium tube and then 200µl aliquots of this suspension was inoculated into shell vials. All cultures were processed by standard laboratory culture technique and examined after a specified incubation time. Organism viability was determined by fluorescing foci counts.

Organisms evaluated were:

Herpes Simplex Virus Type 1 (HSV 1) ATCC® VR-539

Herpes Simplex Virus Type 2 (HSV 2) ATCC® VR-734.

TEST RESULTS

***Table 1 - SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES SWAB ELUTION METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 24 hrs	Average of CFUs recovered at time 48 hrs	Log ₁₀ decline	Interpretation
<i>Streptococcus pyogenes</i> ATCC® 19615	diluted 1:10	MSwab®	2045	3.9 x 10 ⁵	2.7 x 10 ⁵	2.3 x 10 ⁵	-0.23	Acceptable Recovery
			2045/1	4.8 x 10 ⁵	2.5 x 10 ⁵	2.1 x 10 ⁵	-0.36	Acceptable Recovery
			2045/2	4.2 x 10 ⁵	2.3 x 10 ⁵	2.0 x 10 ⁵	-0.32	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC® 6305	diluted 1:10	MSwab®	2045	1.2 x 10 ⁵	1.8 x 10 ⁴	2.1 x 10 ⁴	-1.76	Acceptable Recovery
			2045/1	1.2 x 10 ⁵	1.9 x 10 ⁴	1.7 x 10 ⁴	-1.85	Acceptable Recovery
			2045/2	1.2 x 10 ⁵	2.0 x 10 ⁴	1.9 x 10 ⁴	-1.80	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC® 49136	diluted 1:10	MSwab®	2045	2.0 x 10 ⁵	2.0 x 10 ⁵	2.1 x 10 ⁵	-1.98	Acceptable Recovery
			2045/1	1.8 x 10 ⁵	2.1 x 10 ⁴	1.8 x 10 ⁴	-2.00	Acceptable Recovery
			2045/2	1.9 x 10 ⁵	1.8 x 10 ⁴	1.7 x 10 ⁴	-2.05	Acceptable Recovery
<i>Enterococcus faecalis</i> ATCC® 29212	diluted 1:10	MSwab®	2045	1.1 x 10 ⁵	1.0 x 10 ⁵	9.4 x 10 ⁴	-0.07	Acceptable Recovery
			2045/1	1.1 x 10 ⁵	1.0 x 10 ⁵	9.2 x 10 ⁴	-0.08	Acceptable Recovery
			2045/2	1.1 x 10 ⁵	1.0 x 10 ⁵	9.5 x 10 ⁴	-0.06	Acceptable Recovery
<i>Staphylococcus epidermidis</i> ATCC® 12228	diluted 1:10	MSwab®	2045	1.6 x 10 ⁵	1.2 x 10 ⁵	9.2 x 10 ⁴	-0.24	Acceptable Recovery
			2045/1	1.7 x 10 ⁵	1.3 x 10 ⁵	1.1 x 10 ⁵	-0.19	Acceptable Recovery
			2045/2	1.4 x 10 ⁵	1.2 x 10 ⁵	9.5 x 10 ⁴	-0.17	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 29213	diluted 1:10	MSwab®	2045	1.3 x 10 ⁵	1.0 x 10 ⁵	8.8 x 10 ⁴	-0.17	Acceptable Recovery
			2045/1	1.0 x 10 ⁵	9.5 x 10 ⁴	7.4 x 10 ⁴	-0.13	Acceptable Recovery
			2045/2	1.0 x 10 ⁵	9.5 x 10 ⁴	8.0 x 10 ⁴	-0.10	Acceptable Recovery
<i>Streptococcus agalactiae</i> ATCC® 13813	diluted 1:10	MSwab®	2045	2.1 x 10 ⁵	1.8 x 10 ⁵	1.5 x 10 ⁵	-0.15	Acceptable Recovery
			2045/1	1.9 x 10 ⁵	1.8 x 10 ⁵	1.5 x 10 ⁵	-0.10	Acceptable Recovery
			2045/2	2.2 x 10 ⁵	1.4 x 10 ⁵	1.3 x 10 ⁵	-0.23	Acceptable Recovery
<i>Kocuria rhizophila</i> ATCC® 9341	diluted 1:10	MSwab®	2045	2.0 x 10 ⁵	1.5 x 10 ⁵	1.4 x 10 ⁵	-0.15	Acceptable Recovery
			2045/1	2.1 x 10 ⁵	1.4 x 10 ⁵	1.4 x 10 ⁵	-0.18	Acceptable Recovery
			2045/2	2.2 x 10 ⁵	1.6 x 10 ⁵	1.4 x 10 ⁵	-0.19	Acceptable Recovery
<i>Listeria monocytogenes</i> ATCC® 19114	diluted 1:10	MSwab®	2045	1.4 x 10 ⁵	8.1 x 10 ⁴	3.0 x 10 ⁴	-0.67	Acceptable Recovery
			2045/1	1.3 x 10 ⁵	9.0 x 10 ⁴	3.4 x 10 ⁴	-0.58	Acceptable Recovery
			2045/2	1.5 x 10 ⁵	8.8 x 10 ⁴	4.2 x 10 ⁴	-0.55	Acceptable Recovery
<i>Bacillus cereus</i> ATCC® 10876	diluted 1:10	MSwab®	2045	1.5 x 10 ⁵	6.9 x 10 ⁴	1.6 x 10 ⁴	-0.97	Acceptable Recovery
			2045/1	1.6 x 10 ⁵	7.6 x 10 ⁴	1.8 x 10 ⁴	-0.95	Acceptable Recovery
			2045/2	1.6 x 10 ⁵	8.3 x 10 ⁴	1.4 x 10 ⁴	-1.06	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 43300	diluted 1:10	MSwab®	2045	2.1 x 10 ⁵	1.2 x 10 ⁵	8.0 x 10 ⁴	-0.42	Acceptable Recovery
			2045/1	1.9 x 10 ⁵	1.3 x 10 ⁵	8.9 x 10 ⁴	-0.33	Acceptable Recovery
			2045/2	1.9 x 10 ⁵	1.4 x 10 ⁵	9.4 x 10 ⁴	-0.30	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 6538	diluted 1:10	MSwab®	2045	1.5 x 10 ⁵	9.2 x 10 ⁴	6.9 x 10 ⁴	-0.34	Acceptable Recovery
			2045/1	1.6 x 10 ⁵	1.0 x 10 ⁵	8.3 x 10 ⁴	-0.29	Acceptable Recovery
			2045/2	1.7 x 10 ⁵	1.3 x 10 ⁵	9.5 x 10 ⁴	-0.25	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 700698	diluted 1:10	MSwab®	2045	2.0 x 10 ⁵	1.6 x 10 ⁵	1.4 x 10 ⁵	-0.15	Acceptable Recovery
			2045/1	2.1 X 10 ⁵	1.6 x 10 ⁵	1.4 x 10 ⁵	-0.18	Acceptable Recovery
			2045/2	1.9 X 10 ⁵	1.5 x 10 ⁵	9.5 x 10 ⁴	-0.30	Acceptable Recovery

***Table 2 - SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES SWAB ELUTION METHOD, 20-25°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 24 hrs	Average of CFUs recovered at time 48 hrs	Log ₁₀ decline	Interpretation
<i>Streptococcus pyogenes</i> ATCC® 19615	diluted 1:10	MSwab®	2045	3.9 x 10 ⁵	2.2 x 10 ⁵	1.5 x 10 ⁵	-0.41	Acceptable Recovery
			2045/1	4.8 x 10 ⁵	2.1 x 10 ⁵	1.5 x 10 ⁵	-0.51	Acceptable Recovery
			2045/2	4.2 x 10 ⁵	2.0 x 10 ⁵	1.7 x 10 ⁵	-0.39	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC® 6305	diluted 1:10	MSwab®	2045	1.2 x 10 ⁵	1.0 x 10 ⁴	9.4 x 10 ²	-2.11	Acceptable Recovery
			2045/1	1.2 x 10 ⁵	1.2 x 10 ⁴	1.1 x 10 ³	-2.04	Acceptable Recovery
			2045/2	1.2 x 10 ⁵	1.2 x 10 ⁴	1.1 x 10 ³	-2.04	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC® 49136	diluted 1:10	MSwab®	2045	2.0 x 10 ⁵	1.5 x 10 ⁴	1.6 x 10 ³	-2.10	Acceptable Recovery
			2045/1	1.8 x 10 ⁵	1.5 x 10 ⁴	1.1 x 10 ³	-2.21	Acceptable Recovery
			2045/2	1.9 x 10 ⁵	1.7 x 10 ⁴	1.4 x 10 ³	-2.13	Acceptable Recovery
<i>Enterococcus faecalis</i> ATCC® 29212	diluted 1:10	MSwab®	2045	1.1 x 10 ⁶	1.1 x 10 ⁶	7.7 x 10 ⁵	-0.15	Acceptable Recovery
			2045/1	1.1 x 10 ⁶	9.6 x 10 ⁵	8.1 x 10 ⁵	-0.13	Acceptable Recovery
			2045/2	1.1 x 10 ⁶	9.3 x 10 ⁵	8.6 x 10 ⁵	-0.11	Acceptable Recovery
<i>Staphylococcus epidermidis</i> ATCC® 12228	diluted 1:10	MSwab®	2045	1.6 x 10 ⁶	1.1 x 10 ⁶	5.2 x 10 ⁵	-0.49	Acceptable Recovery
			2045/1	1.7 x 10 ⁶	1.1 x 10 ⁶	6.8 x 10 ⁵	-0.40	Acceptable Recovery
			2045/2	1.4 x 10 ⁶	9.0 x 10 ⁵	5.5 x 10 ⁵	-0.41	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 29213	diluted 1:10	MSwab®	2045	1.3 x 10 ⁶	8.7 x 10 ⁵	5.2 x 10 ⁵	-0.40	Acceptable Recovery
			2045/1	1.0 x 10 ⁶	8.7 x 10 ⁵	5.8 x 10 ⁵	-0.24	Acceptable Recovery
			2045/2	1.0 x 10 ⁶	9.9 x 10 ⁵	6.0 x 10 ⁵	-0.22	Acceptable Recovery
<i>Streptococcus agalactiae</i> ATCC® 13813	diluted 1:10	MSwab®	2045	2.1 x 10 ⁵	1.5 x 10 ⁵	1.3 x 10 ⁵	-0.21	Acceptable Recovery
			2045/1	1.9 x 10 ⁶	1.7 x 10 ⁶	1.3 x 10 ⁵	-0.16	Acceptable Recovery
			2045/2	2.2 x 10 ⁶	1.4 x 10 ⁶	1.0 x 10 ⁶	-0.34	Acceptable Recovery
<i>Kocuria rhizophila</i> ATCC® 9341	diluted 1:10	MSwab®	2045	2.0 x 10 ⁵	1.1 x 10 ⁵	6.0 x 10 ⁴	-0.52	Acceptable Recovery
			2045/1	2.1 x 10 ⁵	1.5 x 10 ⁵	5.4 x 10 ⁴	-0.59	Acceptable Recovery
			2045/2	2.2 x 10 ⁵	1.3 x 10 ⁵	6.5 x 10 ⁴	-0.53	Acceptable Recovery
<i>Listeria monocytogenes</i> ATCC® 19114	Diluted 1:10	MSwab®	2045	1.4 x 10 ⁵	3.8 x 10 ⁵	2.0 x 10 ⁵	-0.84	Acceptable Recovery
			2045/1	1.3 x 10 ⁶	5.3 x 10 ⁵	1.8 x 10 ⁵	-0.86	Acceptable Recovery
			2045/2	1.5 x 10 ⁶	5.3 x 10 ⁵	1.7 x 10 ⁵	-0.94	Acceptable Recovery
<i>Bacillus cereus</i> ATCC® 10876	diluted 1:10	MSwab®	2045	1.5 x 10 ⁶	5.8 x 10 ⁴	2.0 x 10 ³	-1.87	Acceptable Recovery
			2045/1	1.6 x 10 ⁶	7.8 x 10 ⁴	2.1 x 10 ³	-1.88	Acceptable Recovery
			2045/2	1.6 x 10 ⁵	8.1 x 10 ⁴	3.8 x 10 ³	-1.62	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 43300	diluted 1:10	MSwab®	2045	2.1 x 10 ⁶	1.0 x 10 ⁶	5.0 x 10 ⁵	-0.62	Acceptable Recovery
			2045/1	1.9 x 10 ⁶	1.4 x 10 ⁶	4.2 x 10 ⁵	-0.65	Acceptable Recovery
			2045/2	1.9 x 10 ⁶	9.1 x 10 ⁵	4.3 x 10 ⁵	-0.65	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 6538	diluted 1:10	MSwab®	2045	1.5 x 10 ⁶	7.6 x 10 ⁵	5.1 x 10 ⁵	-0.47	Acceptable Recovery
			2045/1	1.6 x 10 ⁶	1.1 x 10 ⁶	6.4 x 10 ⁵	-0.40	Acceptable Recovery
			2045/2	1.7 x 10 ⁶	1.0 x 10 ⁶	7.9 x 10 ⁵	-0.33	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 700698	diluted 1:10	MSwab®	2045	2.0 x 10 ⁶	1.6 x 10 ⁶	1.1 x 10 ⁶	-0.26	Acceptable Recovery
			2045/1	2.1 X 10 ⁶	1.5 x 10 ⁶	1.3 x 10 ⁶	-0.21	Acceptable Recovery
			2045/2	1.9 X 10 ⁶	1.4 x 10 ⁶	8.2 x 10 ⁵	-0.36	Acceptable Recovery

***Table 3 - SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES ROLL-PLATE METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 24 hrs	Average of CFUs recovered at time 48 hrs	Interpretation
<i>Streptococcus pyogenes</i> ATCC® 19615	diluted 10 ⁻³	MSwab®	2045	263.3	172.3	153.0	Acceptable Recovery
			2045/1	266.7	220.0	169.0	Acceptable Recovery
			2045/2	275.7	240.7	160.3	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC® 6305	diluted 10 ^{-1.5}	MSwab®	2045	227.3	153.0	111.7	Acceptable Recovery
			2045/1	215.7	171.0	109.0	Acceptable Recovery
			2045/2	217.3	159.0	104.3	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC® 49136	diluted 10 ^{-1.5}	MSwab®	2045	279.7	206.3	123.7	Acceptable Recovery
			2045/1	256.7	207.7	113.0	Acceptable Recovery
			2045/2	270.0	195.7	123.0	Acceptable Recovery

<i>Enterococcus faecalis</i> ATCC® 29212	diluted $10^{-3.5}$	MSwab®	2045	207.3	172.3	158.3	Acceptable Recovery
			2045/1	214.0	167.7	158.3	Acceptable Recovery
			2045/2	223.3	183.3	140.0	Acceptable Recovery
<i>Staphylococcus epidermidis</i> ATCC® 12228	diluted $10^{-3.5}$	MSwab®	2045	280.7	238.7	123.0	Acceptable Recovery
			2045/1	277.0	238.0	113.0	Acceptable Recovery
			2045/2	273.3	151.3	117.0	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 29213	diluted $10^{-3.5}$	MSwab®	2045	252.7	209.3	179.0	Acceptable Recovery
			2045/1	212.7	188.3	160.0	Acceptable Recovery
			2045/2	217.0	188.7	157.3	Acceptable Recovery
<i>Streptococcus agalactiae</i> ATCC® 13813	diluted 10^{-3}	MSwab®	2045	174.0	121.0	117.7	Acceptable Recovery
			2045/1	162.3	127.0	100.3	Acceptable Recovery
			2045/2	186.0	148.3	124.7	Acceptable Recovery
<i>Kocuria rhizophila</i> ATCC® 9341	diluted 10^{-2}	MSwab®	2045	271.7	237.0	212.0	Acceptable Recovery
			2045/1	267.7	223.3	202.0	Acceptable Recovery
			2045/2	269.3	235.0	179.0	Acceptable Recovery
<i>Listeria monocytogenes</i> ATCC® 19114	diluted 10^{-3}	MSwab®	2045	267.3	243.3	209.7	Acceptable Recovery
			2045/1	275.0	242.0	192.7	Acceptable Recovery
			2045/2	274.3	226.7	198.3	Acceptable Recovery
<i>Bacillus cereus</i> ATCC® 10876	diluted $10^{-1.5}$	MSwab®	2045	234.7	185.3	167.0	Acceptable Recovery
			2045/1	246.0	182.7	151.7	Acceptable Recovery
			2045/2	231.0	179.0	144.7	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 43300	diluted $10^{-3.5}$	MSwab®	2045	204.7	185.0	172.7	Acceptable Recovery
			2045/1	211.7	185.7	177.3	Acceptable Recovery
			2045/2	196.7	183.7	161.0	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 6538	diluted $10^{-3.5}$	MSwab®	2045	149.0	123.0	98.7	Acceptable Recovery
			2045/1	172.0	138.0	108.3	Acceptable Recovery
			2045/2	167.3	142.3	98.3	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 700698	diluted $10^{-3.5}$	MSwab®	2045	257.0	223.3	204.0	Acceptable Recovery
			2045/1	254.3	229.0	203.3	Acceptable Recovery
			2045/2	257.7	238.0	202.7	Acceptable Recovery

***Table 4 - SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES ROLL-PLATE METHOD, 20-25°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 24 hrs	Average of CFUs recovered at time 48 hrs	Interpretation
<i>Streptococcus pyogenes</i> ATCC® 19615	diluted 10^{-3}	MSwab®	2045	283.3	184.0	113.3	Acceptable Recovery
			2045/1	266.7	196.0	123.0	Acceptable Recovery
			2045/2	275.7	211.0	130.0	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC® 6305	diluted $10^{-1.5}$	MSwab®	2045	227.3	131.7	83.0	Acceptable Recovery
			2045/1	215.7	149.3	73.0	Acceptable Recovery
			2045/2	217.3	137.3	80.0	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC® 49136	diluted $10^{-1.5}$	MSwab®	2045	279.7	148.3	71.0	Acceptable Recovery
			2045/1	256.7	171.0	81.7	Acceptable Recovery
			2045/2	270.0	171.0	71.0	Acceptable Recovery
<i>Enterococcus faecalis</i> ATCC® 29212	diluted $10^{-3.5}$	MSwab®	2045	207.3	140.0	117.7	Acceptable Recovery
			2045/1	214.0	140.0	117.7	Acceptable Recovery
			2045/2	223.3	134.0	96.0	Acceptable Recovery
<i>Staphylococcus epidermidis</i> ATCC® 12228	diluted $10^{-3.5}$	MSwab®	2045	280.7	147.0	89.7	Acceptable Recovery
			2045/1	277.0	154.0	93.3	Acceptable Recovery
			2045/2	273.3	167.7	88.0	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 29213	diluted $10^{-3.5}$	MSwab®	2045	252.7	176.7	112.0	Acceptable Recovery
			2045/1	212.7	178.3	126.7	Acceptable Recovery
			2045/2	217.0	131.7	107.3	Acceptable Recovery
<i>Streptococcus agalactiae</i> ATCC® 13813	diluted 10^{-3}	MSwab®	2045	174.0	112.7	80.3	Acceptable Recovery
			2045/1	162.3	119.7	88.0	Acceptable Recovery
			2045/2	186.0	112.7	84.0	Acceptable Recovery
<i>Kocuria rhizophila</i> ATCC® 9341	diluted 10^{-2}	MSwab®	2045	271.7	198.3	165.3	Acceptable Recovery
			2045/1	267.7	196.3	145.7	Acceptable Recovery
			2045/2	269.3	204.7	142.3	Acceptable Recovery
<i>Listeria monocytogenes</i> ATCC® 19114	diluted 10^{-3}	MSwab®	2045	267.3	221.3	167.0	Acceptable Recovery
			2045/1	275.0	221.0	162.0	Acceptable Recovery
			2045/2	274.3	211.3	172.0	Acceptable Recovery

<i>Bacillus cereus</i> ATCC® 10876	diluted $10^{-1.5}$	MSwab®	2045	234.7	156.0	116.0	Acceptable Recovery
			2045/1	246.0	167.0	120.3	Acceptable Recovery
			2045/2	231.0	161.7	117.3	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 43300	diluted $10^{-3.5}$	MSwab®	2045	204.7	173.7	153.0	Acceptable Recovery
			2045/1	211.7	174.3	142.3	Acceptable Recovery
			2045/2	196.7	174.3	137.3	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 6538	diluted $10^{-3.5}$	MSwab®	2045	149.0	117.0	72.7	Acceptable Recovery
			2045/1	172.0	113.3	82.0	Acceptable Recovery
			2045/2	167.3	133.3	86.3	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 700698	diluted $10^{-3.5}$	MSwab®	2045	257.0	205.0	171.7	Acceptable Recovery
			2045/1	254.3	203.0	176.3	Acceptable Recovery
			2045/2	257.7	205.3	170.3	Acceptable Recovery

***Table 5 - SUMMARY OF RESULTS FOR ADDITIONAL BACTERIAL RECOVERY STUDIES ON SPECIFIC STRAINS ROLL-PLATE METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 10 days	Average of CFUs recovered at time 14 days	Interpretation
<i>Staphylococcus aureus</i> ATCC® 29213	diluted $10^{-3.5}$	MSwab®	1210	248	349	401	Acceptable Recovery
			25C080	232	321	497	Acceptable Recovery
			43A133	240	300	521	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 43300	diluted $10^{-3.5}$	MSwab®	1210	231	264	305	Acceptable Recovery
			25C080	220	287	289	Acceptable Recovery
			43A133	212	292	215	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 6538	diluted $10^{-3.5}$	MSwab®	1210	152	175	156	Acceptable Recovery
			25C080	88	139	160	Acceptable Recovery
			43A133	256	147	165	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 700698	diluted $10^{-3.5}$	MSwab®	1210	168	184	225	Acceptable Recovery
			25C080	176	240	209	Acceptable Recovery
			43A133	180	327	247	Acceptable Recovery

***Table 6 - SUMMARY OF RESULTS FOR ADDITIONAL BACTERIAL RECOVERY STUDIES ON SPECIFIC STRAINS ROLL-PLATE METHOD, 20-25°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 72 hrs	Interpretation
<i>Staphylococcus aureus</i> ATCC® 29213	diluted $10^{-3.5}$	MSwab®	1210	248	233	Acceptable Recovery
			25C080	232	303	Acceptable Recovery
			43A133	240	327	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 43300	diluted $10^{-3.5}$	MSwab®	1210	231	60	Acceptable Recovery
			25C080	220	80	Acceptable Recovery
			43A133	212	69	Acceptable Recovery
<i>Staphylococcus aureus</i> ATCC® 6538	diluted $10^{-3.5}$	MSwab®	1210	152	112	Acceptable Recovery
			25C080	88	194	Acceptable Recovery
			43A133	256	125	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC® 700698	diluted $10^{-3.5}$	MSwab®	1210	168	67	Acceptable Recovery
			25C080	176	84	Acceptable Recovery
			43A133	180	75	Acceptable Recovery

***Table 7 - SUMMARY OF RESULTS FOR BACTERIAL OVERGROWTH STUDIES, ROLL-PLATE METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension in saline	Product	Lot Number	Average of CFUs recovered at time 0 hrs	Average of CFUs recovered at time 48 hrs	Interpretation
<i>rfaPseudomonas aeruginosa</i> ATCC® BAA-427	diluted 10^{-4}	MSwab®	2045	39.0	0	Acceptable Recovery
			2045/1	40.3	0	Acceptable Recovery
			2045/2	36.3	0	Acceptable Recovery

***Table 8 - SUMMARY OF RESULTS FOR VIRAL RECOVERY STUDIES, 4-8°C**

Organism	Dilution	Product	Lot Number	Average of foci of infected cells at time 0 hrs	Average of foci of infected cells at time 24 hrs	Average of foci of infected cells at time 48 hrs	Log ₁₀ decline	Interpretation
HSV 1 ATCC® VR539	diluted 10 ⁻²	MSwab®	2045	535.3	302.6	175.6	-0.48	Acceptable Recovery
			2045/ 1	505.7	354.7	195.0	-0.41	Acceptable Recovery
			2045/ 2	501.7	369.3	213.7	-0.37	Acceptable Recovery
HSV 2 ATCC® VR734	diluted 10 ⁻¹	MSwab®	2045	902.0	526.0	209.7	-0.63	Acceptable Recovery
			2045/ 1	889.7	419.0	245.7	-0.56	Acceptable Recovery
			2045/ 2	954.0	486.0	275.7	-0.54	Acceptable Recovery

***Table 9 - SUMMARY OF RESULTS FOR VIRAL RECOVERY STUDIES, 20-25°C**

Organism	Dilution	Product	Lot Number	Average of foci of infected cells at time 0 hrs	Average of foci of infected cells at time 24 hrs	Average of foci of infected cells at time 48 hrs	Log ₁₀ decline	Interpretation
HSV 1 ATCC® VR539	diluted 10 ⁻²	MSwab®	2045	535.3	199.3	142.6	-0.57	Acceptable Recovery
			2045/ 1	505.7	195.7	140.7	-0.55	Acceptable Recovery
			2045/ 2	501.7	152.0	103.3	-0.69	Acceptable Recovery
HSV 2 ATCC® VR734	diluted 10 ⁻¹	MSwab®	2045	902.0	387.7	184.7	-0.69	Acceptable Recovery
			2045/ 1	889.7	414.7	203.7	-0.64	Acceptable Recovery
			2045/ 2	954.0	375.3	187.0	-0.71	Acceptable Recovery

*Translation in other languages of words/sentences used in Table 1, 2, 3, 4, 5, 6, 7, 8, 9:

Organism: Organismo / Organismo / Organismus / Organisme / Organismo / Организм / Organismus / Organisme / Οργανισμός / Organism / Tijelo / Kermenis / Organizmas / Szervezet / Organisme / Organisme / Organism / Organism / Organism / Organismi / Organismus / Organism / Organiza
Dilution 0.5 McFarland bacterial suspension with saline: Diluizione Sospensione batterica 0.5 McFarland con sol. Salina / Dilución: suspensión bacteriana 0.5 McFarland con sol. salina / Verdünnung: 0.5-McFarland-Bakteriensuspension mit Salzlösung / Dilution: Suspension bactérienne 0.5 McFarland avec sol. saline / Diluição em suspensão bacteriana McFarland 0,5 com sol. salina / Разреждение на Бактериална суспензия 0,5 McFarland сс Сольн пазмоор / Zředění bakteriální suspenze 0.5 McFarland fyziologickým roztokem / Fortynding 0.5 McFarland bakteriel suspension med saltvand / Αραίωση Βακτηριακό εναιώρημα 0,5 McFarland με φυσιολογικό ορό / Bakterisuspensiooni lahjendamine 0,5 McFarland füsioloogilise lahusega / Razrzedjenje bakterijske suspenzije 0,5 McFarland sa solozuoinom fiziolozhkom otopinom / Atskaidījums Bakteriju suspensija 0,5 McFarland ar šķīd. Sāls šķīdums / 0,5 McFarland bakterinēs suspensijos atskiedimas druskos tirpalu / 0,5 McFarland baktēriumsuspenzija higitāsa šķīdolatban / Verduunning bacteriesuspensie 0.5 McFarland met zoutoplossing / Fortynding 0.5 McFarland bakteriell suspension med saltvann / Rozcieńczenie Zawiesiny Bakteryjnej 0,5 McFarland z roztworem soli / Dilutie Suspensie bacteriană 0,5 McFarland cu sol. Salină / Redčitev Bakterijske suspenzije 0,5 McFarland s sol. Salino / Dilúcia bakteriálnej suspenzie 0,5 McFarland so soľným roztokom / Laimennus: 0,5 McFarland-bakteerisuspensio suolaliuoksella / Utspadding: 0.5 McFarland bakteriell koksaletslösning / Tuzlu solüsyonla 0.5 McFarland bakteri süspansiyonu seyreltme.

Product: Prodotto / Producto / Produkt / Produit / Produto / Προϊόν / Výrobek / Produkt / Προϊόν / Toode / Proizvod / Produkts / Gaminys / Termék / Product / Produkt / Produkt / Produks / Izdelek / Výrobok / Tuote / Produkt / Ürün

Lot Number: Numero di Lotto / N.° de lote / Chargen-nummer / Numéro de lot / Número do lote Номер пармуда / Číslo šarže / Partinummer / Αριθμός παρτίδας / Partii number / Broj serije (Lot) / Partijas numurs / Partijos numeris / Tételszám / Partijnummer / Partinummer / Numer Partii / Numero di Lotto / Šifra serije / číslo šarže / Eränumero / Parti-nummer / Parti numaras

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In accordance with Clinical and Laboratory Standards Institute M40-A2, viability performance is measured for each test organism at the 48 hrs time point and compared with the acceptance criteria.

In both the Roll-Plate and Swab Elution viability performance studies, Copan MSwab® System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8°C) and room temperature (20 – 25°C). Acceptable recovery for the Roll-Plate Method is defined as ≥5 CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU. Acceptable recovery for the Swab Elution Method is defined as no more than a $3 \log_{10}$ (1×10^3 +/- 10%) decline in CFU between the zero-time CFU count and the CFU of the swabs after the specified holding time.

Additional time-points were tested for *Staphylococcus aureus*, ATCC® 29213 and ATCC® 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC® 43300 and ATCC® 700698.

In the Roll-Plate viability performance studies, Copan MSwab® System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8 °C) for 14 days and room temperature (20 – 25 °C) for 72 hrs. Acceptable recovery for the Roll-Plate Method is defined as ≥ 5 CFU over the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU.

Viability performance studies also include an assessment of bacterial overgrowth at refrigerated temperatures (4 – 8 °C). For the Swab Elution Method, an over- growth assessment is made on all bacteria species tested at the 48 hrs holding time point. Overgrowth assessment using the Swab Elution Method is defined as greater than 1 log₁₀ increase in CFU between the zero-time CFU count and the holding time point. For the Roll-Plate Method, an overgrowth assessment is made with a separate analysis in which swabs are dosed with 100µl containing 10² CFU of *Pseudomonas aeruginosa* culture. Overgrowth under these conditions is defined as greater than 1 log₁₀ increase in CFU between zero-time CFU and the 48 hrs holding time point. Copan MSwab® System demonstrated no overgrowth based on the acceptance criteria described in Clinical and Laboratory Standards Institute M40-A2.

Copan MSwab® System was able to maintain the viability of the following organisms for at least 48 hours at both room temperature (20-25 °C) and in the refrigerator (2-8 °C) under the test conditions described above: Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2.

TABLE OF SYMBOLS

See the table of symbols at the end of the instructions for use.

NOTES FOR THE PROFESSIONAL USER

In the event that a serious incident occurs in relation to this device, it must be reported to the manufacturer (see the contacts at the end of the Instructions for Use) and the competent authority in the country where the user and/or patient is located.

REVISION HISTORY

Last Revision No.*	Release date	Changes made
01	10-2022	Revision of IFU sections (first revision in IVDR)

* Should you need earlier revisions, contact Copan Customer Service.

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BIBLIOGRAPHY

1. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. *Manual of Clinical Microbiology*. 7th edition. Washington, DC: ASM; 1999.
2. Isenberg HD, Schoenkench FD, Von Graevenitz A. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, S. J. Rubin. American Society for Microbiology, Washington, DC, 1979.
3. Isenberg HD. (Editor in Chief). *Clinical Microbiology Procedures Handbook*. American Society for Microbiology, Washington, DC, 1992.
4. Isenberg HD. (Editor in Chief). *Essential Procedures for Clinical Microbiology*. American Society for Microbiology, Washington DC, 1998.
5. Summanen P, Baron EJ, Citron D, Strong C, Wexler HM, Finegold SM. (1993). *Wadsworth Anaerobic Bacteriology Manual*, 5th edn. Star Publishing Company, Belmont, CA.
6. ASM Cumitech publications number 10A, 13A, 17A. American Society for Microbiology, Washington, DC.
7. Miller JM. A Guide to Specimen Management in Clinical Microbiology. Second Edition. American Society for Microbiology. Washington, DC. 1999.
8. Miller JM, Holmes HT. Specimen collection, transport, and storage. In: *Manual of Clinical Microbiology*. 6th ed. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. Washington, DC: ASM; 1995:19-20.
9. Quality Control of Microbiology Transport Systems. M40-A Vol 23 No. 34. December 2003. Clinical and Laboratory Standards Institute (formerly NCCLS).
10. Sng E-H, Rajan VS, Teo K-L, Goh A-J. The recovery of *Neisseria gonorrhoeae* from clinical specimens: effects of different temperatures, transport times, and media. *Sex Trans Dis*. 1982; 9:74-78.
11. Sun Y, Taylor T, Williams L, Sautter RL. Comparison of bacterial viability using both the EZ brand collection and transport system with the Difco swab transport pack. Presented at: 96th ASM General Meeting. 1996; Washington DC. Abstract C35.
12. Arbiqque JC, Forward KR, LeBlanc J. Evaluation of four commercial transport media for the survival of *Neisseria gonorrhoeae*. *Diagnostic Microbiology and Infectious Disease*. 2000; 36:163-168.
13. Perry JL. Effects of temperature on fastidious organism viability during swab transport. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-55.
14. Wilson DA, Tuohy MS, Procop GW, Hall GS. Effects of storage on the recovery of bacteria from three swab transport systems: BD CultureSwab, BD Culturette and Starplex StarSwab II. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-61.
15. Arbiqque J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies described in NCCLS document M40-P Quality Control of Microbiology Transport Devices. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-40.
16. Mitchell E, Berman M, Ginocchio CC. Evaluation of two new Liquid Stuart transport systems: Platinum StarSwab II (Starplex Scientific) and BBL CultureSwab (Becton Dickinson). 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-74.
17. Perry JL, Matthews JS. Compliance of two popular swab transport systems with performance standards detailed by the new NCCLS Proposed Standard, M40-P. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-42.
18. Robinson A, Gruver ML. Comparison of bacterial survival in two transport systems stored at room temperature and refrigerator temperatures. 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-69.
19. Human RP, Jones GA. Evaluation of 4 transport systems against a published standard. 104th General Meeting of the American Society for Microbiology. 2004; New Orleans, LA. Abstract C-161.
20. Human RP, Jones GA. Evaluation of swab transport systems against a published standard. *J Clin Pathol* 2004; 57:762-763.
21. Arbiqque J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies for anaerobic organisms described in NCCLS document M40-P, Quality Control of Microbiology Transport Devices. 13th European Congress of Clinical Microbiology and Infectious Disease (ECCMID). 2003; Glasgow, UK. Abstract P-652.
22. Marler LM, Siders JA, Allen SD. *Direct Smear Atlas, A Monograph of Gram-Stained Preparations of Clinical Specimens*. Lippincott Williams and Wilkins, 2001.
23. Rotimi VO, Yakubu Z, Abudu OO, Banjo TO. Direct Gram's stain of vaginal discharge as a means of diagnosing bacterial vaginosis. *Journal of Medical Microbiology*, 1991 Vol 35, Issue 2 103-106.
24. Spiegel CA, Ansel R, Holmes KK. Diagnosis of bacterial vaginosis by direct gram stain of vaginal fluid. *J Clin Microbiol*. 1983 Jul;18 (1):170-177.
25. Benavides MI, Moncada X, Rodriguez B, Castillo C. Gonococcal urethritis in men: clinical experience in 1978-1988. *Rev Med Chil*. 1992 Oct;120(10):1140-3.
26. Mayaud P, Msuya W, Todd J, Kaatano G, West B, Begkoyian G, Grosskurth H, Mabey D. Rapid assessment in Rwandan refugee camps in Tanzania. *Genitourin Med*, 1997 Feb;73 (1):33-8
27. Lynne S. *Garcia Clinical Microbiology Procedures Handbook*, Third Edition 2010 ASM Press 3.2.1
28. AE Greenberg, LS Clesceri, and AD Eaton. 9215 heterotrophic plate count. In: *Standard Methods for the Examination of Water and Waste Water*. 18th ed. Washington, DC APHA; 1992: 9-33-9-34
29. Henry D. Isenberg *Essential procedures for Clinical Microbiology* 1998 ASM PRESS 477-478
30. Steven Specter, Richard L. Hodinka, Stephen A. Young *Clinical Virology Manual* Third edition, 384-409
31. Fleming D. *Biological Safety: Principles and Practices*. January 2000. ASM, Washington DC.
32. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA. <http://www.cdc.gov/od/ohs/symp5/jyrtxt.htm>.
33. Richardson JH. *Biosafety in Microbiological and Biomedical Laboratories*. December 1994. Diane Publishing Company.
34. Hansen DJ. *Healthcare, Laboratories and Biosafety*. Vol 2., 1992. CRC Press.

Index of Symbols / Tabella dei Simboli / Tabla de símbolos / Symboltabelle / Table des Symboles / Tabela de símbolos / Таблица на Символите / Tabulka symbolů / Tabel med symboler / Πίνακας συμβόλων / Sümbole tabele / Tablica simbola / Simbolu tabula / Simbolių rodyklė / Szimbólumtáblázat / Tabel symbolen / Symboltabelle / Tabela symboli / Tabelul cu Simboluri / Tabela simbolov / Tabuľka symbolov / Symbolitaulukko / Symbolförteckning / Semboller Tablosu

Symbol / Simbolo / Simbolo / Symbole / Symbole / Simbolos / Символ / Symbol / Symbol / Σύμβολο / Sümbole / Simbol / Simbol / Simbols / Simbols / Szimbólum / Symbol / Symbol / Symbol / Symbol / Symbol / Symbol / Symboli / Symbol / Sembol	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Значение / Význam / Betyder / Σημασία / Tähendus / Značenje / Nozime / Reikšmė / Jelenítés / Betekenis / Betydning / Znaczenie / Semnificație / Pomen / Význam / Merkitys / Betyder / Anlami
	Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante / Производител / Výrobce / Producent / Κατασκευαστής / Tootja / Proizvođač / Ražotājs / Gamintojas / Gyártó / Fabrikant / Produzenten / Producent / Producător / Proizvajalec / Výrobca / Valmistaja / Tillverkare / Üretici
IVD	In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro / Dispositivo de diagnóstico in vitro / Устройство за Диагностика инвитро / In vitro diagnostický prostředek / Diagnostisk medicinsk anordning in vitro / In vitro διαγνωστική συσκευή / In vitro diagnostikaseade / In vitro diagnostički proizvod / In vitro diagnostikas ierice / „In vitro“ diagnostikos prietaisai / In Vitro diagnosztikai készülék / Hulpmiddel in-vitrodiagnostiek / In vitro diagnostisk medisinsk utstyr / Urządzenie do diagnostyki in vitro / Dispozitiv pentru diagnostic in vitro/ In vitro diagnostična naprava / Diagnostická pomůcka in vitro / In vitro -diagnostikkalaita / Medicinsk anordning för in vitro-diagnostik / In Vitro Tanılma Cihazı
MD	Medical Device / Dispositivo Medico / Producto sanitario / Medizinprodukt / Dispositif médical / Dispositivo Médico / Медицински Продукт / Zdravotnický prostředek / Medicinsk utstyr / Ιατρική συσκευή / Meditsiiniline seade / Medicinski proizvod / Medicinska ierice / Medicinos prietaisai / Orvosi eszköz / Medisch hulpmiddel / Medisinsk utstyr / Wyrób Medyczny / Dispozitiv Medical / Medicinski pripomoček / Zdravotnícka pomůcka / Lääkinnällinen laite / Medicinteknisk produkt / Tibbi Cihaz
UDI	Unique Device Identifier / Identificativo unico del dispositivo / Identificador único de dispositivo / Einmalige Produktkennung / Identifiant dispositif unique / Identificador de dispositivo único / Уникален идентификатор на устройството / Jednoznačný identifikátor prostředku / Unik udstysidentifikationskode / Μοναδικό αναγνωριστικό συσκευής / Seadme unikaalne identifikaator / Ierices unikālais identifikators / Unikalus įrenginio identifikatorius / Az eszköz egyedi azonosítója / Unieke identificatie hulpmiddel / Identifikator for Unique-enhet / Identificare unică a dispozitivului / Unikatowy numer urządzenia / Enolični identifikator naprave / Jedinečný identifikátor pomůcky / Yksilöivä laitteenunite / Unik produktidentifisering / Benzersiz Cihaz Tanımlayıcı
CE	CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE / Marcação CE / Маркировка CE / Označení CE / CE-mærke / Σήμανση CE / CE-mārgis / Oznaka CE / CE markējums / CE ženklas / CE-jelölés / CE-markering / CE-merke / Oznakowanie CE / Marcajul CE / Oznaka CE / Označenje CE / CE-merkintä / CE-märkning / CE İşareti
CE 0123	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Идентификационный номер на нотифицирания орган / Identifikační číslo oznámeného subjektu / Identifikationsnummer for det bemyndigede organ / Αριθμός αναγνώρισης του κοινοποιημένου οργανισμού / Teavitatud asutuse identifitseerimisnumber / Jedinšteni identifikator proizvoda / Identifikacijski broj obaviještenog tijela / Plinvarotās iestādes identifikācijas numurs / Notifikuotosios įstaigos identifikavimo numeris / A bejelentett szervezet azonosító száma / Identificatienummer aangemelde instantie / Identifikasjon for godkjenningsorganet / Numer identyfikacyjny jednostki notyfikowanej / Numărul de identificare al organismului notificat / Identifikacijska številka priglasenega organa / Identifikačné číslo notifikovaného orgánu / Ilmoitettu laitoksen tunnusnumero / Identifiseringsnummer av certifierad myndighet / Onaylı kuruluş tanımlama numarası
STERILE EO	Sterilized using ethylene oxide / Sterilizzato usando ossido di etilene / Esterilizado por óxido de etileno / Sterilisiert mit Äthylenoxid / Stérilisé à l'aide d'oxyde d'éthylène / Esterilizado utilizando óxido de etileno / Стерилизовано с этилен оксид / Sterilizováno pomocí ethylenoxidu / Steriliseret med ethylenoxid / Αποστειρωμένο με τη χρήση οξειδίου του αιθyleνίου / Steriliseeritud, kasutades etüleenoksiidi / Steriloitu etyleeniksiidillä / Sterilizirano etilen-oksidom / Etilén-oxidall sterilizálva / Sterilizuota etileno oksidu / Sterilizēts, izmantojot etilēnoksidu / Gesteriliseerd met ethyleenoxide / Steriliseret med etylenoksid / Sterylizowany tenkiem etylenu / Sterilizat cu oxid de etilenă / Sterilizované etylenoxidom / Sterilizirano z etilenoksidom / Sterilisano etilen-oksidom / Steriliserad med etylenoxid / Etilen oksit kullanılarak sterilize edilmiştir / Стерилизований завдяки етиленоксиду

	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a usar / Не используйте повторно / Nepoužívejte opakovaně / Nepoužívejte znovu / Μην επαναχρησιμοποιείτε / Mitte taaskasutada / Ne upotrebljavajte ponovno / Nelietot atkārtoti / Nenaudoti pakartotinai / Ne használja fel újra / Niet hergebruiken / Må ikke gjenbrukes / Nie używać ponownie / De unică folosință / Ne uporabljajte ponovno / Nepoužívajte znovu / Ei saa käyttää uudelleen / Nemojte ponovo koristiti / Yeniden kullanmayın
	Single sterile barrier system / Sistema a barriera sterile singola / Sistema de barrera estéril único / Einfaches Sterilbarrieresystem / Système de barrière stérile unique / Sistema de barreira estéril simples / Система с единична стерилна бариера / Systém samostatné sterilní bariéry / Enkelt sterilisert barriersystem / Ένιαίο σύστημα φραγμού αποστείρωσης / Ühekordne steriilne barjäärisüsteem / Jednostruki sterilni sustav / Vienas sterilas barjeras sistēma / Vienkartinė steriliaus barjero sistema / Egyes steril korlátrendszer / Enkel steril barrièresysteem / Sterilt barriersystem til engangsbruk / System pojedynczej bariery sterylnej / Sistem cu barieră sterilă unică / Enojni sterilni pregradni sistem / Systém s jednou sterilnou bariérou / Yksittäinen steriili estejärjestelmä / Enkelt sterilisert barriärsystem / Tek kullanımlık steril ambalaj sistemi
	Do not resterilize / Non risterilizzare / No reesterilizar / Nicht resterilisieren / Ne pas restériliser / Não reesterilizar / Не стерилизуйте повторно / Neprovádějte opakovanou sterilizaci / Må ikke gensteriliseres / Μην επαναστείριζετε / Mitte uesti steriliseerida / Ne sterilizirajte ponovno / Nesterilizēt atkārtoti / Nesterilizuokite pakartotinai / Ne sterilizálja újra / Niet opnieuw steriliseren / Må ikke steriliseres på nytt / Nie sterylizować ponownie / Nu reesterilizați / Ne sterilizirajte ponovno / Nesterilizujte znovu / Ei saa steriloida uudelleen / Får inte omsteriliseras / Tekrar sterilize etmeyin
	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Referência do catálogo / Каталоген номер / Číslo katalogu / Katalognummer / Αριθμός καταλόγου / Katalogi number / Kataloški broj / Kataloga numurs / Katalógusszám / Katalogo numeris / Catalogusnummer / Katalognr / Numer katalogowy / Număr de catalog / Kataloška številka / Katalógové číslo / Luettelonumero / Katalognummer / Katalog numarası
	Temperature limitation / Limiti di temperatura / Límites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Температурни граници / Teplotní limity / Teplotní limity / Ορια θερμοκρασίας / Temperatuuri piirid / Temperaturno ograničenje / Temperatūras ierobežojumi / Temperatūros apribojimas / Hőmérséklethatárok / Temperatuurgrenzen / Temperaturgrenser / Ograniczenie temperatury / Intervale de temperatură / Temperaturne omejitve / Teplotné obmedzenia / Lämpötilarajat / Opseg temperature / Sicaklık kısıtlamaları
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Используйте до / Spotřebujte do / Anvendes før / Χρήση έως / Kasutusaeg / Upotrijebite do / Izlietot līdz / Tinka naudoti iki / Felhasználás a következőn belül / Gebruiken vóór / Må brukes innen / Użyć do / A se folosi până la / Uporabite do / Spotřebujte do / Viimeinen käyttöpäivä / Ska användas innan / Son kullanim tarihi
	Peel / Strappare per aprire / Desprender / Abziehen / Décoller / Destacável / Откъснете за да отворите / Roztrhněte za účelem otevření / Riv opp for å åpne / Σκίσιμο για άνοιγμα / Avamiseks rebida / Oguliti / Noplēst, lai atvērtu / Atplēšti / Openscheuren / A kinyitáshoz vegye le a kupakot / Skall för åpen / Rozerwij, aby otworzyć / Rupeți pentru a deschide / Odtřgajte, da odprete / Roztrhnúťim otvorte / Kuori / Skal för öppen / Soyarak açın
	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Код на партида (партида) / Kód šarže (dávky) / Seriennummer (parti) / Κωδικός παρτίδας (batch) / Partii kood / Broj serije (Lot) / Partijas (sērijas) kods / Partijos kods (partija) / Tételszám (elkönyvelés) / Code partij / Batch-nummer (parti) / Kod partii (seria) / Codul lotului (lotul) / Šifra serije (serija) / Kód šarže / Eräkoodi (erä) / Dra för att öppna / Partii kodu (lot)
	Contains sufficient for <n> tests / Contenu suffisant pour <n> test / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contémo suficiente para <n> testes / Достатъчно съдържание за <n> тестове / Dostatečný obsah pro <n> testů / Inhold tilstrækkeligt til <n> prøver / Περιεχόμενο επαρκές για <n> δοκιμή / Sisaldab varustust <n> prooviks / Sadržaj dovoljno za <n> ispitivanja / Saturs ir pietiekams <n> testiem / Pakanka <n> tyrimų / Elegendő tartalom <n> teszt számára / Inhoud voldoende voor <n> tests / Innhold tilstrekkelig for <n> test / tester - Zawiera materiał pozwalający na wykonanie <n> testów / Conținut suficient pentru <n> teste / Vsebina zadostuje za <n> testov / Obsah dostatočný na <n> testov / Sisältää riittävästi <n> testiin / Innehåller tillräckligt för <n> tester / <n> test için yeterli içerik

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**Copan Italia S.p.A.**

Via F. Perotti, 10
25125 Brescia, Italy
Tel +39 030 2687211
Fax +39 030 2687250

Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:

Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562, USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com