SRK[®]

Instructions for Use

English	3
Italiano	7
Español	11
Deutsch	15
Français	20
Português	24

Intended Use

Surface sampling for monitoring the cleanliness of the environment or investigations for searching for bacterial or fungal contamination.

Summary and principle

Bioburden monitoring in pharmaceutical/food/hospital and similar industry environments plays an important role in internal control for verifying cleanliness, investigating hospital outbreaks and food contamination, verifying surface basal contamination or searching for a single microorganism⁽¹⁾. The broad range of media and conformation proposals available in the SRK® (Hygiene System) family of products represents the ideal tool according to the type of test to be performed. The SRK® medium allows the conservation of the microbial load during the transport phase of the sample up to the moment of its processing, thus allowing a semi-quantitative estimate of the microbial load present on the surface. The flocked swab supplied with some codes is a valid solution for optimising both surface recovery and the release of microorganisms in the transport medium.

Reagents/ Medium Formulation

Sodium Chloride
Potassium Chloride
Calcium Chloride
Lecithin
Sodium thiosulphate
Sodium Thioglycollate
Sodium Pyruvate
Sodium Bicarbonate
Polysorbate 80

Product codes

Preferred application	Code	Description
Food and cosmetics	902C	16X100 TUBE 10 ML SRK® SOLUTION BLUE PLASTIC APPLICATOR RAYON TIPPED
	904C	16X100 TUBE 10 ML SRK® SOLUTION LARGE PLASTIC APPLICATOR RAYON TIPPED
	905C	SQUARE SAMPLING KIT STERILE SWABS + FLEXIBLE MASKS TUBES WITH NEUTRALIZING SOLUT.
	906C	SRK® 2,5ML 12X86 FLOC.AP.200PKG
	906C.CR	SRK® 2,5ML 12X86 FLOCKED APPLICATOR TRIPLE PACK
	908C	16X100 TUBE 2,5 ML SRK® SOLUTION BLUE PLASTIC APPLICATOR RAYON TIPPED
	9U017S	SRK® SURFACE RECOVERY KIT 1 REGULAR FLOQSwabs® APPLICATOR 10 ML COPAN NEUTRALIZING TRANSPORT MEDIUM IN TUBE INDIVIDUALLY WRAPPED
	9U025S01	SRK® SURFACE RECOVERY KIT 1 REGULAR FLOQSwabs® APPLICATOR, 20 MM BREAKING POINT 1 ML SRK® MEDIUM IN SCREW CAP TUBE INDIVIDUALLY WRAPPED
	922C	SWAB RINSE KIT(SRK®) 10 ML SINGLE PLASTIC APPLICATOR VISCOSA TIPPED
	924C	SWAB RINSE KIT(SRK®) 2,5 ML SINGLE APPLICATOR VISCOSA TIPPED
	926C	SWAB RINSE KIT(SRK®) 10 ML IN BULK SINGLE PLASTIC APPLICATOR VISCOSA TIPPED
Pharmaceutical	940C	SRK® SPONGE DEVICE WITH 10 ML SRK® SOLUTION
	9U004S	HYGIENE PHARMA KIT 10ML SRK® SOLUTION WITH REGULAR FLOQSWAB IN TUBE IN TRIPLE PACKAGING
	952C	SRK® PHARMA 1ML W/REG.APP.
	958C	SRK® PHARMA 1 TUBE, TRIPLE PACK REGULAR FLOCKED SWAB
Hospital	9U011S01	SRK® HYGIENE MONITORING SYSTEM – SINGLE WRAPPING REGULAR FLOQSwabs® IN TUBE WITH 1ML SRK® PREMOISTENING SOLUTION
	9U012S	SRK® HYGIENE MONITORING KIT – SINGLE WRAPPING 1 REGULAR FLOQSwabs® IN TUBE WITH 1 ML SRK® PREMOISTENING SOLUTION 1 ML SRK® SOLUTION IN SCREW CAP TUBE

Technical Notes

SRK® formula contains Sodium Thioglycollate, an important component for the performance of the product and the support of organism viability. Sodium Thioglycollate has a natural sulfur-like odor. It may be possible to detect this odor momentarily when first opening the SRK® peel pouch. This odor is a perfectly normal and completely harmless characteristic.

Direct contact with Sodium Thioglycollate may cause skin irritation.

Foam during vortexing may occur. Foaming enables a better re-suspension of the sample for an effective neutralization action.

Precautions, warnings and limitations

1. This product is not intended to be used for transport and/or enrichment of clinical specimens in diagnostic procedures.
2. This product is for single use only; its reuse could lead to a risk of cross-contamination and/or to inaccurate results.
3. Adopt approved biohazard precautions and use aseptic techniques as appropriate.
4. The use of the product is reserved exclusively to specially trained and qualified personnel.
5. After use, all samples and the materials used to process them must be handled with appropriate precautions and disposed of according to applicable regulations.
6. Read and follow the instructions carefully.
7. Do not resterilise unused products before use.
8. Do not tamper with the product packaging until use.
9. The product is not suitable for applications other than its intended use.

10. The product is not suitable for the collection, transport and conservation of viruses for their processing with culture techniques.
11. The product can be used for molecular biology investigations after internal laboratory validation. The product is not suitable for chemical analysis.
12. Do not use if the product is visibly damaged (e.g., if the swab tip or applicator is broken or if the tube is leaking).
13. The use of this product combined with rapid kits or other equipment must be validated by the user before use.
14. Do not ingest the medium.
15. For codes 906C.CR, 958C and 9U004S: do not use the swab to directly plate the medium.
16. Not for use on humans
17. Not for use for human diagnostics

The product is intended for surface sampling only. Do not use for sampling of live or dead people or animals and of food intended for consumption.

Storage

The product is ready to use and requires no further preparation. The product must be stored in the original packaging until use at the temperature indicated on the label. Do not incubate or freeze before use. In case of incorrect storage conditions, the effectiveness could be compromised.

Product deterioration

Do not use if: (1) the product shows visible marks of damage or contamination; (2) there is evidence of leakage; (3) the expiration date has passed; (4) there are other signs of deterioration (i.e., medium is turbid).

Specimen storage and transportation

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations. Transport of specimens should comply with internal guidelines. All specimens should be processed as soon as they are received in the laboratory.

Code	Sample stability
902C	4°C – 25°C / 24h
904C	
905C	
906C	
906C.CR	
908C	
9U017S	
9U025S01	
922C	
924C	
926C	
940C	
9U004S	
958C	
9U012S	
952C	Swab shall immediately streak the onto the plate and incubate as per internal procedures.
9U011S01	

Materials required but not supplied

Appropriate materials for isolating and culturing bacteria, yeasts, and moulds. These materials include culture media plates or tubes and incubation systems, gas jars or anaerobic workstations. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques.

Appropriate material for delimiting the surface to be sampled is available from Copan. Please refer to COPAN S.p.A. customer service or to your local representative for specific product codes.

Directions for use

The SRK® line has a dedicated pharma products line which have a triple wrapping. This packaging configuration allow to qualify the use of the product with isolator systems. The triple wrapping is validated to be resistant to VHP (vaporized hydrogen peroxide) treatment. The inner layers are already sterilized by gamma radiation 10^{-6} and does not require VHP treatment for sterilization.

Wear gloves and prepare all the necessary materials for the collection. Specific precautionary measures may be required in order to avoid contamination of both the sampling area and the environment. Please refer to your standard operational procedures.

Specimen collection

For correct use, refer to the instructions in the specific sections according to the codes indicated.

Devices with separate premoistening

Codes **902C, 904C, 905C, 906C, 906C.CR, 908C, 9U004S, 9U017S, 9U025S01.**

1. In case of triple wrapped product (i.e. 906C.CR, 9U004S), open the external bag and remove the content. Open the second bag and remove the peel pouch.
2. Open the peel-pouch and take out all the components. Place the SRK® tube in a suitable rack. Discard the peel-pouch.
3. Pretwetting step: aseptically unscrew the tube and moisten ONLY the tip of the swab with the medium by immersing it in the tube, making sure not to spill the medium. Close the tube with the cap.
Note: If pretwetting is not necessary, proceed with step 4
4. Proceed with sampling as follows:
 - A) Flat surfaces: sampling in 3 directions (vertical, horizontal and diagonal).

NOTE: sampled surface shall not exceed 100 cm²

B) Articulated surfaces: sampling all possible hidden areas (corners, angles etc...)

5. Take from the rack the SRK® tube while holding the swab, unscrew the cap from SRK® tube making sure not to spill the medium. Insert the swab inside the tube with the breaking point against the rim of the tube. Break the shaft and discard the upper part of the swab. Close the cap of the tube and secure it tightly.
6. Identify the sample and transport the swab within the sample stability conditions described in the paragraph Specimen storage and transportation.

Devices with premoistening tube

Codes **9U012S, 958C**

1. In case of triple wrapped product (958C), open the external bag and remove the content. Open the second bag and remove the peel pouch.
2. Open the peel-pouch and take out all the components. Place the SRK® tube in a suitable rack. discard the peel-pouch.
3. Prewetting step: squeeze the sponge at the bottom of the tube with thumb and index finger until the tip of the swab is wetted.
Note: If prewetting is not necessary, proceed with step 4
4. Open the tube and take out the swab. Discard the pre-moistening tube with the sponge
5. Proceed with sampling as follows:
A) Flat surfaces: sampling in 3 directions (vertical, horizontal and diagonal).
NOTE: sampled surface shall not exceed 100 cm²
B) Articulated surfaces: sampling all possible hidden areas (corners, angles etc...)
6. Take from the rack the SRK® tube while holding the swab, unscrew the cap from SRK® tube making sure not to spill the medium. Insert the swab inside the tube with the breaking point against the rim of the tube. Break the shaft and discard the upper part of the swab. Close the cap of the tube and secure it tightly.
7. Identify the sample and transport the swab within the sample stability conditions described in the paragraph Specimen storage and transportation.

Codici **952C e 9U011S01**

1. In case of triple wrapped product (952C) open the external bag and remove the content. Open the second bag and remove the peel pouch.
2. Open the peel-pouch and take out all the components. Remember to discard the peel-pouch.
3. Prewetting step: squeeze the sponge at the bottom of the tube with thumb and index finger until the tip of the swab is wetted.
Note: If prewetting is not necessary, proceed with step 4
4. Open the tube and take out the swab. Place the tube in a suitable rack.
5. Proceed with sampling as following:
A) Flat surfaces: sampling in 3 directions (vertical, horizontal and diagonal).
NOTE: sampled surface shall not exceed 100 cm²
B) Articulated surfaces: sampling all possible hidden areas (corners, angles etc...)
6. Immediately streak the swab and incubate the plate as per internal procedures.
NOTE: If transport of the swab is necessary, replace the swab inside its tube previously identified. Identify the sample and transport the swab within the sample stability conditions described in the paragraph Specimen storage and transportation.

Applicator immersed in the medium glued to the cap

Codes **922C, 924C, 926C**

1. Open the peel-pouch and take out all the components (not necessary for products in bulk, i.e. 926C). Remember to discard the peel-pouch.
2. Unscrew cap from the tube. Place the tube in a suitable rack.
3. Proceed with sampling as following:
A) Flat surfaces: sampling in 3 directions (vertical, horizontal and diagonal).
NOTE: sampled surface shall not exceed 100 cm²
B) Articulated surfaces: sampling all possible hidden areas (corners, angles etc...)
4. Close the tube.
5. Identify the sample and transport the swab within the sample stability conditions described in the paragraph Specimen storage and transportation.

Sponge and premoistening

Code **940C**

1. Open the peel-pouch and take out all the components. Remember to discard the peel-pouch.
2. Place the tube in a suitable rack.
3. Use the cap with the sponge and proceed with sampling as following:
A) Flat surfaces: sampling in 3 directions (vertical, horizontal and diagonal).
NOTE: sampled surface shall not exceed 100 cm²
B) Articulated surfaces: sampling all possible hidden areas (corners, angles etc...)
4. Take from the rack the SRK® tube while holding the cap with the sponge, unscrew the cap from SRK® tube making sure not to spill the medium. Discard the green cap and use the cap with the sponge to close the tube.

Processing in the laboratory

For all codes except 952C, 9U011S01.

- Do not use the swab to inoculate plates;
 - Shake the tube briefly to homogenise the sample in the medium before performing any type of analysis;
 - Process the sample by plating 100ul directly from the tube onto non-selective medium.
NOTE: for samples where a high bioburden is expected, it is possible to perform serial dilutions of the original sample in a suitable diluent (also SRK®).
- NOTE: Analyses on selective media or passages in enrichment broths must first be validated internally.

For codes 952C, 9U011S01.

- Use the swab to directly inoculate a non-selective culture plate.
NOTE: for samples where a high bioburden is expected, it is possible to dilute the original sample in a suitable diluent
NOTE: Analyses on selective media or passages in enrichment broths must first be validated internally.

Quality control and performance characteristics

The product has been tested for its ability to keep a panel of strains representative of the environmental flora viable (gram+, gram-, aerobes, anaerobes, moulds and yeasts). The product has been tested for sampling capability in accordance with ISO 18593. The pH is tested for each batch produced.

The test aims to demonstrate the retention capacity of the device, indicated as the percentage loss of the CFU value between those found by plating at zero time with respect to the moment of inoculation and those found after a simulated transport time of 24h at room temperature (25°C) and of 24h and 48h at refrigerated temperatures. Specifically, the device is inoculated with a bacterial suspension which guarantees a growth of between 10 and 120 CFU per 100 µL of liquid in a plate with non-selective medium. The device is then maintained for 24h at room temperature (25°C) and for 24h and 48h at refrigerated temperatures as a simulation of transport and after the established time, aliquots are plated on a plate. The value expressed as CFU obtained from the plates after transport and that of the plates at time 0 is entered in the formula for calculating the percentage loss $[(T_{24/T0} - 1) \times 100]$.

The following strains have been verified:

Aspergillus brasiliensis (ATCC® 16404)
Candida albicans (ATCC 10231)
Escherichia coli (ATCC 25922)
Staphylococcus aureus (ATCC 6538)
Staphylococcus epidermidis (ATCC 12228)
Salmonella Typhimurium (ATCC 14028, ATCC 13311, ATCC 13076)
Listeria monocytogenes (ATCC 13932, ATCC 19114, ATCC 19111)
Enterococcus faecalis (ATCC 29212)
Bacillus subtilis ATCC 6633
Candida tropicalis (ATCC 750)
Saccharomyces cerevisiae (ATCC 9763)
Pseudomonas aeruginosa (ATCC 9027, ATCC 27853)*
Campylobacter jejuni (ATCC 33291)**
Staphylococcus aureus MRSA (ATCC 43300, ATCC 33591, ATCC 700698)
Enterococcus faecalis (ATCC 51299, ATCC 51575, ATCC 700221)
Escherichia coli (ATCC BAA-2469)
Klebsiella pneumoniae (ATCC BAA-1705)
Enterobacter cloacae (BAA-2341)

*Tested only at room temperature for 24h

**Tested only at refrigerated temperature for 24h.

The results obtained largely depend on correct and adequate sampling, the type of surface sampled and any disinfectant residues, as well as on the timeliness with which transport and laboratory analyses are performed.

The device has been validated for its neutralising capacity. Each ml of medium is able to neutralise the following classes of disinfectants at the concentrations shown in the following table:

Antimicrobial Family	Tested Chemical	Maximum Concentration Neutralized (ppm)
Oxidizing compounds	Sodium Hypochlorite	12.000
	Peracetic Acid	35.000
Quaternary ammonium compounds	Benzalkonium Chloride Solution	500
Biguanides	Chlorhexidine Digluconate solution	200
Alcohols	Isopropyl Alcohol	700.000
Mix of Phenolic	o-phenylphenol o-benzyl-p-chlorophenol p-tert-amylphenol	10.000

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Copan Italia S.p.A.
Copan Diagnostics Inc.