
LBM[®] CABroth

Instructions for Use

English	3
Italiano	5
Español	7
Deutsch	10
Français	13
Português	15
БЪЛГАРСКИ	18
Česky	20
Dansk	23
ΕΛΛΗΝΙΚΑ	25
Eesti Keel	28
Hrvatski	30
Latviešu	32
Lietuvių K.	35
Magyar	37
Nederlands	40
Norsk	42
Polski	45
Română	47
Slovenčina	49
Slovenščina	52
Suomi	54
Svenska	56
Türkçe	59

Candida Auris Broth (CABroth)

Instructions for use

INTENDED USE

Selective enrichment medium for *Candida auris*.

The CABroth is a prepared medium for the selective enrichment of *Candida auris* for subsequent culture or molecular analysis.

The CABroth is intended for use with compatible molecular assays.

SUMMARY AND PRINCIPLES

Candida auris is a globally emerging, multidrug-resistant fungal pathogen associated with severe nosocomial infections and high mortality rates, particularly in critically ill and immunocompromised patients. Its rapid spread across continents and frequent involvement in healthcare-associated outbreaks have made *C. auris* a critical priority for infection prevention and control efforts ^(1,2).

Effective detection and surveillance of *C. auris* colonization are essential to prevent transmission and guide timely clinical and epidemiological responses. However, early detection remains challenging due to limitations in current diagnostic workflows. One of the key barriers is the low fungal burden in colonized individuals and environmental reservoirs, which may go undetected using standard culture techniques ⁽³⁾.

The implementation of selective and enrichment media specifically designed for *C. auris* represents a crucial step forward in enhancing cultural and molecular-based detection. Traditional media often lack the selectivity required to isolate *C. auris*. ⁽⁴⁾ Non-specific growth conditions can allow competing microbial flora to overgrow the target organism, leading to false-negative results or delayed identification ⁽⁵⁾.

CABroth is a selective enrichment medium specifically formulated to promote the growth of *Candida auris* while suppressing microbial contaminants. Designed for use prior to culture or molecular testing, it enables targeted enrichment of *C. auris* from clinical samples.

REAGENTS

Aminoacids
Inorganic Salts
Organic Salts
Carbohydrates
Antibiotic
Selective Reagent

PRODUCT DESCRIPTION

The product CABroth is ready to use and does not require additional preparation. The product is composed by a polymeric cap and tube containing **2mL of lighter yellow medium.**

The **CABroth product is not sterile.**

The CABroth is available in the following product configurations:

Ref. Number	Product Description	Package	Capturability
4U161N	Selective enrichment medium for <i>Candida auris</i> , 2 mL, in 12x80 mm tube with yellow capture screw-cap.	50 devices per package 50 x 6 devices per box	NOT APPLICABLE The capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service.

PRODUCT STORAGE

The product is ready to use and does not require additional preparation. **The product must be transported and stored in the original container at a temperature of 5-25°C until use. Do not overheat. Do not incubate or freeze prior to use.**

If incorrectly stored, its effectiveness will be compromised. Do not use after the expiry date, clearly printed on the inner and outer boxes and on the label of each tube.

REQUIRED MATERIAL BUT NOT PROVIDED

Suitable material for the specimen collection (eSwab® or FLOQSwabs®) and material for the cultivation and isolation of Fungi, alternatively material for extraction, purification and amplification of nucleic acid. Refer to the laboratory reference manuals for the recommended culture techniques, molecular diagnostics techniques, and identification procedure.

SPECIMEN PRESERVATION

The CABroth must be processed promptly in the laboratory after clinical sample inoculation from eSwab® or FLOQSwabs®.

Analytical results depend on correct sample collection, transportation and storage conditions as well as correct laboratory processing. In case of delay of the processing, refer to internal laboratory reference manuals.

LIMITATIONS

- The product must be used by healthcare professional only.
- Transfer of the clinical sample into the product should be done only through the eSwab® or FLOQSwabs® devices. For information regarding FLOQSwabs® codes suitable to be used with CABroth, please contact customer service.
- The product is not intended for sample transport.**

WARNINGS AND PRECAUTION

- Carefully read and follow the instructions.
- The results obtained depend largely on the correct and adequate collection of the sample, as well as the timeliness of analysis in the laboratory.**

3. For in vitro diagnostic use only.
4. The product is not suitable for any other application different from its intended use. The manufacturer cannot be held responsible for any improper or unqualified use of the product.
5. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
6. The medium is not intended to be used for pre-moisten or pre-wet of swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
7. This product is ready for use and requires no further preparation.
8. **Do not sterilize.**
9. Do not re-pack.
10. Do not freeze.
11. Do not ingest the product.
12. Adopt the approved precautions on biological risks and use aseptic techniques. **The device is to be used by trained and qualified personnel only.** The manufacturer cannot be held responsible for any use by unauthorized or unqualified persons.
13. All clinical samples are considered potentially infectious and must be handled with appropriate precautions. Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens.
14. The use of this product in association with any diagnostic tests or any diagnostic instrumentation must be validated before use, according to the laboratory's internal procedures.
15. The product should not be used if (1) there is evidence of damage or contamination, (2) there is evidence of leakage, (3) the expiration date has passed or (4) there are other signs of deterioration.
16. Check for the presence of medium inside the product before use. Do not use the device if there is no medium inside the product.
17. **Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format and can be identified by the e-IFU indicator on the package label.**
18. Open the test tube only at the time of use.
19. Unused reagents, waste and specimens must be disposed of in accordance with local regulations.
20. The turbidity does not necessarily identify the presence of *Candida auris*. Use diagnostic technologies that can identify the various types of microorganisms.
21. The proper functioning of the swab capture feature by the cap is strictly dependent on the geometry of the swab. Using an incompatible swab may pose contamination risks for both the user and the environment. For information about the correct component combination, please contact Customer Service.

INSTRUCTIONS FOR USE – SPECIMEN COLLECTION

A correct specimen collection represents the most critical phase for the isolation and identification of the pathogenic microorganism. Copan eSwab® and FLOQSwabs® are compatible devices that can be used for sampling. Refer to the instructions for use of the products.

INSTRUCTION FOR USE – LABORATORY PROCESSING

1. Identified a tube of CABroth with the patient ID, unscrew and remove the cap, taking care not to spill the medium. (Fig. 1)

FOR SPECIMEN COLLECTED AND TRANSPORTED TO LABORATORY with FLOQSwabs®:

- 2.a Inoculate the broth by transferring the specimen on swab into the opened tube.
 - **for FLOQSwabs® with break points:** insert the swab into the test tube until the break point is at the level of the test tube opening (Fig. 2). Bend and break the swab at the break point, turning the test tube away from your face. If necessary, gently bend the swab shaft up to 180° and rotate the swab shaft to complete the break (Fig. 3a and Fig. 3b). Discard the upper part of the swab shaft.
 - **for FLOQSwabs® without break points:** insert the swab into the test tube and cut off with a sterile instrument the excess part of the shaft.

FOR SPECIMEN COLLECTED AND TRANSPORTED TO LABORATORY IN Copan eSwab®:

- 2.b Vortex the eSwab® tube for 5-10 seconds at 2000-2500 rpm to homogenize the sample. From eSwab® tube transfer an aliquot between 30 µL and 200 µL into the CABroth.

The best ratio between sample and medium is 1:10.
NOTE: Specimens should be collected early in the course of the disease and specimens give better results if transported to laboratory with Copan eSwab® to maintain microorganisms viability.
3. Reposition the cap on the test tube, hermetically seal it (Fig. 4).
4. Vortex the tube for 5-10 seconds at 2000-2500 rpm in order to mix tube contents.
5. Incubate inoculated CABroth tubes at a temperature between 35 °C and 40 °C (preferably at 40 °C).
6. **After 48 hours**, inoculate 1 to 100 µl of CABroth onto appropriate culture medium (i.e. Sabourad Dextrose Agar, Cromogenic Candida Agar Bioline).


Fig. 1


Fig. 2


Fig. 3.a


Fig. 3.b


Fig. 4

QUALITY CONTROL

CABroth is verified for pH value and its ability to enrich *Candida auris* while inhibiting the growth of unwanted flora.

PERFORMANCE CHARACTERISTICS

CABroth was qualified for its ability in the enrichment of *Candida auris* strains **using both culture and molecular methods**. Enrichment was qualified at both 35 °C and 40 °C. CABroth was qualified for its ability of selective action toward a panel of *Candida* strains and bacteria at both 35 °C and 40 °C. CABroth was qualified for its compatibility with eSwab® transport medium, for transfer of primary sample from liquid Amies medium. CABroth was qualified for its absence of *Candida auris* amplifiable DNA, and its absence of interference with MDx.

Molecular performance tests were performed with a molecular platform (kit. AurisID, IMMY; Fungiplex, Bruker) after extraction with automated system (i.e. Magnetic beads extraction).

STRAIN	ZERO-TIME: CFU/PLATE	TIME 48H CFU/PLATE
<i>Candida auris</i> ATCC [®] , CDC B11903	179	CONFLUENT GROWTH
<i>Candida parapsilosis</i> ATCC [®] 22019	SEMI-CONFLUENT	GROWTH FROM PARTIAL TO TOTAL INHIBITION
<i>Candida albicans</i> ATCC [®] 10231	SEMI-CONFLUENT	GROWTH FROM PARTIAL TO TOTAL INHIBITION
<i>Stafilococcus aureus</i> ATCC [®] 6538	SEMI-CONFLUENT	GROWTH FROM PARTIAL TO TOTAL INHIBITION

NOTE FOR THE PROFESSIONAL USER

In the event that a serious incident occurs in relation to this device, it must be reported to the manufacturer (see the contacts at the end of the Instructions for Use) and the competent authority in the country where the user and/or patient is located.

REVISION HISTORY

Last Revision No.*	Release date	Changes made
00	06-2025	First issue

* Should you need earlier revisions, contact Copan Customer Service.

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BIBLIOGRAPHY

1. CDC (Centers for Disease Control and Prevention). (2023). *Candida auris*. <https://www.cdc.gov/fungal/candida-auris>
2. WHO (World Health Organization). (2022). Fungal Priority Pathogens List to Guide Research, Development and Public Health Action.
3. Kordalewska, M., & Perlin, D. S. (2019). Identification of drug-resistant *Candida auris*. *Frontiers in Microbiology*, 10, 1918. <https://doi.org/10.3389/fmicb.2019.01918>
4. Welsh RM, Bentz ML, Shams A, Houston H, Lyons A, Rose LJ, Litvintseva AP. Survival, Persistence, and Isolation of the Emerging Multidrug-Resistant Pathogenic Yeast *Candida auris* on a Plastic Health Care Surface. *J Clin Microbiol*. 2017 Oct;55(10):2996-3005. doi: 10.1128/JCM.00921-17. Epub 2017 Jul 26. PMID: 28747370; PMCID: PMC5625385.
5. Jeffery-Smith, A., et al. (2018). *Candida auris*: a review of the literature. *Clinical Microbiology Reviews*, 31(1), e00029-17. <https://doi.org/10.1128/CMR.00029-17>

INDEX OF SYMBOLS / INDICE DEI SIMBOLI / ÍNDICE DE SÍMBOLOS / SYMBOLVERZEICHNIS / INDEX DES SYMBOLES / ÍNDICE DE SÍMBOLOS / ИНДЕКС НА СИМВОЛИТЕ / REJSTŘÍK SYMBOLŮ / INDEKS OVER SYMBOLER / ΕΥΡΕΤΗΡΙΟ ΣΥΜΒΟΛΩΝ / TINGMÄRKIDE KIRJELDUSED / KAZALO SIMBOLA / IZMANTOTIE SIMBOLI / SIMBOLIŲ RODYKLĖ / SZIMBÓLUMOK JEGYZÉKE / INDEX VAN SYMBOLEN / SYMBOLINDEKS / INDEKS SYMBOLI / INDEXT SIMBOLURI / ZOZNAM SYMBOLOV / KAZALO SIMBOLOV / SYMBOLIAKEMISTO / LISTA ÖVER SYMBOLENER / SEMBOLLER DIZÍNÍ

<p>Symbol / Símbolo / Símbolo / Symbol / Symbole / Símbolos / Символ / Symbol / Symbol / Σύμβολο / Tingmärk / Symboli / Símbol / Szimbólum / Simbols / Símbols / Symbol / Símbol / Symbol / Símbol / Symbol / Sembol</p>	<p>Meaning / Significato / Signification / Bedeutung / Sens / Significado / Значение / Význam / Betydning / Σημασία / Tähendus / Merkitys / Značenje / Jelentés/ Reikšmė / Nozime / Betekenis / Betydning / Znaczenie / Semnificație / Význam / Pomen / Betydelse / Anlam</p>
	<p>Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante / Производител / Výrobce / Producent / Κατασκευαστής / Toolija / Valmistaja / Proizvodáč / Gyártó/ Gamintojas / Razotājs / Fabrikant / Produzent / Producent / Producător / Výrobca / Proizvajalec / Tilverkare / Üretici</p>
	<p>In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro/ Dispositivo de diagnóstico in vitro / ИнвИТРО диагностично изделие / Diagnostický prostředek pro in vitro diagnostiku / Udstyr til in vitro-diagnostik / Διαγνωστικό προϊόν in vitro / In vitro diagnostika seade / In vitro -diagnostikkalaite / In vitro dijagnostiikki proizvod / In vitro diagnostikai eszköz / „In vitro“ diagnostikos prietaisais / In vitro diagnostikas ierīce / Hulpmiddel voor in-vitrodiagnostiek / In vitro-diagnostisk enhet / Urządzenie do diagnostyki in vitro / Dispozitiv de diagnostic in vitro / Diagnostická pomůcka in vitro / In vitro diagnostični pripomoček / Enhet för in vitro-diagnostik / In vitro tani cihazı</p>
	<p>Unique Device Identifier / Identificativo unico del dispositivo / Identificador único de dispositivo / Einmalige Produktkennung / Identifiant dispositif univoque / Identificador de dispositivo único / Унікален ідентифікатор на изделиє / Jedinečný identifikátor zařízení / Unik enheds-id / Αποκλειστική ταυτοποίηση ιατροτεχνολογικού προϊόντος / Seadme kordumatu identifitseerimistunnus / Yksilöivä laitetunniste / Jedinstveni identifikator proizvoda / Egyedi eszközazonosító / Unikalus įrenginio identifikatorius / Unikālais ierīces identifikators / Unieke hulpmiddelidentificatie / Identifikator for Unique-enhet / Unikatowy numer urządzenia / Identificator unic dispozitiv / Unikátny identifikátor pomůcky / Edinstveni identifikator naprave / Unik enhetsidentifiarare / Benzersiz Cihaz Tanımlayıcı</p>
	<p>CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE/ Marcação CE / CE маркировка / Označení CE / CE-mærkning / Σήμανση CE / CE-mārgis / CE-merkintä / Oznaka CE / CE-jelölés/ CE ženklas / CE marķējums / CE-markering / CE-merking / Oznakowanie CE / Marcaj CE / Oznaczenie CE / Oznaka CE / CE-märkning / CE işareti</p>
	<p>Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Не използвайте повторно / Výrobek nepoužívejte opakovaně / Må ikke genbruges / Μην επαναχρησιμοποιείτε / Æhekordne kasutus / Ei saa käyttää uudelleen / Ne upotrebļavajate ponovno / Ne használja újra / Nenaudoti pakartotini / Nelietot atkārtoti / Niet hergebruiken / Skal ikke brukes flere ganger / Nie używać ponownie / A nu se reutiliza / Nepoužívajte opakovane / Ne uporabljajte ponovno / Får ej återanvändas / Yeniden kullanmayın</p>
	<p>Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Réference du catalogue / Número do catálogo / Каталоген номер / Katalogové číslo / Katalognummer / Αριθμός καταλόγου / Katalognummer / Luettelonumero / Kataloški broj / Katalógusszám/ Katalogo numeris / Kataloga numurs / Catalogusnummer / Katalognummer / Numer katalogowy / Număr catalog / Katalogové číslo / Številka kataloga / Katalognummer / Katalog numarasi</p>
	<p>Temperature limitation / Limiti di temperatura / Límites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Температурно ограничение / Teplotní omezení / Temperaturbegrænsning / Περιορισμός θερμοκρασίας / Temperaturuppiirang / Lämpötilaraja / Temperaturu ograniczenie / Hőmérsékletkorlátozás/ Temperatūras apribojimas / Temperatūras ierobežojums / Temperaturlimit / Temperaturbegrænsning / Ograniczenie temperatury / Limite de temperatură / Teplotné limity / Temperatura omejitve / Temperaturbegrænsning / Sıcaklık kısıtlamaları</p>
	<p>Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Използвайте преди / Spotřebujte do / Udløbsdato / Χρήση έως / Kasuta kuni kuuräevani / Viimeinen käyttöpäivä / Upotrijebite do / Felhasználhatósági idő/ Tinka naudoti iki / Izljetot līdz / Te gebruiken voor / Utløpsdato / Użyć do / A se utiliza înainte de / Použití do / Rok uporabe / Används före / Son kullanın tarihi</p>
	<p>Peel / Strappare per aprire / Desprender / Abziehen / Décoller / Destacar / Отрпепете тук / Odloupněte / Riv op / Avoiyōmo / Eemaldamine / Kuori / Oguliti / Húzza le / Atplēsti / Nolobit / Open trekken / Rives av / Peel / Desprindeți / Odtrhnúť / Sterilizacijska vrečka / Riv av / Soyarak Açın</p>
	<p>Batch code (Lot) / Codice del lotto (partija) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Код на партидата (Серия) / Kód dávky (šarža) / Batchcode (lot) / Κωδικός παρτίδας (Παρτίδα) / Partii tāhis (Lot) / Eräkoodi (erä) / Broj serije (Lot) / Tételkód (tétel) / Partijos kods (partija) / Partijas kods (partija) / Batchcode (partij) / Batchkode (lot) / Kod partii (Seria) / Cod lot (lot) / Kód dávky (šarža) / Serijska oznaka (šarža) / Satskod (Lot) / Parti kodu (Lot)</p>

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Email: info@copangroup.com
Website: www.copangroup.com

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Copan Italia S.p.A.