

Virus, Chlamydia, Mycoplasma and Ureaplasma

UTM

UNIVERSAL TRANSPORT MEDIUM

• *Room Temperature Stable* •

Product Features

UTM System is available as filled tubes in bulk or as patient sample collection packs that incorporate various combinations of polyester fiber or flocked swabs and media tubes.

Safe, shatterproof polypropylene skirted tube with distinctive internal conical shape enabling centrifugation of samples and allowing tubes to stand upright on the laboratory bench

Product can be conveniently stored at room temperature prior to use for collecting samples for all four organism groups; Viruses, Chlamydia, Mycoplasma and Ureaplasma

Choice of regular UTM and mini UTM tubes to suit different work patterns




Large 3ml medium fill volume enough for multiple tests on the same specimen or 1ml medium fill volume to prevent dilution effect.

Media formulation includes antibiotics to inhibit bacterial and fungal flora in patient samples

Three glass beads in each tube facilitates the release and dispersion of patient sample material and virus particles from the swab during vortexing

Wide selection of regular and minitip specimen collection applicators available

All swabs feature unique molded breakpoint enabling safe, reliable break off and deposit of swab sticks into transport tubes



FLOQSwabsTM
by COPAN



*For collection, transport,
maintenance and long
term freeze storage*

UTM
UNIVERSAL TRANSPORT MEDIUM

Copan UTM Ordering Information



Cat. No.	Kits of UTM with Flocked Swabs	Packaging
346C	16x100mm size tube, 3ml UTM medium, 1 regular flocked swab, molded breakpoint	6 boxes of 50 kits
307C	16x100mm size tube, 3ml UTM medium, 1 minitip flocked swab, molded breakpoint	6 boxes of 50 kits
305C	16x100mm tube, 3 ml UTM medium, 1 flexible minitip flocked swab, molded breakpoint	6 boxes of 50 kits
359C	12x80mm size tube, 1ml UTM medium, 1 regular flocked swab, molded breakpoint	6 boxes of 50 kits
355C	12x80mm size tube, 3ml UTM medium, 1 regular flocked swab, molded breakpoint	6 boxes of 50 kits
361C	12x80mm size tube, 1ml UTM medium, 1 minitip flocked swab, molded breakpoint	6 boxes of 50 kits
360C	12x80mm size tube, 1 ml UTM medium, flexible minitip flocked swab, molded breakpoint	6 boxes of 50 kits

Additional packaging configurations may be available. Contact Copan or your authorized Copan dealer.



Cat. No.	Kits of UTM with Polyester Swabs	Packaging
302C	16x100mm size tube, 3ml UTM medium 2 regular size plastic/polyester applicators	6 boxes of 50 kits
303C	16x100mm size tube, 3ml UTM medium, 2 plastic/polyester applicators (1 regular, 1 minitip)	6 boxes of 50 kits
328C	16x100mm size tube, 3ml UTM medium, 1 regular size plastic/polyester applicator	6 boxes of 50 kits
340C	16x100mm size tube, 3ml UTM medium, 1 sterile applicator swab, 1 sterile combo wire - plastic minitip swab	6 boxes of 50 kits
343C	16x100mm size tube, 1.5ml UTM medium, no glass beads, 2 regular size plastic/polyester applicator	6 boxes of 50 kits

Additional packaging configurations may be available. Contact Copan or your authorized Copan dealer.



Cat. No.	UTM Tubes in Bulk	Packaging
330C	16x100mm size tube with 3ml UTM medium	6 boxes of 50 kits
331C	30ml tube with 10ml UTM medium	6 boxes of 50 kits
350C	12x80mm size tube, 1mL UTM medium	6 boxes of 50 kits
353C	12x80mm size tube, 3mL UTM medium	10 boxes of 100 pc

Collection of Virus, Chlamydia, Mycoplasma and Ureaplasma



Copan Diagnostics Inc.
26055 Jefferson Ave.
Murrieta, California 92562 USA

Toll Free: (800) 216-4016 (US & Canada)
Phone: (951) 696-6957 Fax: (951) 600-1832
E-Mail: info@copanusa.com

*For collection, transport,
maintenance and long
term freeze storage*



012610VDP

www.copanusa.com