



COPAN

SELECTIVE ENRICHMENT BROTHS

COPAN's Selective Enrichment Broths are easily managed on automated specimen processors, decreasing manual handling and maximizing automation utilization to bring further efficiencies to specimen processing in clinical Microbiology.



Shatterproof plastic tube compatible with automation



Easy storage of room temperature stable enrichment broths



Convenient skirted, flat-bottomed tube stands on the benchtop reducing risk of biohazardous spills



Selenite

- Isolation medium for Salmonella
- Compatible with COPAN FecalSwab™ or original stool sample



LIM Broth

- Modified Todd Hewitt Broth intended for the isolation and preservation of *Streptococcus agalactiae* (Group B Strep)
- Contains antibiotics to inhibit the growth of normal bacterial flora



TSB Salt

- Isolation medium for Methicillin-resistant *Staphylococcus aureus* (MRSA)
- Intended to enrich samples for subculture on selective MRSA agar plates



Thio Broth

- Intended for cultivation of fastidious or slow growing, anaerobic microorganisms
- Contains growth supplements to enrich many fastidious anaerobic bacteria
- Vacuum sealed in PET tube to maintain anaerobic conditions during incubation



BHI Broth

- Enrichment medium for fastidious bacteria, including *Streptococcus spp.* and *Pneumoniae spp.*



GN Broth

- Enrichment medium for Gram-negative bacteria such as *Salmonella* and *Shigella*
- Inhibits Gram-positive bacteria

Ordering Information

CAT. NO.	PRODUCT DESCRIPTION
085CU.A	GN Broth 4mL GN Broth Medium in Polypropylene 12x80mm Screw-Cap Tube with Internal Conical Shape with Black Cap - 1 Case; 6 Packs/Case; 50 Tubes/Pack
474CE.A	BHI Enrichment Broth 3 mL BHI Broth Medium in Polypropylene 12x80mm Screw-Cap Tube with Internal Conical Shape with Purple Cap - 1 Case; 6 Packs/Case; 50 Tubes/Pack
475CE.A	Selenite Enrichment Broth 2 mL Selenite Medium in Polypropylene 12x80mm Screw-Cap Tube with Internal Conical Shape with White Cap - 1 Case; 6 Packs/Case; 50 Tubes/Pack
476CE.A	LIM Enrichment Broth 2 mL LIM Broth Medium in Polypropylene 12x80mm Screw-Cap Tube with Internal Conical Shape with Light Blue Cap - 1 Case; 6 Packs/Case; 50 Tubes/Pack
4C058N.A	TSB Salt Enrichment Broth 2 mL TSB Medium with 2.5% NaCl in Polypropylene 12x80mm Screw-Cap Tube with Internal Conical Shape with Dark Blue Cap - 1 Case; 6 Packs/Case; 50 Tubes/Pack
4U002N	Thioglycollate Broth 4mL Thioglycollate Broth Medium in Polypropylene 12x80mm Screw-Cap Tube with Internal Conical Shape with Dark Blue Cap - 1 Case; 6 Packs/Case; 50 Tubes/Pack



COPAN'S FULL LABORATORY AUTOMATION LINE is a modular, scalable, and customizable solution for automated specimen processing and full specimen management, automated incubation and Digital Microbiology



WASP[®] DT
Daniela Triva

CONTACT COPAN TODAY to streamline and standardize your specimen collection and to learn more about automation in Clinical Microbiology.



Innovating Together™

🌐 COPANUSA.com

☎ [800.216.4016](tel:800.216.4016)

🐦 [@COPANUSA](https://twitter.com/COPANUSA)

✉ info@copanusa.net

📘 facebook.com/COPAN.USA

COPAN Diagnostics, Inc.
26055 Jefferson Avenue ■ Murrieta, CA 92562 USA