

Enrichment and Selective Broths in Automation Ready Tubes



LBM® Enrichment Broths

Copan's enrichment broth line enhances sample processing efficiency, featuring a range of enrichment and selective broths tailored for both aerobic and anaerobic bacteria. Each LBM enrichment broth tube is optimally designed for integration with the WASP® Walk-Away Specimen Processor for streamlined specimen processing.

Automation Ready

Copan enrichment and selective broths are designed for seamless compatibility with automated platforms, reducing manual handling and enhancing efficiency in clinical Microbiology workflows.

Safe and Stable Design

The shatterproof plastic tubes are compatible with automated systems and feature a skirted, flat-bottom design that stands securely on benchtops, minimizing spills.

Convenient Storage






Copan enrichment broths are stable at room temperature, making them easy to store without special temperature controls, simplifying inventory management.

Compliance and Standardization

Fully compliant with ISO 13485-2016 and ISO 9001-2015 standards, these ready-to-use broths require no pretreatment and fit a wide range of laboratory devices, ensuring consistent quality and compatibility.

Ordering information

Select the LBM® Enrichment and Selective Broths you need, in bulk packs

Cat. No.	Description	Pack size
4U161N Learn More	CABroth 2ml in Screw Cap Tube in Bulk  NEW	300 pieces 6 boxes of 50 pieces
474CE.A Learn More	BHI Broth 3ml in Screw Cap Tube in Bulk 	300 pieces 6 boxes of 50 pieces
4U002N Learn More	THIOL Broth 4ml in PET Tube in Bulk 	300 pieces 6 boxes of 50 pieces
477CE.A Learn More	TSB Broth plus 2.5% NaCl, 2ml in Screw Cap Tube in Bulk 	300 pieces 6 boxes of 50 pieces
475CE.A Learn More	SELENITE Broth 2ml in Screw Cap Tube in Bulk 	300 pieces 6 boxes of 50 pieces
476CE.A Learn More	LIM Broth 2ml in Screw Cap Tube in Bulk 	300 pieces 6 boxes of 50 pieces
085CU.A Learn More	4ml GN in Screw Cap Tube in Bulk 	300 pieces 6 boxes of 50 pieces

All above products are suitable for WASP® automation system. Please contact Copan for additional information.





These products may not be available in all countries. Please refer to www.copanusa.com for the most recent version of this document. Always consult product technical documentation for complete information. The use of this product in association with diagnostic kits or instrumentation should be validated by the end user. All trademarks, logos and brand names are the property of their respective owners. ©2026 COPAN Diagnostics Inc. All rights reserved.
Code: Enrichment_Broths_Brochure_CD_05_2026

COPANUSA.com [800.216.4016](tel:800.216.4016) [@COPANUSA](https://twitter.com/COPANUSA) CommunicationsUS@copangroup.com [/COPAN.USA](https://facebook.com/COPAN.USA)

COPAN Diagnostics, Inc.

26055 Jefferson Avenue | Murrieta, CA 92562 USA